
● CHROMagar™
COL-APSE

**For detection of Colistin resistant
gram-negative bacteria**




Plate Reading

-  • COL-R *E.coli*
→ Dark pink to reddish
-  • COL-R *Klebsiella*, *Enterobacter*, *Citrobacter*
→ Metallic blue
-  • COL-R *Pseudomonas*
→ translucent cream to green
-  • COL-R *Acinetobacter*
→ cream, opaque

For detection of Colistin resistant gram-negative bacteria

Background

Polymyxin E (colistin) and B are increasingly used as antimicrobials in the treatment of multi-drug resistant bacterial infections. Polymyxin resistance, although intrinsic in Gram-positive and some Gram-negative species (*Proteus*, *Morganella*, *Serratia*), is now a problem in a number of other pathogens (*Acinetobacter baumannii*, *Pseudomonas aeruginosa*, *Escherichia coli*, *Salmonella enterica*, *Klebsiella pneumoniae*).¹

Resistance arises due to mutations / insertions in genes involved in LPS biosynthesis (lpx, pmrA/B, mgrB, phoP/Q) and / or the acquisition of phosphoethanolamine transferases (PEtN). Of great concern is the recently described plasmid-encoded PEtN, MCR-1, now found worldwide in a range of animal, human and environmental bacterial isolates.¹

CHROMagar™ COL-APSE is a sensitive and specific medium for the growth of Colistin resistant bacterial pathogens with a lower limit of detection of 10 CFU/ml. This new chromogenic medium may be useful as a primary isolation medium in the surveillance and recovery of Colistin resistant bacteria from complex human, veterinary and environmental samples especially those with plasmid mediated MCR-1 or novel mechanisms of polymyxin resistance.¹

1- Novel Chromogenic Culture Media (CHROMagar™ COL-APSE) for the Isolation and Differentiation of Colistin Resistant Gram-negative Pathogens. Muhd Haziq F Abdul Momin, David W Wareham. ECCMID 2017.

Medium Performance

1 COLOURFUL DIFFERENTIATION OF COLONIES WITH ACQUIRED COLISTIN RESISTANCE

Distinguishing between *E.coli*, Coliforms, *Pseudomonas* and *Acinetobacter*.

2 IMPRESSIVE LIMIT OF DETECTION

10cfu/ml

Medium Description

Powder Base	Total	42.5 g/L
	Agar	15.0
	Peptones	20.0
	Salt	5.0
	Chromogenic and selective mix	0.8
	Growth factors	1.7
	Storage at 15/30 °C - pH: 7.1 +/-0.2	
	Shelf Life	2 years
+		
Supplement (included in the pack)	Liquid form.....	2 ml/L
	Storage at 15/30°C	
	Shelf Life	3 years

Usual Samples	stools, rectal swabs, veterinary, environmental samples
Procedure	Direct Streaking. Incubation 18-24h at 37°C Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Order References


Please use these product references when contacting your local distributor:

5000 ml pack CO262

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com

Find your nearest distributor on
www.CHROMagar.com/contact

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