
## CHROMagar™ Acinetobacter


Click below:

EN

FR

ES


# ● CHROMagar™ Acinetobacter

**For detection of *Acinetobacter*  
and MDR *Acinetobacter* species**

# ● CHROMagar™ Acinetobacter

www.CHROMagar.com


## Plate Reading

For detection of *Acinetobacter* spp.:

- *Acinetobacter* spp.  
→ red
- Other Gram (-)  
→ blue or mostly inhibited
- Gram (+) bacteria and yeasts  
→ inhibited

For detection of MDR *Acinetobacter* spp. (if using the optional supplement CR102):

- MDR *Acinetobacter*  
→ red
- Non-MDR *Acinetobacter*  
→ inhibited

## Medium Description

<b>Powder Base</b>	Total ..... 32.8 g/L Agar ..... 15.0 Peptone & Yeast extract ..... 12.0 Salt ..... 4.0 Chromogenic mix ..... 1.8 Storage at 15/30 °C - pH: 7.0 ± 0.2 Shelf Life ..... > 18 months
<b>Supplement</b> (included in the pack)	Growth and regulator factors ..... 4 mL/L Storage at 15/30 °C Aspect: Liquid Form Shelf Life ..... > 18 months
<b>CHROMagar™ MDR Selective Suppl.: CR102 and CR103-25</b> Order separately	Selective mix ..... 5 doses Storage at 2/8 °C Shelf Life ..... > 12 months

Usual Samples	rectal, nares and wounds swabs, stools, urine and surface samples.
Procedure	Direct Streaking. Incubation 18-24 h, 35-37 °C. Aerobic condition.

Scientific Publications on this product: available on [www.CHROMagar.com](http://www.CHROMagar.com)  
Please read carefully the instructions for use (IFU document) available on [www.CHROMagar.com](http://www.CHROMagar.com)

## For detection of *Acinetobacter* and MDR *Acinetobacter* spp.

### Background

Common bacteria widely spread in the nature, *Acinetobacter* has the capacity to survive in dry as well as moist environments. It becomes a source of infection in hospital environment when colonizing medical equipments, human skin and sometimes foodstuff. *Acinetobacter* species are generally not pathogenic for healthy people but are life threatening in compromised patients. It is often isolated in nosocomial infections cases, intensive care units, and can for instance cause nosocomial pneumonia, bacteraemia, and meningitis.

Especially, *Acinetobacter baumannii* is becoming a major hospital-acquired infection issue because of its often multi-drug resistance (MDR : resistance to C3G, quinolones, carbapenems, etc). This contributes to the increase of morbidity and mortality.

Active surveillance is necessary to control its spread in the facilities, to reduce the risk of cross-contamination, and to identify the carriers. Rapid identification of patients that are colonized with *Acinetobacter* would lead to infection control practices aimed at preventing spread of the organisms.

### Intended Use

CHROMagar™ *Acinetobacter* is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of colonization with *Acinetobacter* to aid in the prevention and control of *Acinetobacter*, drug-susceptible or multi-drug resistant (MDR), in healthcare settings. The test is performed with rectal swabs, nare swabs, wound swabs, stools and urine samples from patients to screen for *Acinetobacter* colonization. It can also be used in hygiene monitoring in the clinical environment with surface sampling. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

The medium can also be used as an early warning indicator for diagnostic tests of infections to signal the possible presence of multi drug-resistant bacteria. This use does not replace the institution's protocols.

CHROMagar™ *Acinetobacter* is not intended to diagnose *Acinetobacter* infection nor to guide nor monitor treatment for infections.

A lack of growth or the absence of colonies on CHROMagar™ *Acinetobacter* does not preclude the presence of *Acinetobacter*. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

### Medium Performance

1 **ONE UNIQUE RED COLOUR:** Detection of *A. baumannii* from traditional culture media might be a difficult and tedious task due to the abundance of background flora found in collected specimens, especially when using media based on differentiation by the lactose/non-lactose fermentation ability. To overcome these difficulties, CHROMagar™ *Acinetobacter* was designed as a highly selective medium, allowing the growth of *Acinetobacter* in conspicuously red colonies, after overnight incubation.

2 **UNIQUE chromogenic medium for *Acinetobacter* detection.**

3 **Screening of MDR *Acinetobacter*:** This medium can be supplemented to enhance MDR specificity allowing the growth of carbapenem-resistant strains.

	Analytical data	Clinical data
<b>Sensitivity</b>	100 % *	100 % ** compared to 46 % with Drigalski
<b>Specificity</b>	100 % *	99.9 % **

\* Data obtained after a 18-20 h incubation at 35-37 °C in aerobic conditions in the study «Laboratory evaluation of different agar media for isolation of carbapenem-resistant *Acinetobacter* spp.» Moran-Gilad *et al.*, 2014, *Eur. J. Clin. Microbiol. Infect. Dis.*

\*\* Data obtained after a 18 h incubation at 37 °C in aerobic conditions with rectal and nasal swabs from 1022 patients in the study «Overnight identification of imipenem-resistant *Acinetobacter baumannii* carriage in hospitalized patients». Gaillot *et al.*, Poster ICAAC 2010.


## Order References

Please use these product references when contacting your local distributor:

- 5000 mL pack ..... AC092 (included in this reference: powder base AC092(B) + supplement AC092(S))
- 25 L pack ..... AC093-25 (included in this reference: powder base AC093-25(B) + supplement 2xAC093-13(S))
- Optional:
- 5000 mL pack ..... CR102 (MDR Selective suppl.)
- 25 L pack ..... CR103-25 (MDR Selective suppl.)

Manufacturer: CHROMagar, 29 avenue George Sand, 93210 La Plaine Saint-Denis - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Website: [www.CHROMagar.com](http://www.CHROMagar.com)  
Find your nearest distributor on [www.CHROMagar.com/contact](http://www.CHROMagar.com/contact)

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