## CHROMagar™ *C.difficile*

Click below:

[EN](#)

[FR](#)

[ES](#)

# CHROMagar<sup>TM</sup> *C.difficile*

For detection of *Clostridioides difficile*


# CHROMagar™ C.*difficile*

For detection of *Clostridioides difficile*

## Plate Reading

- *C. difficile*  
→ colourless and fluorescent under UV light at 365 nm
- Other bacteria  
→ colourless, not fluorescent or inhibited

## Medium Description

<b>Powder Base</b>	Total ..... 54.7 g/L Agar ..... 15.0 Peptone & Yeast extract ..... 25.0 Salts ..... 9.0 Growth factors ..... 4.0 Chromogenic mix ..... 1.7 Storage at 15/30 °C - pH: 7.8 ± 0.2 Shelf Life ..... > 18 months
<b>Supplement</b> (Included in the pack)	Powder form ..... 325 mg/L Storage at 2/8 °C ..... Shelf Life ..... > 12 months
Usual Samples	Stools

## Order References

Please use these product references when contacting your local distributor:  
5000 mL pack ..... CD122  
(included in this reference: powder base CD122(B) + supplement CD122(S))

## Background

*Clostridioides difficile* is the leading cause of nosocomial infectious diarrhea in adults. These infections occur mostly in patients who have both medical care and antibiotic treatment. The symptoms of *C. difficile* infection are fever, abdominal cramps and severe diarrhea leading to death. In the United States, nearly 250,000 people each year develop *C. difficile* infections with at least 14 000 deaths (CDC estimate, 2013). Due to the emergence of highly toxigenic *C. difficile* strains, these infections have become more frequent and more difficult to treat in the last years.

Although PCR has become the leading *C. difficile* detection technique, culture is essential for strain typing and antimicrobial susceptibility testing. CHROMagar™ C.difficile is a fluorogenic culture medium, extremely sensitive and selective, especially designed to simplify and speed up (24 h) the culture of *C. difficile*.

## Intended Use

CHROMagar™ C.difficile is a selective and differential fluorogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with *Clostridioides difficile* to aid in the prevention and control in healthcare settings. The test is performed with stools. Results can be interpreted under UV light after 24 h of anaerobic incubation at 35-37 °C.

CHROMagar™ C.difficile is not intended to diagnose an infection nor to guide nor monitor treatment for infections. A lack of growth or the absence of fluorescent colonies on CHROMagar™ C.difficile does not preclude the presence of *C. difficile*. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

## Medium Performance

### 1 RAPID DETECTION COMPARED TO TRADITIONAL MEDIA

Big colonies (around 2 mm) of *C. difficile* after only 24 h of incubation in anaerobic atmosphere, contrary to traditional media requiring 48 h.

### 2 HIGH SENSITIVITY AND SPECIFICITY

	Analytical data *		Clinical data**	
	CHROMagar™ C. difficile	Reference medium (Taurocholate-CCFA)		
Sensitivity	-	95.4 %	70 %	
Specificity	100 %	88.8 %	97 %	

\* Data obtained after a 24 h incubation at 35 °C under anaerobic atmosphere in the study «Comparison of CHROMagar C.difficile and taurocholate-CCFA media for isolation of toxigenic *Clostridium difficile* from stools». CCFA = Cefoxitin-Cycloserine-Fructose Agar. Roux et al. Poster ASM 2014.

\*\* Data obtained after a 24 h incubation at 35 °C under anaerobic atmosphere with 594 stool samples, being positive 174, in the study «Comparison of CHROMagar C.difficile and taurocholate-CCFA media for isolation of toxigenic *Clostridium difficile* from stools». Roux et al. Poster ASM 2014.

Scientific Publications on this product: available on [www.CHROMagar.com](http://www.CHROMagar.com)  
Please read carefully the instructions for use (IFU document)  
available on [www.CHROMagar.com](http://www.CHROMagar.com)

Manufacturer: CHROMagar  
4 place du 18 juin 1940 75006 Paris - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Website: [www.CHROMagar.com](http://www.CHROMagar.com)  
Find your nearest distributor on  
[www.CHROMagar.com/contact](http://www.CHROMagar.com/contact)

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