● CHROMagar™ C.difficile

Click below:

● CHROMagar™
C.difficile

For detection of *Clostridium difficile*

● CHROMagar™ C.difficile

www.CHROMagar.com

For detection of *Clostridium difficile*

Background

Clostridium difficile is the leading cause of nosocomial infectious diarrhea in adults. These infections occur mostly in patients who have both medical care and antibiotic treatment. The symptoms of *C. difficile* infection are fever, abdominal cramps and severe diarrhea leading to death. In the United States, nearly 250,000 people each year develop *C. difficile* infections with at least 14 000 deaths (CDC estimate, 2013). Due to the emergence of highly toxigenic *C. difficile* strains, these infections have become more frequent and more difficult to treat in the last years.

Although PCR has become the leading *C. difficile* detection technique, culture is essential for strain typing and antimicrobial susceptibility testing. CHROMagar™ *C.difficile* is a fluorogenic culture medium, extremely sensitive and selective, especially designed to simplify and speed up (24 h) the culture of *C. difficile*.

Medium Performance

1 RAPID DETECTION COMPARED TO TRADITIONAL MEDIA
Big colonies (around 2 mm) of *C. difficile* after only 24 h of incubation in anaerobic atmosphere, contrary to traditional media requiring 48 h.

2 HIGH SENSITIVITY AND SPECIFICITY
C. difficile is detected by characteristic fluorescent colonies (under UV light at 365 nm) and the specimen's flora largely inhibited

Specificity / Sensitivity: ≈ 100 %*

* Sensitivity from scientific studies: «Comparison of CHROMagar™ *C.difficile* and taurochocolate-CCFA media for isolation of toxigenic *Clostridium difficile* from stools» Gaillot O. et al. ASM 2014.

3 POLYVALENCE
This medium can be used for clinical specimens as well as environmental samples.

Medium Description

Powder Base	Total	54.7 g/L
	Agar	15.0
	Peptone & Yeast extract	25.0
	Salts	9.0
	Growth factors	4.0
	Chromogenic mix.....	1.7
	Storage at 15/30 °C - pH: 7.8 ± 0.2	
	Shelf Life	> 18 months
Supplement (Included in the pack)	Powder form.....	325 mg/L
	Storage at 2/8 °C Shelf Life	> 12 months



Usual Samples	Stools
Procedure	Direct streaking or after an appropriate enrichment step of the sample. Incubation at 35-37 °C for 24 h. Anaerobic condition.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Plate Reading

- *C. difficile*
→ colourless and fluorescent under UV light at 365 nm
- Other bacteria
→ colourless, no fluorescent or inhibited

Order References

Please use these product references when contacting your local distributor:

5000 mL pack CD122
(included in this reference: powder base CD122(B) + supplement CD122(S))

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on
www.CHROMagar.com/contact

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