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

Figure 1: *Streptococcus agalactiae* Group B *Streptococcus* Figure 2: Results of a rapid agglutination test for Group B *Streptococcus*

CHROMagar

- CHROMagar StrepB: for detection of Group B *Streptococcus* (= *S. agalactiae*) by colony colour.
- *Streptococcus* B, A and *S. porcinus* (b) = mauve to pink
- other Gram + and other bacteria = blue, colourless or inhibited.

Purpose No. 1

Results No. 1 : Sensitivity, reading after one night incubation (<24h)

- Evaluate the technical sensitivity of CHROMagar StrepB media to detect *Streptococcus agalactiae* (Group B Strep).

- Genital samples to detect GBS (N=143) on pregnant women.
- Genital samples on non-pregnant women (N=178)

		CHROMagar StrepB			
N= 321		nég	Positive +	Positive ++	Positive +++
Blood agar plate	nég	Sc ≈ 97,4%			
	nég	277	4**	1**	
	Positive +	1*	10		
	Positive ++			13	
	Positive +++				15

Purpose No. 2

Results No. 2 : Specificity, reading after one night incubation (<24h)

- Evaluate the technical specificity of CHROMagar StrepB media.
- On a 66 collection strains of *Streptococcus* and *Enterococcus* spp.

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Figures

Figure 3 : Culture of Group B *Streptococcus* on CHROMagar StrepB medium


Figure 4: Culture of Group D *Streptococcus* on CHROMagar StrepB medium

Figure 5 : Comparison of culture: - left side : *Streptococcus agalactiae* (B) - right side : *Streptococcus pyogenes* (A)

Figure 6 : Culture of *Streptococcus porcinus* (b) on CHROMagar StrepB medium

Conclusion

We can conclude from this study that :

The CHROMagar StrepB medium is particularly rapid (<24h) for detecting *Streptococcus* involved in materno-faetal infections community associated or health-care associated (eventually nosocomial infections):

- *S. agalactiae* (B)
- *S. pyogenes* (A)
- *S. porcinus* (b)

For more information on CHROMagar StrepB :

CHROMagar Company:

Email: chromagar@chromagar.com

Other Data

References

- The reading is not influenced by incubation with or without CO₂.
- The suspect colonies have to be confirmed by serogrouping or identification (gallery or mass spectroscopy)
- Rapid Detection after overnight incubation, <24 hours.

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