CHROMagar™ Orientation

Click below:

EN

FR

ES

● CHROMagar™ Orientation

**For isolation and differentiation
of urinary tract pathogens**

For isolation and differentiation of urinary tract pathogens

Background

Urinary tract infections (UTIs) are a severe public health problem. Up to 40 % of women will develop UTI at least once during their lives, and a significant number of these women will have recurrent urinary tract infections. Thus, any workload reduction related to this analysis will dramatically improve the efficiency of the laboratory.

Intended Use

CHROMagar™ Orientation is a non-selective chromogenic culture medium intended for use in the qualitative direct detection, differentiation and presumptive identification of uropathogens to aid in the diagnosis of urine tract infections. The test is performed with urine samples. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

Concomitant cultures are necessary to recover organisms for further microbiological testing or epidemiological typing. A lack of growth or the absence of colonies on CHROMagar™ Orientation does not preclude the presence of bacteria. CHROMagar™ Orientation is not intended to diagnose infection nor to guide nor monitor treatment for infections.

CHROMagar™ Orientation can also be used with veterinary samples.

Medium Performance

1 INSTANT PALETTE OF COLOURS TO OBTAIN A LARGE SPECTRUM OF SPECIES DIFFERENTIATION

CHROMagar™ Orientation has several advantages over traditional media:

- allows in most cases full differentiation of the pathogens
- allows a semi-quantitative assessment and presumptive identification of urinary tract pathogens
- allows an easier recognition of mixed growth
- provides higher detection rates

2 HIGH DETECTION OF MINOR POPULATION

The proper use of CHROMagar™ Orientation will correctly pinpoint the presence of a minor population and will help to establish the right diagnosis and therapy.

3 SAVE TIME AND REDUCE WORKLOAD

The most common UTI pathogen is *E. coli*, found in 40-70 % of infections. CHROMagar™ Orientation has a specificity of 99,3 %* for *E. coli*. One plate of CHROMagar™ Orientation will give the same information as the combination of the 3 classical plates used for UTI analysis (blood agar, CLED and MacConkey agar). Moreover, since it is easy to differentiate mixed flora on CHROMagar™ Orientation, antimicrobial susceptibility tests can be performed directly from primary isolates without the need of subcultures**

	Analytical data *	Clinical data**
Sensitivity	100 %	100 %
Specificity	98 % For <i>E. coli</i> 99.3 %	-

* Data obtained after 16-24 h incubation at 35-37 °C in aerobic conditions with 1478 isolates were tested, being 74 enterococcal isolates and 1404 Gram-negative bacilli in the study «Evaluation of CHROMagar Orientation for Differentiation and Presumptive Identification of Gram-Negative Bacilli and Enterococcus Species». Merlino et al. 1996. *J. Clin. Microbiol.*

** Data obtained after 24 h incubation at 35 ± 2 °C in aerobic conditions with 900 urine samples in the study «Evaluation of Use of a New Chromogenic Agar in Detection of Urinary Tract Pathogens». Samra et al. 1998. *J. Clin. Microbiol.*

Medium Description

Powder Base	
Total	33 g/L
Agar	15.0
Peptone and yeast extract	17.0
Chromogenic mix	1.0
Storage at 15/30 °C - pH: 7.0 ± 0.2	
Shelf Life	> 18 months

Usual Samples	urine.
Procedure	Direct Streaking. Incubation at 35-37 °C, 18-24 h. Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Plate Reading

- *E. coli*
→ dark pink to reddish
- *Enterococcus*
→ turquoise blue
- *Proteus*
→ brown halo
- *Klebsiella, Enterobacter, Serratia*
→ metallic blue
- *S. aureus*
→ golden, opaque, small
- *Citrobacter*
→ metallic blue with red halo
- *S. saprophyticus*
→ pink, opaque, small
- *Candida albicans*
→ colourless
- *Streptococcus agalactiae*
→ Light blue
- *Pseudomonas aeruginosa*
→ Translucent, cream to blue

Order References

Please use these product references when contacting your local distributor:

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