
# Leaflet

LF-EXT-010

Version 8.1

## CHROMagar™ Salmonella

Click below:


EN

FR

ES

CHR  Magar™  
The Chromogenic Media Pioneer

# ● CHROMagar™ Salmonella

**For detection and isolation of *Salmonella***

## For detection and isolation of *Salmonella* species, including *S. Typhi* and *S. Paratyphi* in clinical specimens

### Background

Infections caused by *Salmonella* spp, including *Salmonella* Typhi, remain a major worldwide health problem:

- In the US, *Salmonella* has an incidence rate of 16.47 cases per 100,000 (CDC estimation, 2010).
  - In Europe, it is reported as the first bacterial cause of food outbreaks (EFSA/ECDC 2011 report, 2009 figures)
  - In developing countries, *Salmonella* Typhi and Paratyphi are commonly encountered with an estimated annual incidence of about 17 million cases. (2007 EFSA report)
- Moreover, according to a recent WHO report, *Salmonella* infections are responsible for 2 million deaths per year from diarrhoea. *Salmonella* is the second most reported zoonotic infection in humans (EFSA/ECDC 2011 report, 2009 figures).

Mainly due to contamination in the food chain and/or during food-production processes, *Salmonella* commonly induces enteric illness whose major symptoms are abdominal cramps, diarrhea, nausea, vomiting. More severe cases, for instance typhoid cases or infections in immuno-depressed patients, can lead to body dehydration with renal failure or bacteraemia.

### Intended Use

CHROMagar™ Salmonella is a selective chromogenic culture medium intended for use in the qualitative direct detection, differentiation, and presumptive identification of *Salmonella*. The test is performed with rectal swabs and stools, to aid in the diagnosis of *Salmonella* infections. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.


Concomitant cultures are necessary to recover organisms for further microbiological testing or epidemiological typing. A lack of growth or the absence of mauve colonies on CHROMagar™ Salmonella does not preclude the presence of *Salmonella*.

CHROMagar™ Salmonella is not intended to diagnose infection nor to guide nor monitor treatment for infections.

CHROMagar™ Salmonella can also be used in the detection of *Salmonella* in the analyses of food products for human consumption, animal feed and in environmental samples.


### Medium Performance

- EASY READING**  
intense mauve colony colours for better identification and partial inhibition of *E. coli* and coliforms.
- GREATER SPECIFICITY / LESS WORKLOAD**  
conventional media for the detection of *Salmonella* by H<sub>2</sub>S character have very poor specificity resulting in numerous false positives (*Citrobacter*, *Proteus*, etc.) among the rare, real positive *Salmonella*. The workload for unnecessary examination of suspect colonies is so heavy that real positive *Salmonella* colonies might often be overlooked in routine testing. Because of their poor specificity, conventional media require a tedious examination of at least 10 colonies per suspected sample. On the contrary, CHROMagar™ Salmonella eliminates most of those false positives and allows technicians to focus on the real contaminated samples.
- HIGH SENSITIVITY AND SPECIFICITY**  
leading to a higher detection rate of *Salmonella*
- DRAMATIC REDUCTION OF THE WORKLOAD**  
Number of useless confirmatory tests is minimized since there is no need of duplicating them.

### Plate Reading

- *Salmonella* including *S. Typhi*  
→ mauve
- Other bacteria  
→ blue, colourless or inhibited

### Medium Description

Powder Base	Total .....	34.9 g/L
	Agar .....	15.0
	Peptone & Yeast extract .....	7.0
	Chromogenic and selective mix .....	12.9
	Storage at 2/30 °C - pH: 7.6 ± 0.2	
Shelf Life .....	> 12 months	

Usual Samples	Rectal swabs and stools
Procedure	Direct Streaking. Incubation 18-24 h, 35-37 °C. Aerobic condition.

Scientific Publications on this product: available on [www.CHROMagar.com](http://www.CHROMagar.com)  
Please read carefully the instructions for use (IFU document) available on [www.CHROMagar.com](http://www.CHROMagar.com)

	Analytical data	Clinical data
Sensitivity	(81 %) and 93 % *	95 % **
Specificity	100 % *	88.9 % ** compared to 78.5 % with Hektoen agar

\* In-house data obtained after a 24-48 h incubation at 37 °C in aerobic conditions. Sensitivity % in parenthesis includes lactose positive *Salmonella* species growing in blue. 2012.

\*\* Data obtained after a 18-24 h incubation at 37 °C in aerobic conditions with 508 stool samples analyzed in the study "Comparison of CHROMagar™ Salmonella medium and Hektoen Enteric Agar for isolation of Salmonellae from stool samples". Gaillot *et al.*, 1998. *J. Clin. Microbiol.*


### Order References

Please use these product references when contacting your local distributor:

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