

CHROMagar™ Malassezia

Click below:

CHROMagar™
Malassezia

For detection of *Malassezia* spp.

CHRO Magar™
The Chromogenic Media Pioneer



Plate Reading

- *Malassezia furfur*
→ large, pale pink and wrinkled
- Other *Malassezia* spp including *M. globosa* & *M. restricta*
→ mostly pink to purple

For detection of *Malassezia* spp.

Background

Malassezia is a fungi naturally found on animal skins. On normal healthy skin it does not cause infections, but when the environment of the skin is altered, *Malassezia* species are able to cause several cutaneous diseases as severe dermatitis or otitis (inflammation of the skin or ears respectively).

1 Medium Performance

EASY AND COLOURFUL DETECTION

Since members of the genus *Malassezia* share similar morphological and biochemical characteristics, the use of traditional culture media for differentiating them based on phenotypic features is not suitable. CHROMagar™ Malassezia was developed with the goal of facilitating not only their detection, but also to improve an algorithm for the differentiation of the most common species (see: Revised Culture-Based System for Identification of *Malassezia* Species, by Takamasa et al. JCM No-2007)

Typical Appearance of the colonies reported in the study:

M. pachydermatis CBS 1879 → large, pale pink & smooth
M. restricta CBS 7877 → small, pink & smooth
M. dermatis JCM11348 and JCM11470 → small, pale pink & smooth
M. slooffiae CBS 7956 → large, pale pink & smooth
M. obtusa CBS 7876 → medium, pink & rough
M. globosa CBS 7966 → small, purple & smooth
M. sympodialis CBS 7222 → large, pale pink & smooth
M. furfur CBS 1878 → large, pale pink & wrinkled

Extracted from: «Revised Culture-Based System for Identification of *Malassezia* Species», by Takamasa et al. JCM No-2007)

Medium Description

Powder Base	Total 56.3 g/L Agar 15.0 Peptones and extracts 38.0 Chromogenic mix 2.8 Chloramphenicol 0.5 Storage at 15/30 °C - pH: 6.3 +/- 0.3 Shelf Life > 18 months
Mandatory Supplements to the Powder Base (Not provided by CHROMagar™):	
Glycerol	2 g
Tween 40	10 g

Usual Samples	Veterinary samples (skin, ear canals etc.)
Procedure	Direct Streaking. Incubation 72 h at 30-37 °C Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available
on www.CHROMagar.com

Order References


Please use these product references when contacting your local distributor:

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