CHROMagar™ C3G®

C3G R

Click below:

EN

FR

ES

CHRO Magar™
The Chromogenic Media Pioneer

CHROMagar™
C3G®

C3G® R

For overnight detection of beta-lactamase producing
Enterobacteria achieving third generation cephalosporin
resistance

CHRO Magar™
The Chromogenic Media Pioneer


Plate Reading

- *E. coli*
→ dark pink to reddish
- *Klebsiella, Enterobacter, Citrobacter*
→ metallic blue (+/- red halo)
- *Proteus*
→ brown halo
- *Pseudomonas*
→ translucent cream to blue
- *Acinetobacter*
→ cream, opaque

For overnight detection of beta-lactamase producing Enterobacteria achieving third generation cephalosporin resistance

Background

β -Lactamase production (ESBL, AmpC, carbapenemase) is the most common mechanism of β -lactam drug resistance in Gram negative bacteria. Each of the involved mechanisms can be plasmid-encoded and transfer horizontally. Moreover, plasmids continuously gather other antibiotic resistance determinants, ultimately closing to pan-resistance. Their spread leads to increased colonization pressure in hospital and community as well. It is crucial to ensure that proper surveillance is in place to help implement appropriate guidelines and policies for their spread control. Rapid detection of bacteria producing these enzymes also allows for de-escalation to more targeted therapy, to conserve carbapenem antibiotics for more serious infections.

Intended Use

CHROMagar™ C3G® is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with 3rd generation cephalosporin-resistant Enterobacteria (C3G®-E) to aid in the prevention and control of C3G®-E in healthcare settings. The test is performed with rectal swab and stools from patients to screen for C3G®-E colonization. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

The medium can also be used as an early warning indicator for diagnostic tests of infections to signal the possible presence of multi drug-resistant bacteria. This use does not replace the institution's protocols. CHROMagar™ C3G® is not intended to diagnose C3G®-E infection nor to guide nor monitor treatment for infections. A lack of growth or the absence of colonies on CHROMagar™ C3G® does not preclude the presence of C3G®-E. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

Medium Performance

1 FAST RESULTS

Detection after overnight incubation.

2 SPECIES DIFFERENTIATION

Thanks to the chromogenic performances of supplemented CHROMagar™ Orientation. Indeed, the product is composed of a powder base CHROMagar™ Orientation and a supplement to select β -Lactamase producing bacteria.

3 HIGH SENSITIVITY

Unique medium not inhibiting plasmid-mediated AmpC-producing bacteria.

4 TIME AND WORKLOAD SAVINGS

Direct culture from specimen. There is no need of a selective pre-enrichment.

	Analytical data *	Clinical data**	
		CHROMagar™ C3G®	Reference medium (ESBL Isolation Agar Non chromogenic media)
Powder Base (CHROMagar™ Orientation)	Total 33 g/L Agar 15.0 Peptone and yeast extract 17.0 Chromogenic mix 1.0 Storage at 15/30 °C - pH: 7.0 +/- 0.2 Shelf Life > 18 months		
CHROMagar™ C3G® Supplement (included in the pack)	Selective mix (Powder form) 0.37 g/L Storage at 2/8 °C Shelf Life > 18 months	Sensitivity 100 % Specificity -	77 % - 67 % -

* Data obtained after 24 h incubation at 37 °C in aerobic conditions in the study «Evaluation of Four Media for the Isolation of Resistant Gram Negative Organisms from Surveillance Rectal Swabs». Joshi-Caeser et al. Poster ASM 2012.

** Data obtained after 24 h incubation at 37 °C in aerobic conditions with 96 rectal samples in the study «Evaluation of Four Media for the Isolation of Resistant Gram Negative Organisms from Surveillance Rectal Swabs». Joshi-Caeser et al. Poster ASM 2012.

Medium Description

Powder Base (CHROMagar™ Orientation)	Total 33 g/L Agar 15.0 Peptone and yeast extract 17.0 Chromogenic mix 1.0 Storage at 15/30 °C - pH: 7.0 +/- 0.2 Shelf Life > 18 months
CHROMagar™ C3G® Supplement (included in the pack)	Selective mix (Powder form) 0.37 g/L Storage at 2/8 °C Shelf Life > 18 months

Usual Samples stools and rectal specimens.

Procedure Direct Streaking. Incubation 18-24 h at 35-37 °C.
Aerobic conditions

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on
www.CHROMagar.com

Order References

Please use these product references when contacting your local distributor:

5000 mL pack CGRT2

(included in this reference: powder base RT412 + supplement CG632)

25 L pack CGRT3-25

(included in this reference: powder base RT413-25 + supplement CG633-25)

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
CHR MagarTM
The Chromogenic Media Pioneer

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Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

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