

CHROMagarTM Orientation

Click below:

CHROMagarTM Orientation


For isolation and differentiation
of urinary tract pathogens

CHROMagar™ Orientation

www.CHROMagar.com

Plate Reading

For isolation and differentiation of urinary tract pathogens

Background

Urine analysis is the most common clinical microbial test.

For instance, in France in 2007, out of 10 million bacteriology tests carried out, 6 million (60 %) were urinalyses. Thus, any workload reduction related to this analysis will dramatically improve the efficiency of the laboratory.

Medium Performance

1 INSTANT PALETTE OF COLOURS TO OBTAIN A LARGE SPECTRUM OF SPECIES DIFFERENTIATION

CHROMagar™ Orientation has several advantages over traditional media:

- allows in most cases full differentiation of the pathogens
- allows for reliable detection semi-quantitative and presumptive identification of urinary tract pathogens
- easier recognition of mixed growth
- provides higher detection rates

2 HIGH DETECTION OF MINOR POPULATION

The proper use of CHROMagar™ Orientation will correctly pinpoint the presence of a minor population and will help to establish the right diagnosis and therapy.

3 SAVE TIME AND REDUCE WORKLOAD

The most common UTI pathogen is *E. coli*, found in 40-70 % of infections.

CHROMagar™ Orientation has a specificity of **99,3 %*** for *E. coli*, rendering the species confirmatory test largely unnecessary.

*Merlino, J. et al. 1996. Evaluation of CHROMagar™ Orientation for Differentiation and Presumptive Identification of Gram-Negative Bacilli and Enterococcus Species, J.C.M. 34: 1788-1793.

One plate of CHROMagar™ Orientation will give the same information as the combination of the 3 classical plates used for UTI analysis (blood agar, CLED and MacConkey agar). Moreover, since it is easy to differentiate mixed flora on CHROMagar™ Orientation, antimicrobial susceptibility tests can be performed directly from primary isolates without the need of subcultures*.

*Samra, Z et al. 1998. Evaluation of use of a new chromogenic agar in detection of urinary tract pathogens, J.C.M. 36: 990-994

4 ISOLATION OF A VARIETY OF MICROORGANISMS

The major target of this medium is the detection of urinary tract pathogens but CHROMagar™ Orientation has a broader application when supplemented with various antibiotics in detecting increasingly important nosocomial and multidrug resistant microorganisms (See CHROMagar™ ESBL and CHROMagar™ KPC).

CHROMagar™ can also be used to differentiate various microorganisms in environmental field.

Medium Description

Powder Base	Total 33 g/L Agar 15.0 Peptone and yeast extract 17.0 Chromogenic mix 1.0 Storage at 15/30 °C - pH: 7.0 ± 0.2 Shelf Life > 18 months
--------------------	---

Usual Samples	urine, surfaces, air, clinical and other materials.
Procedure	Direct Streaking, Incubation at 35-37 °C, 18-24 h. Aerobic condition.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Order References

Please use these product references when contacting your local distributor:

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

CHRO MagarTM
The Chromogenic Media Pioneer

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