● CHROMagar™ Staph aureus

Click below:

● CHROMagar™ Staph aureus

**For isolation and direct differentiation
of *Staphylococcus aureus***

For isolation and direct differentiation of *Staphylococcus aureus* in clinical and industrial samples

Background

Food Industry: Human beings are the main reservoir of *S. aureus*. A carrier contaminates the surrounding environment when coughing, sneezing and by touching food with a hand. It is often found in the environment and on food preparation surfaces and also in certain uncooked foods (dairy products, salads, sandwiches...). It is important to check the presence of *S. aureus* before and after the foodstuff sterilisation process.

Clinical relevance: *S. aureus* is the leading cause of skin and soft tissue infections and can also cause serious infections such as bloodstream infections, pneumonia, or bone and joint infections.

Medium Performance

Clinical application

1 EASY TO READ
compared to Blood Agar or Mannitol Salt Agar. CHROMagar™ Staph aureus allows easier differentiation of *S. aureus* colonies enhanced by a mauve colour and is of considerable help in identifying suspect colonies. Thus, it reduces the confirmatory workload.

2 HIGH SPECIFICITY **exceeds 99 %***

*Specificity from scientific study: "Evaluation of CHROMagar Staph aureus, a new chromogenic medium, for isolation and presumptive identification of *Staphylococcus aureus* from human clinical specimens." Gaillot O. et al. ASM 2000.

Food and environmental quality control

1 EASY TO PREPARE
The conventional medium for *S. aureus* is the Baird-Parker which has to be supplemented with RPF (Rabbit Plasma Fibrinogen), rendering the plate manufacturing delicate and complex, and also reducing the shelf life of the poured plates to a couple of weeks. On the contrary, CHROMagar™ Staph aureus comes with all the compounds already in the agar (no need of any supplement) and remains stable.

2 FAST
The results on Baird Parker have to be read after 48 h of incubation while with CHROMagar™ Staph aureus the results are available after only 24 h.

Medium Description

| | | |
|--------------------|--------------------------------------|-------------|
| Powder Base | Total | 82.5 g/L |
| | Agar | 15.0 |
| | Peptone and yeast extract | 40.0 |
| | Salts | 25.0 |
| | Chromogenic mix | 2.5 |
| | Storage at 15/30 °C - pH 6.9 +/- 0.2 | |
| | Shelf Life | > 18 months |

| | |
|---------------|--|
| Usual Samples | Clinical: stools, nasal swab, skin, sputum, throat and wounds. Industrial: Food stuff |
| Procedure | Direct streaking. Incubate at 35-37 °C for 18-24 h. Aerobic conditions. |

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com



Plate Reading

- *S. aureus*
→ pink to mauve
- Other bacteria
→ colourless, blue or inhibited

Order References

Please use these product references when contacting your local distributor:

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