

● CHROMagar™ MRSA

Click below:

● CHROMagar™ MRSA


For isolation and differentiation of
Methicillin Resistant *Staphylococcus aureus*
(MRSA) including low level MRSA

Plate Reading

- Methicillin Resistant *Staphylococcus aureus* (MRSA)
→ rose to mauve
- Methicillin Susceptible *Staphylococcus aureus* (MSSA)
→ inhibited
- Other bacteria
→ blue, colourless or inhibited

For isolation and differentiation of Methicillin Resistant *Staphylococcus aureus* (MRSA) including low level MRSA

Background

Leading cause of nosocomial infections, especially in intensive care units, the MRSA sources are either endogenous (the patient) or through cross contamination (environmental or by person to person contact).

The major issue with this pathogen is its resistance to a large panel of antibiotics, among them beta-lactam antibiotics, limiting the therapeutic options for clinicians.

Early detection is essential for controlling the spread of MRSA, providing appropriate care, and avoiding complex and expensive treatments. Pre-admission screening for MRSA has proved to be an effective method for reducing the hospital burden of MRSA-colonised patients. The savings due to consistent decolonisation before elective admission outweigh the costs of screening. Today, in the US, the extra-expenses linked to difficult treatments of MRSA infections are estimated at \$2.4 billion for about 370,000 hospital stays. (Genetic Engineering and Biotechnology News, August 2009).

In the UK, the estimation of the additional cost of discharging every hospital patient who acquires MRSA is £9,000.

Medium Performance

1

ABSOLUTELY RELIABLE

CHROMagar™ MRSA, introduced in 2002, was the first chromogenic medium for MRSA detection. It led to such significant reductions in both, the response time and laboratory workload, that it allowed an absolutely necessary wide-scale patient screening.

2

EFFICIENT

The medium exhibits sensitivity and specificity values close to 100 %. CHROMagar™ MRSA allows an accurate detection of MRSA with a higher level of sensitivity than oxacillin containing media.

3

FAST & EASY INTERPRETATION

Intense mauve colony colour in 18-24 h.

Medium Description

Powder Base	Total	82.5 g/L
	Agar	15.0
	Peptones and yeast extract.....	40.0
	Salts	25.0
	Chromogenic mix	2.5
	Storage at 15/30 °C - pH: 6.9 +/- 0.2	
+ Lyophilized supplement (included in the pack)	Shelf Life	> 18 months
	Powder form QSF 20 L	20 mL
	Storage at 2/8 °C	Shelf Life > 18 months

Usual Samples	nasal, perineal, throat, rectal specimens
Procedure	Direct Streaking. Incubation 18-24 h at 35-37 °C. Aerobic conditions
Scientific Publications on this product: available on www.CHROMagar.com Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com	


Order References

Please use these product references when contacting your local distributor:

5000 mL pack MR502 (included in this reference: powder base MR502 + supplement SU620)

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on
www.CHROMagar.com/contact

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