

● CHROMagar™ STEC

Click below:

● CHROMagar™ STEC


For detection of
Shiga-Toxin producing *E. coli* (STEC)

Plate Reading

- Most common Shiga-Toxin *E. coli* serotypes
→ mauve
- Other Enterobacteriaceae
→ colourless, blue or inhibited
- Gram (+) bacteria
→ inhibited

DETECTION OF THE BIG SIX

Serotype	isolates	Sensitivity
O157	394/410	96 %
O26	152/155	98 %
O45	1/1	100 %
O103	62/90	69 %
O111	64/66	97 %
O121	31/36	86 %
O145	45/49	92 %

Internal figures

More Scientific content
available online

For detection of Shiga-Toxin producing *E. coli* (STEC)

Background

An increasing and worrisome number of studies have lately shown that, non-O157 ShigaToxin-producing *E. coli* (STEC) have been responsible for foodborne poisoning outbreaks. The CDC has also reported warnings about this potential risk:

“ Disease caused by Shiga toxin-producing *Escherichia coli* (STEC) ranges from self-limiting diarrhea to hemorrhagic colitis and hemolytic uremic syndrome (HUS). Serotype O157:H7, the most frequently implicated STEC causing hemorrhagic colitis and HUS, has been isolated from large foodborne outbreaks, as well as sporadic cases, in North America and abroad. However, 60 STEC serotypes have been implicated in diarrheal disease, and several non-O157:H7 serotypes have been implicated as the cause of foodborne outbreaks and HUS in the United States, Europe, and Australia. Studies from Canada, Europe, Argentina, and Australia suggest that non-O157:H7 STEC infections are as prevalent, or more so, than O157:H7 infection.”

CDC report « Prevalence of Non-O157:H7 Shiga Toxin-Producing *Escherichia coli* in Diarrheal Stool Samples [...] »

Therefore, several regulatory authorities urge worldwide food industry to implement measures to control the absence of such organisms in their production. In many cases, laboratories have limited their search for pathogenic *E. coli* to the common O157 serotype. This is due, among other reasons, to the fact that there were no available selective culture media for non-O157 *E. coli*. CHROMagar™ STEC is designed to fill this gap: detection, as mauve colonies, of not only the classical STEC O157, but also many other serotypes.

Medium Performance

- 1 **Easy reading :** A majority of STEC strains grow in mauve colony color, while other bacteria grow in blue, colourless or are inhibited.
- 2 **Highly STEC selective medium:** Excellent tool for large number of samples screening procedures.
- 3 **Worldwide premiere:** Unique medium in the market for STEC detection.
- 4 **Flexibility:** It can be supplemented with additional compounds to render it even more selective for the strain causing an outbreak.

Medium Description

Powder Base CHROMagar™ STEC base	Total 30.8 g/L Agar 15.0 Peptone and yeast extract 8.0 Salts 5.2 Chromogenic mix 2.6 Storage at 15/30 °C - pH: 6.9 +/-0.2 Shelf Life > 18 months
CHROMagar™ STEC Supplement (included in the pack)	Freeze-dried vials 5 x 1000 mL Storage at 15/30 °C Shelf Life > 18 months

Usual Samples	Food, environmental, faeces.
Procedure	Direct Streaking. Incubation 18-24 h at 35-37 °C. +/- prior enrichment step. Aerobic conditions.
Scientific Publications on this product: available on www.CHROMagar.com Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com	


Order References

Please use these product references when contacting your local distributor:

5000 mL pack ST162
(included in this reference: powder base base ST162(B) + 5 supplements ST162(S))

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on
www.CHROMagar.com/contact

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