
# Leaflet

LF-EXT-050

Version 4.0

## ● CHROMagar™ Streptococcus

Click below:


[EN](#)

[FR](#)

[ES](#)

CHR  Magar™  
The Chromogenic Media Pioneer

# ● CHROMagar™ Streptococcus

**For isolation and differentiation  
of *Streptococcus* spp.**

## For isolation and differentiation of *Streptococcus* spp.

### Background

Environmental Streptococci are considered to be major mastitis pathogens because they typically cause high somatic cell counts, persistent infections and milk bacterial contamination. CHROMagar™ Streptococcus was designed to help in the Mastitis diagnosis, particularly as an aid to differentiate species involved in the udder infection.


There are other applications where CHROMagar™ Streptococcus has also shown to be a useful tool to differentiate Streptococci in samples heavily loaded with a complex bacterial flora.

### Medium Performance

- 1 **HIGHLY SELECTIVE**  
for Streptococci group.
- 2 **COLOR DIFFERENTIATION**  
Allows a clear color differentiation of various Streptococci species.
- 3 **FAST RESULT**  
Results within 24 h of incubation.
- 4 **EASY INTERPRETATION**  
CHROMagar™ Streptococcus is suitable for specimens heavily loaded with microbial flora.
- 5 Possibility to carry additional tests directly from a colony.

### Plate Reading

- *S. agalactiae* → blue
- *Enterococcus* → mauve to violet-mauve
- *S. uberis* → blue to metallic blue
- Other *Streptococcus* → various

### Medium Description

<b>Powder Base</b>	Total ..... 44.9 g/L Agar ..... 15.0 Peptones and growth factors ..... 20.0 Salts ..... 7.5 Chromogenic and selective mix ..... 2.4 Storage at 15/30 °C - pH: 7.3 +/- 0.2 Shelf Life ..... > 18 months
<b>+</b> <b>Supplement</b> (included in the pack)	Growth factors ..... 8.0 mL/L Storage at 15/30 °C Aspect: Liquid Form Shelf Life ..... > 18 months

Usual Samples	Milk
Procedure	Direct Streaking. Incubation 24 h at 37 °C. Aerobic conditions.
Scientific Publications on this product: available on <a href="http://www.CHROMagar.com">www.CHROMagar.com</a> Please read carefully the instructions for use (IFU document) available on <a href="http://www.CHROMagar.com">www.CHROMagar.com</a>	

### Order References


Please use these product references when contacting your local distributor:

500 g pack .....CQ39-500  
(included in this reference: powder base CQ39B-500G + supplement CQ39S-05)

Manufacturer: CHROMagar, 29 avenue George Sand, 93210 La Plaine Saint-Denis - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Website: [www.CHROMagar.com](http://www.CHROMagar.com)

Find your nearest distributor on [www.CHROMagar.com/contact](http://www.CHROMagar.com/contact)

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