Rambach™ Agar

Instructions For Use

Available in several languages

NT-EXT-006

Version **10.0**

Click below for:

ENGLISH


FRANCAIS

ESPAÑOL

DEUTSCH

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Rambach™ Agar plate

MEDIUM PURPOSE

Chromogenic medium for detection and isolation of *Salmonella* species in food samples.

Infections caused by *Salmonella* spp, including *Salmonella typhi*, remain a major worldwide health problem:

- In the US, *Salmonella* has an incidence rate of 16.47 cases per 100,000 (CDC estimation, 2010).
- In Europe, it is reported as the first cause of food outbreaks (EFSA/ECDC 2011 report, 2009 figures).
- In developing countries, *Salmonella typhi* and *paratyphi* are commonly encountered with an estimated annual incidence of about 17 million cases (2007 EFSA report).

Moreover, according to a recent WHO report, *Salmonella* infections are responsible for 2 million deaths per year from diarrhea. *Salmonella* is the second most reported zoonotic infection in humans (EFSA/ECDC 2011 report, 2009 figures).

COMPOSITION

The product is composed of a powder base (B) and one liquid supplement (S).

Product	=	Base (B)	Supplement (S)
Total g/L		30.7 g/L	10 mL/L
Composition g/L		Opaque Agar 20.0 Peptone and yeast extracts 8.0 Chromogenic mix 2.7	Propylene glycol 10.4
Aspect		Powder Form	Liquid Form
STORAGE		15/30 °C	15/30 °C
FINAL MEDIA pH		7.1 +/- 0.2	

PREPARATION (Calculation for 1 L)

Step 1

Preparation

- Disperse slowly 30.7 g of powder base in 1 L of purified water.
- Stir until agar is well thickened.
- Add 10 mL of the LIQUID content of the supplement vial.
- Swirl for mixing.
- Heat and bring to boil (100 °C) while swirling or stirring regularly.

DO NOT HEAT TO MORE THAN 100 °C. DO NOT AUTOCLAVE AT 121 °C.

Warning 1: If using an autoclave, do so without pressure.

Advice 1: For the 100 °C heating step, mixture may also be brought to a boil in a microwave oven: after initial boiling, remove from oven, stir gently, then return to oven for short repeated bursts of heating until complete fusion of the agar grains has taken place (large bubbles replacing foam).

Final media HELPING CALCULATION (Supplement)

1 L	Add 10 mL of supplement to the prepared base
5 L	Add 50 mL of supplement to the prepared base
25 L	Add 250 mL of supplement to the prepared base

Step 2

Pour plates

- Cool down in a water bath to 45-50 °C.
- Swirl or stir gently to homogenize.
- Pour into sterile Petri dishes.
- Let it solidify and dry.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to two months under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.

INOCULATION

- Related samples can be processed by direct streaking on the plate, as well as prior appropriate enrichment step.
- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate and incubate at 37 °C for 24 hours in aerobic condition.

Typical Samples


e.g. Food and environmental samples

Possible enrichment step
Direct streaking or spreading technique

INTERPRETATION

Microorganism	Typical colony appearance
<i>Salmonella</i>	→ red
Many coliforms	→ blue or violet
<i>Proteus</i> , etc.	→ colourless
Gram positive	→ inhibited

Typical colony appearance

The pictures shown are not contractual.

PERFORMANCE & LIMITATIONS

- Various studies have shown that 97 % to 99 % *Salmonella* strains give the typical red colour colonies excluding *S. paratyphi* A and *S. typhi* (colourless colonies or no colonies), and Lactose positive *Salmonella* (blue colonies).
- This selective medium is inhibitory for many microorganisms other than *Salmonella*.
- Rare strains of *Pseudomonas* can appear positive and can be eliminated by an Oxydase test.
- Final identification must be done by biochemistry and serology (e.g. Latex agglutination test from Microgen), and can be performed directly from the plates on suspected colonies.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the ATCC strains below:

Microorganism	Typical colony appearance
<i>S. enteritidis</i> ATCC® 13076	→ red
<i>S. abaeetuba</i> ATCC® 35640	→ red
<i>E. coli</i> ATCC® 8739	→ metallic blue
<i>C. freundii</i> ATCC® 8090	→ purple blue
<i>S. aureus</i> ATCC® 25923	→ inhibited

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by propriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

- Catalogue reference
- Consult instructions for use
- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity
- Protect from light
- Manufacturer

Need some Technical Documents?

Available for download on www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

Pack Size

5000 mL

250 Tests of 20mL

=

Ordering References

RR702

=

Base

RR702B
Weight: 153.5 g

+

Supplement

RR702S
Quantity: 50 mL

25 L

1250 Tests of 20mL

=

RR703-25

=

RR703-25B
Weight: 767.5 g

+

RR703-25S
Quantity: 250 mL

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach
ATCC® is a registered trademark of the American Type Culture Collection
NT-EXT-006 V10.0 / EN 26-Nov-18

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