CHROMagar™ **Vibrio**

Instructions For Use

Available in several languages

NT-EXT-017

Version **8.1**

ENGLISH

English Version

FRANCAIS


Version Française

ESPAÑOL

Version Español

DEUTSCH

Deutsch Version

MEDIUM PURPOSE

Chromogenic medium for isolation and detection of *V. parahaemolyticus*, *V. vulnificus* and *V. cholerae*.

Naturally present on marine plants and animals, *Vibrio* genus counts over 20 species among which four (*V. parahaemolyticus*, *V. vulnificus*, *V. alginolyticus* and *V. cholerae*) represent a serious public health hazard.

COMPOSITION

The product is composed of a single powder base.

| Product | = | Pack |
|-----------------|---|---|
| Total g/L | | 74.7 g/L |
| Composition g/L | | Agar 15.0 Peptone and yeast extract 8.0 Salts 51.4 Chromogenic mix 0.3 |
| Aspect | | Powder Form |
| STORAGE | | 15/30 °C |
| FINAL MEDIA pH | | 9.0 +/- 0.2 |

PREPARATION (Calculation for 1 L)

Step 1

Preparation of the mix

- Disperse slowly 74.7 g of powder base in 1 L of purified water.
 - Stir until agar is well thickened.
 - Heat and bring to boil (100 °C) while swirling or stirring regularly.
- DO NOT HEAT TO MORE THAN 100 °C. DO NOT AUTOCLAVE AT 121 °C.

Warning 1: If using an autoclave, do so without pressure.

Advice 1: For the 100 °C heating step, mixture may also be brought to a boil in a microwave oven: after initial boiling, remove from oven, stir gently, then return to oven for short repeated bursts of heating until complete fusion of the agar grains has taken place (large bubbles replacing foam).

Step 2

Pouring

- Cool in a water bath at 45-50 °C, swirling or stirring gently.
- Advice 2: In case of samples with a high presence of *Aeromonas*, 50 mg of cefsulodin can be added to the mix once cooled down at 45-50 °C (50 mg/L).**
- Pour medium into sterile Petri dishes.
 - Let it dry and gel.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to 1 month under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.

INOCULATION

Related samples can be processed by direct streaking on the plate, as well as prior appropriate enrichment step.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate.
- Incubate at 37 °C for 24 h in aerobic conditions.

Typical Samples


e.g. environmental, water,
seafood, surfaces

Direct streaking
or spreading technique

INTERPRETATION

| Microorganism | Typical colony appearance |
|---|--------------------------------|
| <i>V. parahaemolyticus</i> | → mauve |
| <i>V. vulnificus</i> / <i>V. cholerae</i> | → green blue to turquoise blue |
| <i>V. alginolyticus</i> | → creamy |

Typical colony appearance

The pictures shown are not contractual.

PERFORMANCE & LIMITATIONS

- The final identification must be confirmed by biochemical tests or by mass spectrophotometry (eg. MALDI-TOF). They can be done directly from the suspicious colonies observed on the medium.
- For the Oxydase test of blue colonies, we suggest the use of normal Oxydase test.
- For the Oxydase test of mauve colonies, we suggest the use of a reagent giving a blue colour with oxydase positive bacteria (tetramethyl-p-phenylenediamine solution at 10 mg/mL).

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the ATCC strains below:

| Microorganism | Typical colony appearance |
|--|---------------------------|
| <i>V. vulnificus</i> ATCC® 27562 * | → green |
| <i>V. parahaemolyticus</i> ATCC® 33845 | → mauve |
| <i>V. alginolyticus</i> ATCC® 33839 | → creamy |
| <i>S. aureus</i> ATCC® 25923 | → inhibited |
| <i>E. coli</i> ATCC® 25922 | → inhibited |

* It's highly recommended to use fresh strains for this isolate.

| Pack Size | Ordering References | Weight |
|-------------------------------|---------------------|----------|
| 5000 mL 250 Tests of 20 mL | VB912 | 373.5 g |
| 25L 1250 Tests of 20 mL | VB913-25 | 1867.5 g |

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Inappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

- Catalogue reference
- Consult instructions for use
- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity
- Protect from light
- Manufacturer

Need some Technical Documents?

Available for download on www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach
ATCC® is a registered trademark of the American Type Culture Collection
NT-EXT-017 V8.1 / 07-May-24

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