CHROMagar™ Streptococcus

Instructions For Use

Available in several languages

NT-EXT-137

Version 1.0

Click below for:


ENGLISH

FRANCAIS

ESPAÑOL

DEUTSCH

CHROMagar™ Streptococcus plate

MEDIUM PURPOSE

Chromogenic medium for the isolation and the differentiation of *Streptococcus* spp.

Environmental *Streptococci* are considered to be major mastitis pathogens because they typically cause high somatic cell counts, persistent infections and milk bacterial contamination. CHROMagar™ Streptococcus was designed to help in the Mastitis diagnosis, particularly as an aid to differentiate species involved in the udder infection.

There are other applications where CHROMagar™ Streptococcus has also shown to be a useful tool to differentiate *Streptococci* in samples heavily loaded with a complex bacterial flora.

COMPOSITION

The product is composed of a powder base (B) and a supplement (S).

Product	=	Base (B)	Supplement (S)
Total g/L		44.9 g/L	8.0 mL/L
Composition g/L		Agar 15.0 Peptones and growth factors 20.0 Salts 7.5 Chromogenic and selective mix 2.4	Growth factors
Aspect		Powder Form	Liquid Form
STORAGE		15-30 °C	15-30 °C
FINAL MEDIA pH		7.3 +/- 0.2	

PREPARATION (Calculation for 1 L)

Step 1

Preparation

- Disperse slowly 44.9 g of powder base in 1 L of purified water.
- Add 8.0 mL of the CHROMagar™ Streptococcus supplement.
- Stir until agar is well thickened.
- Autoclave at 121 °C during 15 min.

Step 2

Pouring

- Cool in a water bath to 45-50 °C.
- Swirl or stir gently to homogenize.
- Pour medium into sterile Petri dishes.
- Let it solidify and dry.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to one month at 2-8 °C if properly prepared and protected from light and dehydration.

Warning: Sometimes the background can be heterogeneous but do not interfere with the performance of the media.

INOCULATION

Related samples can be processed by direct streaking on the plate.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate.
- Incubate at 37 °C in aerobic conditions and read at 24 h.

Typical Samples

e.g. Milk


Direct streaking
or spreading technique

CHROMagar™ Streptococcus

INTERPRETATION

Microorganism	Typical colony appearance
<i>S. agalactiae</i>	→ blue
<i>S. uberis</i>	→ blue to metallic blue
<i>Enterococcus</i>	→ mauve to violet-mauve
Other <i>Streptococcus</i>	→ various
Other Gram (+) bacteria	→ inhibited
Gram (-) bacteria	→ inhibited

Typical colony appearance

The pictures shown are not contractual.

PERFORMANCE & LIMITATIONS

- Rare *Staphylococcus* can grow on the medium but they can be differentiated with a catalase test.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms.

Good preparation of the medium can be tested, isolating the ATCC strains below:

Microorganism	Typical colony appearance
<i>S. agalactiae</i> ATCC® 13813	→ blue
<i>S. uberis</i> ATCC® 700407	→ steel blue
<i>E. faecalis</i> ATCC® 29212	→ mauve
<i>C. albicans</i> ATCC® 10231	→ inhibited
<i>E. coli</i> ATCC® 25922	→ inhibited
<i>S. aureus</i> ATCC® 25923	→ inhibited

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Inappropriate storage may affect the shelf life of the product.
- Recap the bottles tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

-  Catalogue reference
-  Consult instructions for use
-  Quantity sufficient for X liters of media
-  Expiry date
-  Required storage temperature
-  Store away from humidity
-  Protect from light
-  Manufacturer

Need some Technical Documents?

Available for download on www.CHROMagar.com

• Certificate of Analysis (CoA) --> One per Lot

• Material Safety Data Sheet (MSDS)

Pack Size

500 g

≈550 Tests
of 20 mL

=

Ordering Reference

CQ39-500

Base (B)

CQ39B-500G
Weight: 500 g
▽ 11.1 L

+

Supplement (S)

CQ39S-05
Volume: 100 mL
▽ 12.5 L

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach. ATCC® is a registered trademark of the American Type Culture Collection.
NT-EXT-137 V1.0 / 22-January-2026

CHROMagar
The Chromogenic Media Pioneer

 CHROMagar 29 Avenue George Sand,
93210 La Plaine Saint-Denis - France
Email: CHROMagar@CHROMagar.com
Tel +33 (0)1.45.48.05.05. Website: www.CHROMagar.com

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