
# Leaflet

LF-EXT-077

Version 2.0

## CHROMagar™ B.cepacia

Click below:


EN

FR

ES

CHROMagar™  
The Chromogenic Media Pioneer

# ● CHROMagar™ B.cepacia

**For detection and numeration  
of *Burkholderia cepacia* complex**

## For detection and enumeration of *Burkholderia cepacia* complex

### Intended Use

*Burkholderia cepacia* complex (BCC), is a group of lactose-nonfermenting, Gram-negative bacteria composed of at least 20 different species, including *B. cepacia*, *B. multivorans*, *B. cenocepacia* and *B. contaminans*. BCC species can survive for long periods and multiply in hostile environments such as antiseptic and disinfectant solutions, distilled water, whirlpool baths, nebulizers and urinary catheters. Moreover, bacteria from the BCC show multidrug resistance.

BCC is among the most important pathogens isolated from cystic fibrosis patients and in hospital acquired infections. The pathogen causes life-threatening infections in cystic fibrosis patients which is associated with pulmonary function decline. It is also a cause of fatal disease in patient with chronic granulomatous disease. Several BCC outbreaks have been described due to contaminated medications, medical products and equipments. Therefore, a rapid detection and identification of bacteria from BCC are critical point for cystic fibrosis patients, but also for various non-sterile products used in hospitals and in environment.

The slowest growing species of BCC can be missed on conventional media such as Drigalski or MacConkey Agar due to the overgrowth of other organisms. CHROMagar™ B.cepacia is a strongly selective chromogenic medium which will detect most of the bacteria from the BCC within 24 - 36h.

### Medium Performance

#### 1 HIGH INTENSITY

*Burkholderia cepacia* complex (BCC) colonies develop with an intense green-blue colony colour, clearly visible to the naked eye.

#### 2 LIMITS CONTAMINANTS

Mold is largely inhibited.

#### 3 ALL IN ONE

Save time, the pre-weighed medium powder can be added to distilled water. No supplement is required.

#### 4 STRAINS OF INTEREST

A few non-fermenting strains can form green-blue colonies. Some of these strains detected were reported as of diagnostics interest.

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### Medium Description

|                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Powder Base</b>                                                                                                                                                                                                                                    | Total ..... 32.9 g/L<br>Agar ..... 15.0<br>Peptones ..... 11.0<br>Growth factors ..... 6.6<br>Chromogenic and selective mix ..... 0.3<br>Storage at 15/30 °C - pH: 6.2 +/- 0.2<br><b>Shelf Life ..... &gt; 18 months</b> |
| Usual Samples                                                                                                                                                                                                                                         | Clinical: broncho-alveolar lavages, sputum, nasopharyngeal aspirations, and oropharyngeal swabs.<br>Pharma: Non-sterile products and purified water.                                                                     |
| Procedure                                                                                                                                                                                                                                             | Direct streaking or inoculated membrane. Incubate at 35-37 °C during 36-72h in aerobic conditions for clinical samples.                                                                                                  |
| Scientific Publications on this product: available on <a href="http://www.CHROMagar.com">www.CHROMagar.com</a><br>Please read carefully the instructions for use (IFU document) available on <a href="http://www.CHROMagar.com">www.CHROMagar.com</a> |                                                                                                                                                                                                                          |

### Plate Reading

- *Burkholderia*  
→ green-blue +/- blue halo

### Order References


Please use these product references when contacting your local distributor:

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