CHROMagar™ Mycoplasma

Instructions For Use

Available in several languages

NT-EXT-123

Version 1.0

Click below for:


ENGLISH

FRANÇAIS

ESPAÑOL

DEUTSCH

CHROMagar™ Mycoplasma plate

CHROMagar™ Mycoplasma

Chromogenic medium for the detection of *Mycoplasma bovis*

INTENDED USE

Bovine infections can result in a variety of clinical signs, such as pneumonia mainly but also arthritis, mastitis and keratoconjunctivitis, none of which are *M. bovis*-specific. Laboratory diagnosis is therefore important.

Identification of *M. bovis* from nose swabs, broncho-alveolar lavage samples can be done by streaking the sample directly onto CHROMagar™ Mycoplasma agar plates. The detection of *M. bovis* by red coloration allows an easy detection with an fried egg aspect under binocular loupe after 3-7 days of incubation at 37 °C under CO₂ atmosphere.

COMPOSITION

The product is composed of a powder base (B) and 3 supplements (S1) + (S2) + (S3).

| Product | = | Base (B) | + | Supplement (S1) | + | Supplement (S2) | + | Supplement (S3) |
|-----------------|---|---|---|-----------------|---|-------------------------------|---|-----------------|
| Total g/L | | 43.8 g/L | | 1 g/L | | 1 g/L | | 150 mL/L |
| Composition g/L | | Agar 10.0 Peptones and yeast extracts 25.0 Mineral salts 6.3 Growth factors 2.2 Chromogenic and selective mix 0.3 | | Mineral salts | | Chromogenic and selective mix | | Growth factors |
| Aspect | | Powder Form | | Powder Form | | Powder Form | | Lyophilized |
| STORAGE | | 15-30 °C | | 15-30 °C | | 2-8 °C | | 2-8 °C |
| FINAL MEDIA pH | | 7.6 +/- 0.2 | | | | | | |

PREPARATION (Calculation for 1 L)

Step 1

Preparation of the base (B) CHROMagar™ Mycoplasma

- Suspend the powder base CHROMagar™ Mycoplasma in the proportion of 43.8 g into 850 mL of purified water.
- Stir until agar is well thickened.
- Heat and bring to boil (100 °C) while swirling or stirring regularly.
- Autoclave at 121 °C during 15 min.
- Cool in a water bath to 45-50°C, swirling or stirring gently.

Step 2

Preparation of supplements (S1), (S2) and (S3)

- Resuspend 1 g of CHROMagar™ Mycoplasma supplement S1 in 10 mL purified water. Swirl well until complete dissolution and sterilize by filtration at 0.45 µm.
- **WARNING: The supplement S1 has a pH 12 in solution. Take the appropriate precautions.**
- Resuspend 1 g of CHROMagar™ Mycoplasma supplement S2 in 5 mL purified water. Swirl well until complete dissolution and sterilize by filtration at 0.45 µm. The solution is yellow with foam.
- **Advice : To avoid foam you can resuspend 1 g of CHROMagar™ Mycoplasma supplement S2 in a mix of 5 mL of purified water and 50 µl of NaOH at a concentration of 100 g/L. Swirl well until complete dissolution and sterilize by filtration at 0.45 µm. The solution is yellow.**
- Resuspend sterily 10 vials of CHROMagar™ Mycoplasma supplement S3 with 5 mL purified water each. Swirl well and add another 10 mL purified water sterily to each vial.
- Swirl well until completely dissolved.

Step 3

Base (B) + supplements (S1), (S2) and (S3)

- Add 10 mL of supplement S1, 5 mL of supplement S2 and 150 mL of supplement S3 (10 vials) into slurry. Mix well.

Step 4

Pouring

- Pour into sterile Petri dishes with a diameter of 55 mm (usually used for the research of *Mycoplasma*), but other plate formats can also be used.
- Let it solidify and dry.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to 4 weeks under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.
- If not fully used, rehydrated CHROMagar™ Mycoplasma supplements S1 and S2 can be stored up to 1 month at -20°C.

CHROMagar™ Mycoplasma

INOCULATION

Related samples are inoculated by direct streaking on the plate.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate.
- Incubate at 37 °C for 3-7 days in a CO₂ atmosphere.

Typical samples


e.g. nasal swab, broncho-alveolar lavage samples

INTERPRETATION

Qualitative reading and interpretation of the Petri dishes.

| Microorganism | Typical colony appearance |
|-------------------------|-----------------------------------|
| <i>Mycoplasma bovis</i> | → red with a fried egg morphology |
| Gram (+) bacteria | → mostly inhibited |
| Gram (-) bacteria | → inhibited |
| Yeast and mould | → inhibited |

Typical colony appearance

The pictures shown are not contractual.

LIMITATIONS AND COMPLEMENTARY TESTS

- Some particular resistant strains from annexe flora may grow as metallic blue colonies.
- Some strains of *Mycoplasma* spp. may appear red with a fried egg morphology.
- The fried egg morphology is visible under a microscope.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the following ATCC strains:

| Microorganism | Typical colony appearance |
|---|-----------------------------------|
| <i>Mycoplasma bovis</i> ATCC® 25523 | → red with a fried egg morphology |
| <i>Mycoplasma bovis</i> ATCC® 25025 | → red with a fried egg morphology |
| <i>Staphylococcus aureus</i> ATCC® 25923 | → inhibited |
| <i>Klebsiella</i> ATCC® 13883 | → inhibited |

REFERENCES

Σ Pack Size

1000 mL

=

125 Tests
of 8 mL

Ordering Reference

MB640

=

Base (B) Supplement (S1) Supplement (S2) Supplement (S3)

MB640(B) + MB640(S1) + MB640(S2) + MB640(S3)
 Weight: 43.8 g Weight: 1 g Weight: 1 g 10 vials (1 vial for 100 mL of final media)

WARNINGS AND PRECAUTIONS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

-  Catalogue reference
-  Consult instructions for use
-  Quantity of powder sufficient for X liters of media
-  Expiry date
-  Required storage temperature
-  Store away from humidity
-  Protect from light
-  Manufacturer

Need some Technical Documents?

Available for download on www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach
 ATCC® is a registered trademark of the American Type Culture Collection

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