CHROMagar™ **P.aeruginosa**

Instructions For Use

Available in several languages

NT-EXT-163

Version 1.0


Click below for:

ENGLISH

FRANCAIS

ESPAÑOL

DEUTSCH

CHROMagar™ P.aeruginosa

Instructions For Use
NT-EXT-163 V1.0 / 18-Jun-26

ENGLISH

Instructions For Use

Chromogenic medium for the detection and enumeration of *Pseudomonas aeruginosa*.

REFERENCES

Pack size

100 g

≈440 tests
of 7 mL

=

Order references

PA27-100G

▽ 3.1 L

INTENDED USE

CHROMagar™ P.aeruginosa is a selective chromogenic culture medium intended for the direct enumeration, differentiation, and presumptive identification of *Pseudomonas aeruginosa*. The test is performed using water samples filtered through membranes according to the method specified in ISO 16266. Results can be interpreted after 18h to 24 h of incubation at 36 ± 2 °C under aerobic conditions.

COMPOSITION

Product	=	Base
Total g/L		32.0 g/L
Composition g/L		Agar 12.0 Peptones 8.0 Salts 9.0 Growth factors 2.2 Chromogenic and selective mix 0.8
Aspect		Powder Form*
STORAGE		2/8 °C
MEDIA pH		7.4 ± 0.2

Need some
Technical Documents?

Available
for download on
www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

* The color may vary from yellow to pink, but this has no effect on the medium's performance.

PREPARATION (Calculation for 1 L)

Step 1

Preparation

- Let the powder return to room temperature.
- Gently dissolve 32 g of in 1 L of purified water.
- Stir until agar is well thickened.
- Heat to 100 °C until the agar is completely dissolved.

DO NOT HEAT TO 121 °C.

Advice for microwaving: After the initial boil, remove from the microwave, stir gently, then return to the microwave for repeated short heating intervals until the agar granules have completely melted (large bubbles replacing the foam).

Advice for autoclave heating: Heat for 10 minutes at 100°C.

Step 2

Pouring

- Cool in a water bath at 45–50 °C, stirring gently.
- Pour into sterile Petri dishes.
- Allow to solidify and dry.

Storage

- Store away from light.
- Prepared plates can be kept at room temperature for one day.
- Prepared plates can be stored for up to 2 months in the refrigerator (2–8 °C) if they are properly prepared and protected from light and drying out.

Note: The agar base may sometimes appear heterogenous, but this does not affect the medium's performance.

SPECIMEN COLLECTION AND HANDLING

CHROMagar™ P.aeruginosa can be used with the following samples: drinking water and recreational water, treated and untreated.

INOCULATION

Water samples can be used for the direct enumeration of *Pseudomonas aeruginosa* after membrane filtration.

Warning: Do not use polycarbonate nor regenerated cellulose membranes.

Advice: For the filtration method, we recommend using cellulose nitrate membranes for optimal performance.

If the prepared plate has been refrigerated, allow it to return to room temperature before use.

- Place the membrane directly on the agar.
- Incubate at 36 ± 2 °C during 18h to 24h under aerobic conditions.


INTERPRETATION

Microorganism	Typical colony appearance
<i>P. aeruginosa</i>	→ red *
<i>Pseudomonas</i> spp.	→ inhibited or colourless
Other Gram (-) bacteria	→ mainly inhibited
Gram (+) bacteria	→ inhibited
Yeast and fungi	→ inhibited


*Some *P. aeruginosa* strains may appear dark red to brown due to natural pigmentation.

Typical colony appearance

The pictures shown are not contractual.

P. aeruginosa

P. aeruginosa

Pseudomonas spp.

LIMITATIONS AND COMPLEMENTARY TESTS

Some rare strains of *Aeromonas caviae/hydrophila* may appear red on CHROMagar™ P.aeruginosa. The natural pigmentation of *Pseudomonas aeruginosa* may stain the filter.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms.

The good preparation of the medium can be tested by streaking the following ATCC strains:

Microorganism	Typical colony appearance
<i>P. aeruginosa</i> WDCM 00025 / ATCC® 27853	→ red
<i>C. tropicalis</i> ATCC® 1369	→ inhibited
<i>E. faecalis</i> WDCM 00087 / ATCC® 29212	→ inhibited
<i>E. coli</i> WDCM 00013 / ATCC® 25922	→ inhibited

WARNINGS AND PRECAUTIONS

- This laboratory product should be used only by trained personnel (healthcare professional, etc). Handle appropriately with procedures and good laboratory practices.
- Use of the medium may be difficult for people who have problems recognising colours.
- For optimal microbial detection, water samples must be collected and transported properly and in a manner appropriate for each sample, in accordance with good laboratory practices.
- Culture media should not be used as manufacturing material or components.
- Do not ingest or inhale the product.
- Do not use the product after the expiry date.
- Do not use the product if it show any evidence of contamination or any sign of deterioration.
- Do not use the product if the packaging is damaged.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Be sure to tightly reseal the bottles after each use and store them in a dry place, away from moisture and light.
- Reading and interpretation should be performed using isolated colonies.
- Interpretation of the test results should be made taking into consideration colonial and microscopic morphology and if necessary, the results of any other tests performed.
- Laboratory, chemical or biohazardous wastes must be handled and discarded in accordance with all local and national regulations.
- For hazard and precaution recommendations related to some chemical components in this medium, please refer to the pictogram(s) mentioned on the labels. The Safety Data Sheet (SDS) is available on www.chromagar.com

CHROMagar™ *P.aeruginosa*

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 min.

LITERATURE REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

-  Catalogue reference
-  Consult instructions for use
-  Sufficient quantity for X liters of media
-  Expiry date
-  Required storage temperature
-  Store away from humidity
-  Protect from light
-  Manufacturer

CHROMagar™ is a trademark created by Dr A. Rambach
ATCC® is a registered trademark of the American Type Culture Collection

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