CHROMagar™ **Listeria**

Instructions For Use

Available in several languages

NT-EXT-167

Version 1.0


Click below for:

ENGLISH

FRANCAIS

ESPAÑOL

DEUTSCH

MEDIUM PURPOSE

Chromogenic medium for detection, isolation and enumeration of *L. monocytogenes* and *Listeria* spp.

This medium meets the ISO 11290-1: 2017 and ISO 11290-2: 2017 norms.

Listeria monocytogenes is a widespread bacteria that can survive and multiply at low temperatures (refrigerators), and high salt concentrations, forming biofilms that are difficult to remove. These biofilms can constitute permanent reservoirs of this pathogen for processed foods. The presence of *Listeria* spp. in the industrial environment is an indicator of *L. monocytogenes* contamination. Detection, prevention and control of *Listeria* spp. and *L. monocytogenes* are essential in the food industry.

COMPOSITION

The product is composed of a powder base (B) and a supplement (S).

Product	=	Base (B)	+	Supplement (S)
Total g/L		51.5 g/L		9.0 g/L
Composition g/L		Agar 15.0 Peptone and yeast extract 23.0 NaCl 5.0 Chromogenic mix 8.5		Selective and enrichment mix 9.0
Aspect		Powder Form		Powder Form
STORAGE		2/30 °C		2/8 °C
FINAL MEDIA pH		7.0 +/- 0.5		

Typical Samples

e.g. meat samples, dairy products, vegetables, sea-food and environmental.

PREPARATION (Calculation for 1 L)

Step 1

Preparation of the base
CHROMagar™
Listeria base (B)

- Disperse slowly 51.5 g of powder base in 1 L of purified water.
- Stir until agar is well thickened.
- Heat to 121 °C +/- 1 °C for 15 minutes.
- Cool in a water bath to 45-50 °C.

Step 2

Preparation of the
CHROMagar™
Listeria supplement (S)


- Add 9 g to 40 mL of STERILE purified water.
- Add a magnetic bar and stir at high speed (700-1000 rpm) at least 30 minutes WITHOUT HEATING until getting a creamy homogeneous suspension as below:

Final medium HELPING CALCULATION

- 1 L --> Add 9 g of supplement in 40 mL of sterile water
- 5 L --> Add 45 g of supplement in 200 mL of sterile water

STAGE 1:

Possible presence of grains on the surface of the vial and in the mix

Aspect at this step:

- Liquid with undissolved grains

-->

STAGE 2: Formation of white dense foam on the surface and rise of foam

Aspect at this step:

- Dense/compact foam (SURFACE)
- Foamy mix

-->

STAGE 3: less dense/compact Foam, more creamy mix

Aspect at this step:

- No Rise of foam
- More ventilated foam
- Liquid

-->

STAGE 4:
VERY GOOD MIX

Aspect at this step:

- Homogeneous
- Fluid (liquid)
- Without grains nor foam

GOOD PREPARATION

Step 3

Mixing of the prepared base (B) and the prepared supplement (S)

- Place the melted 45-50 °C cooled CHROMagar™ *Listeria* base under gentle stirring.
- Add the homogeneous reconstituted supplement, keeping the gentle stirring during 1 or 2 minutes until complete homogenisation.
- Pour IMMEDIATELY into sterile Petri dishes.

Warning: DO NOT STACK THE PETRI DISHES

- Let them cool down. Let dry.

Storage

- Store in the dark before use.
- Prepared medium plates can be kept for one day at room temperature.
- Plates can be stored for up to one month under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.

INOCULATION

Please refer to Part **Procedures**

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Collection and transport of the samples should be well handled and adapted to the particular specimen according to good laboratory practices.

PROCEDURES

CHROMagar™ *Listeria* Method


• For the Detection of Presence / Absence of *Listeria* spp. and *L. monocytogenes* in all human food products and environmental samples (Refer to Fig. 1)

- 1 Perform an enrichment step of the sample, mixing x g or x mL in 9 mL of half Fraser broth.
→ 24 ± 2 h at 30 ± 1 °C
- 2 Spread 0.1 mL on a well dried CHROMagar™ *Listeria* plate, by streaking on 1/3 of the plate until liquid is absorbed. Isolate on the rest of the plate and incubate.
→ 24 ± 2 h at 37 ± 1 °C
- 3 Typical colonies:
NO → Absence of *Listeria* spp. and *L. monocytogenes*
YES → All samples identified as positive by the alternative method must be confirmed in one of the following ways:
1/ Using conventional tests described in the standardized methods by CEN or ISO (including the purification step).
- 4 2/ Using CHROMagar™ Identification *Listeria* directly from the suspect colony of *L. monocytogenes* without purification step.

In the event of discordant results (positives presumptive by the alternative method, non-confirmed by one of the means described above) the laboratory must follow the necessary steps to ensure validity of the results observed. Test portions weighing more than 25 g have not been tested.

For environmental samples, typical colonies should be absent after 24 h, incubate for additional 24 ± 2 h.

Figure 1: Detection Method

* Typical colonies must be tested for confirmation


CHROMagar™ *Listeria* Enumeration Method

• For the Enumeration of *Listeria* spp. and *L. monocytogenes* in all human food products and environmental samples (Refer to Fig. 2)

- 1 Perform a revivification step to the sample and mix x g or x mL in 9 mL of BPW or non-selective half Fraser broth.
→ 1 h ± 5 min at 20 ± 1 °C
- 2 Spread 0.1 mL on a well dried CHROMagar™ *Listeria* plate (in case of estimation of small numbers spread 1 mL on 3 plates) and incubate.
→ 48 ± 2 h at 37 ± 1 °C
A first reading at 24 h allow a more rapid detection of heavily contaminated samples. However, the final result count is reached after 48 h ± 2 h. Refer to ISO 7218 standard for calculation and interpretation of results. The spreading of a single dilution may present a risk of gap counts compared to the reference method.
- 3 Typical colonies:
NO → Absence of *Listeria* spp. and *L. monocytogenes*
YES → Enumerate colonies and confirm suspect colonies directly from CHROMagar™ *Listeria*. If you do not confirm 5 counted colonies, it can create a risk of making an overestimated result because of the possible presence of typical colonies that would not be *L. monocytogenes*. If confirmation has already been done during the detection phase with CHROMagar™ *Listeria* method, a new confirmation is not necessary after the CHROMagar™ *Listeria* enumeration method. Otherwise, all samples identified as positive by the alternative method must be confirmed in one of the following ways:
1/ Using conventional tests described in the standardized methods by CEN or ISO (including the purification step).
- 4 2/ Using CHROMagar™ Identification *Listeria* directly from the suspect colony without purification step.

In the event of discordant results (positives presumptive by the alternative method, non-confirmed by one of the means described above) the laboratory must follow the necessary steps to ensure validity of the results observed. Test portions weighing more than 25 g have not been tested. For environmental samples, typical colonies should be absent after 24 h, incubate for additional 24 ± 2 h.

Figure 2: Enumeration Method


* Typical colonies must be tested for confirmation

INTERPRETATION

Microorganism	Typical colony appearance
<i>L. monocytogenes</i>	→ blue, diameter less than 3 mm, regular and white halo
<i>Listeria</i> spp.	→ blue, diameter less than 3 mm
Other microorganisms	→ blue, colourless, other colour, inhibited

Typical colony appearance on CHROMagar™ *Listeria*

Non-contractual photo.

PERFORMANCE & LIMITATIONS

- Some strains of *L. ivanovii* may produce blue colonies with a white halo.
- Some strains of *Bacillus cereus* may also grow blue colonies, recognizable by their larger size and irregular shape surrounded by a very large white halo.
- Some recently described *Listeria* species, e.g. *L. fleischmannii* sp. nov., *L. marthii* sp. nov., *L. weihenstephanensis* sp. nov., and *L. rocourtiae* sp. nov. may not be detected or counted.
- Some organisms other than *Listeria* spp. can produce blue colonies on this medium. See Annex C of NF EN ISO 11290-1

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the ATCC strains below:

Microorganism	Typical colony appearance
<i>L. monocytogenes</i> ATCC® 13932 / WDCM 00021	→ blue with halo
<i>L. innocua</i> ATCC® 33091	→ blue without halo
<i>E. faecalis</i> ATCC® 29212 / WDCM 00087	→ inhibited
<i>E. coli</i> ATCC® 25922 / WDCM 00013	→ inhibited

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

- Catalogue reference
- Consult instructions for use
- Quantity sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity
- Protect from light
- Manufacturer

Need some Technical Documents?

Available for download on www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

Pack Size

Ordering References

Base (B)

Supplement (S)

500 g

≈485 Tests
of 20 mL

=

LM85-500

=

LM85B-500G

Weight: 500 g

▽ 9.7 L

+

2 x LM85S-05

Weight: 2 x 45 g

▽ 2 x 5 L

NT-EXT-167 V1.0 / 01-Sep-2026

CHROMagar is a trademark created by Dr A. Rambach. ATCC® is a registered trademark of the American Type Culture Collection. WDCM : World Data Center for Microorganism.

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