


|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA base</b><br><b>Ref : AO883/B</b> | AD-ASQ-265  |
|   |  | Version 1.0 |
| Page 1/7  |  | 04-Oct-2016 |

## SECTION 1: Identification of the mixture and of the company / undertaking

### 1.1 Product identifier

· Trade name: CHROMagar™ AOLA base

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

· Application of the substance / the preparation : Reagent for Laboratory Use

### 1.3 Details of the supplier of the safety data sheet

· Manufacturer/Supplier :

CHROMagar  
 4, place du 18 Juin 1940 - 75006 Paris  
 Tel: +33 1 45 48 05 05  
 Email : chromagar@chromagar.com

### 1.4 Emergency telephone number :

France (ORFILA 24h/24) - Tel : +33 (0)1 45 42 59 59  
 Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
 Belgique - Tel : 32 070/245 245  
 UK : +44 (0)-870-190-6777 (CARECHEM24)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

· Classification according to Regulation (EC) No 1272/2008 :



GHS07 : Eye Irrit. 2 H319 Causes serious eye irritation.

### 2.2 Label elements

· Labelling according to Regulation (EC) No 1272/2008 :

The product is classified and labelled according to the CLP regulation

· Hazard pictograms :



GHS07

Warning

· Signal word :

· Hazard statements :

H319 : Causes serious eye irritation.

· Precautionary statements :

P280 : Wear protective gloves/protective clothing/eye protection/face protection.

P264 : Wash hands thoroughly after handling.

P305 + P351 + P338 : If in eyes rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 : If eye irritation persists: Get medical advice/attention.

· Labelling of packages where the contents do not exceed 125 ml

· Hazard pictograms :




GHS07

Warning

· Signal word :

· Hazard statements :

Not classified as hazardous.

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA base</b><br><b>Ref : AO883/B</b> | AD-ASQ-265  |
|   |  | Version 1.0 |
| Page 2/7  |  | 04-Oct-2016 |

### 2.3 Other hazards


- **Results of PBT and vPvB assessment**

- **PBT :** Not applicable.
- **vPvB :** Not applicable.

## SECTION 3: Composition / information on ingredients

### 3.2 Mixtures

- **Description :** Mixture: consisting of the following components.
- **Dangerous components :**

|                   |   |        |
|-------------------|---|--------|
| CAS: 7447-41-8    | Lithium chloride  | 10-25% |
| EINECS: 231-212-3 |  Acute Tox. 4, H302; Acute Tox. 4, H312; Eye Irrit. 2, H319; STOT SE 3, H335 |        |

RTECS: OJ 5950000  
Reg.nr.: 01-2119560574-35

- **SVHC :** No
- **Additional information :** For the wording of the listed risk phrases, refer to section 16.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

- **General information :** Position and transport stably in side position.  
In case of irregular breathing or respiratory arrest provide artificial respiration.
- **After excessive inhalation :** Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact :** Generally the product does not irritate the skin.  
In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.  
If skin irritation continues, consult a doctor.
- **After eye contact :** Rinse opened eye for several minutes under running water.
- **After swallowing :** Do not induce vomiting; call for medical help immediately.

### 4.2 Most important symptoms and effects, both acute and delayed :

No further relevant information available.

### 4.3 Indication of any immediate medical attention and special treatment needed :

No further relevant information available.

## SECTION 5: Firefighting measures


### 5.1 Suitable extinguishing agents :

CO<sub>2</sub>, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

**5.2 Special hazards arising from the substance or mixture :** No further relevant information available.

### 5.3 Advice for firefighters

- **Protective equipment :** No special measures required.

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA base</b><br><b>Ref : AO883/B</b> | AD-ASQ-265  |
|   |  | Version 1.0 |
| Page 3/7  |  | 04-Oct-2016 |

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures :

No special measures required. Avoid formation of dust.

### 6.2 Environmental precautions :

Do not allow product to reach sewage system or any water course. Dilute with plenty of water.

### 6.3 Methods and material for containment and cleaning up : Pick up mechanically.

### 6.4 Reference to other sections :

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling :

Keep away from heat and direct sunlight.  
Thorough dedusting.

### · Information about fire - and explosion protection:

No special measures required.

### 7.2 Conditions for safe storage, including any incompatibilities :

#### · Storage :

· **Requirements to be met by storerooms and receptacles :** Store away from light and moisture.  
No special requirements.

· **Information about storage in one common storage facility :** Not required.

#### · Further information about storage conditions :

Store in cool,  
dry conditions in well sealed receptacles.

#### · Recommended storage temperature :

15-30°C

### 7.3 Specific end use(s) :

No further relevant information available.

## SECTION 8: Exposure controls / personal protection

· **Additional information about design of technical facilities :** No further data; see section 7.

### 8.1 Control parameters

#### · Ingredients with limit values that require monitoring at the workplace :

##### 7447-41-8 lithium chloride (10-25%)

MAK (Switzerland) Short-term value: 0.2 e mg/m<sup>3</sup>

Long-term value: 0.2 e mg/m<sup>3</sup>

SSc;als Li

#### · Additional information :

The most current valid lists have been used as a basis for the production of this document.

### 8.2 Exposure controls

#### · Personal protective equipment :


#### · General protective and hygienic measures :

The usual precautionary measures are to be adhered to when handling chemicals.

#### Respiratory protection :

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.



|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA base</b><br><b>Ref : AO883/B</b> | AD-ASQ-265  |
|   |  | Version 1.0 |
| Page 4/7  |  | 04-Oct-2016 |

**Protection of hands :**

Protective gloves.  
EN374 gloves for chemicals.



**Gloves material :**

Butyl rubber, BR

**Eye protection :**

Safety glasses.



**Body protection :**

Protective work clothing.



## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

· General Information

· Appearance :

**Form :** Powder

**Colour :** Beige

· Odour : Characteristic

· Odour threshold : Not determined.

· pH-value : Not applicable.

· Change in condition

**Melting point/Melting range :** Undetermined.

**Boiling point/Boiling range :** Undetermined.

· Flash point : Not applicable.

· Flammability (solid, gaseous) : Not determined.

· Ignition temperature : Not applicable.

· Decomposition temperature : Not determined.

· Self-igniting : Product is not self igniting.

· Danger of explosion : Product does not present an explosion hazard.

· Explosion limits :

**Lower :** Not determined.

**Upper :** Not determined.

· Vapour pressure : Not applicable.

· Density at 20°C : 0.5 +/- 0.1 g/cm<sup>3</sup>

· Relative density : Not determined.

· Vapour density : Not applicable.

· Evaporation rate : Not applicable.

· Solubility in / Miscibility with water : Partly soluble.

· Segregation coefficient (n-octanol/water) : Not determined.


· Viscosity :

**Dynamic :** Not applicable.

**Kinematic :** Not applicable.

**9.2 Other information :**

No further relevant information available.

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA base</b><br><b>Ref : AO883/B</b> | AD-ASQ-265  |
|   | Page 5/7   | Version 1.0 |
|   |  | 04-Oct-2016 |

## SECTION 10: Stability and reactivity

- 10.1 Reactivity :** No further relevant information available.
- 10.2 Chemical stability :**
- **Thermal decomposition / conditions to be avoided :** No decomposition if used according to specifications.
- 10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- 10.4 Conditions to avoid :** No further relevant information available.
- 10.5 Incompatible materials :** No further relevant information available.
- 10.6 Hazardous decomposition products :** No dangerous decomposition products known.


## SECTION 11: Toxicological information

- 11.1 Information on toxicological effects**
- **Acute toxicity :** Based on available data, the classification criteria are not met.
  - **LD/LC50 values relevant for classification :**

|           |                       |
|-----------|-----------------------|
| 7447-41-8 | Lithium chloride      |
| Oral      | LD50 526 mg/kg (rat)  |
| Dermal    | LD50 1100 mg/kg (ATE) |
  - **Primary irritant effect :**
  - **On the skin :** Based on available data, the classification criteria are not met.
  - **On the eyes :** Causes serious eye irritation.
  - **Respiratory or skin sensitisation :** Based on available data, the classification criteria are not met.
  - **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction) :**
  - **Germ cell mutagenicity :** Based on available data, the classification criteria are not met.
  - **Reproductive toxicity :** Based on available data, the classification criteria are not met.
  - **STOT-single exposure :** Based on available data, the classification criteria are not met.
  - **STOT-repeated exposure :** Based on available data, the classification criteria are not met.
  - **Aspiration hazard :** Based on available data, the classification criteria are not met.

## SECTION 12: Ecological information

- 12.1 Toxicity**
- **Aquatic toxicity :** No further relevant information available.
- 12.2 Persistence and degradability :** No further relevant information available.
- 12.3 Bioaccumulative potential :** No further relevant information available.
- 12.4 Mobility in soil :** No further relevant information available.
- **Additional ecological information :**
  - **General notes :** Generally not hazardous for water.
- 12.5 Results of PBT and vPvB assessment**
- **PBT :** Not applicable.
  - **vPvB :** Not applicable.
- 12.6 Other adverse effects :** No further relevant information available.

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA base</b><br><b>Ref : AO883/B</b> | AD-ASQ-265  |
|   |  | Version 1.0 |
| Page 6/7  |  | 04-Oct-2016 |

### SECTION 13: Disposal considerations

#### 13.1 Waste treatment methods

- **Recommendation :** Must be specially treated according to official regulations.
- **Uncleaned packaging :**
- **Recommendation :** Disposal must be made according to official regulations.

### SECTION 14: Transport information

#### 14.1 UN-Number

- **ADR, ADN, IMDG, IATA :** Void

#### 14.2 UN proper shipping name

- **ADR, ADN, IMDG, IATA :** Void

#### 14.3 Transport hazard class(es)

- **ADR, ADN, IMDG, IATA**
- **Class :** Void

#### 14.4 Packing group

- **ADR, IMDG, IATA :** Void

**14.5 Environmental hazards :** Not applicable.

**14.6 Special precautions for user :** Not applicable.

#### 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code :

- **UN "Model Regulation":** Not applicable.
- **UN "Model Regulation":** Void

### SECTION 15: Regulatory information

#### 15.1 Safety, health and environmental regulations/legislation specific for the mixture :

- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I :** None of the ingredients is listed.

**15.2 Chemical safety assessment :** A Chemical Safety Assessment has not been carried out.

### SECTION 16: Other information


This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

#### · **Relevant phrases :**

- H302 Harmful if swallowed.
- H312 Harmful in contact with skin.
- H319 Causes serious eye irritation.
- H335 May cause respiratory irritation.

#### · **Abbreviations and acronyms :**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
 IMDG: International Maritime Code for Dangerous Goods

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA base</b><br><b>Ref : AO883/B</b> | AD-ASQ-265  |
|   |  | Version 1.0 |
| Page 7/7  |  | 04-Oct-2016 |

IATA: International Air Transport Association.

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic


SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity, Hazard Category 4

Eye Irrit. 2: Serious eye damage/eye irritation, Hazard Category 2

STOT SE 3: Specific target organ toxicity - Single exposure, Hazard Category 3

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-215  |
|   | <b>CHROMagar™ AOLA supplement E1</b><br><b>Ref : AO882/E1 – AO883/E1</b> | Version 1.0 |
| Page 1/6  |  | 27-May-2016 |

## SECTION 1: Identification of the mixture and of the company / undertaking

### 1.1 Product identifier

- Trade name: CHROMagar™ AOLA supplement E1

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

- Application of the substance / the preparation : Reagent for Laboratory Use

### 1.3 Details of the supplier of the safety data sheet

- Manufacturer/Supplier:

CHROMagar  
 4, place du 18 Juin 1940 - 75006 Paris  
 Tel: +33 1 45 48 05 05  
 Email : chromagar@chromagar.com

### 1.4 Emergency telephone number:

France (ORFILA 24h/24) - Tel : +33 (0)1 45 42 59 59

Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112

Belgique - Tel : 32 070/245 245

UK : +44 (0)-870-190-6777 (CARECHEM24)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

- Classification according to Regulation (EC) No 1272/2008 :

The product is not classified according to the CLP regulation.

### 2.2 Label elements

- Labelling according to Regulation (EC) No 1272/2008 : Not classified as hazardous
- Hazard pictograms : Not classified as hazardous
- Signal word : Not classified as hazardous
- Hazard statements : Not classified as hazardous

### 2.3 Other hazards

- Results of PBT and vPvB assessment

• PBT: Not applicable.

• vPvB: Not applicable.

## SECTION 3: Composition / information on ingredients

### 3.2 Mixtures


- Description : Mixture: consisting of the following components.
- Dangerous components : Not classified as hazardous.
- SVHC : No
- Additional information : For the wording of the listed risk phrases, refer to section 16.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

- General information : Position and transport stably in side position.  
In case of irregular breathing or respiratory arrest provide artificial respiration.
- After excessive inhalation : Provide fresh air; consult a doctor in case of medical complaints.



|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-215  |
|   | <b>CHROMagar™ AOLA supplement E1</b><br><b>Ref : AO882/E1 – AO883/E1</b> | Version 1.0 |
| Page 2/6  |  | 27-May-2016 |

- **After skin contact :** In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.  
If skin irritation continues, consult a doctor.
- **After eye contact :** Rinse opened eye for several minutes under running water.
- **After swallowing :** Do not induce vomiting; call for medical help immediately.

**4.2 Most important symptoms and effects, both acute and delayed :**

No further relevant information available.

**4.3 Indication of any immediate medical attention and special treatment needed :**

No further relevant information available.

## SECTION 5: Firefighting measures

**5.1 Suitable extinguishing agents:**

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

**5.2 Special hazards arising from the substance or mixture :** No further relevant information available.

**5.3 Advice for firefighters**

- **Protective equipment :** No special measures required.

## SECTION 6: Accidental release measures

**6.1 Personal precautions, protective equipment and emergency procedures :** No special measures required. Avoid formation of dust.

**6.2 Environmental precautions :** No special measures required.

**6.3 Methods and material for containment and cleaning up :** Pick up mechanically.

**6.4 Reference to other sections :**

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

## SECTION 7: Handling and storage

**7.1 Precautions for safe handling :**

No special measures required.

- **Information about fire - and explosion protection:**

No special measures required.

**7.2 Conditions for safe storage, including any incompatibilities :**

- **Storage :**

- **Requirements to be met by storerooms and receptacles :** Store away from light and moisture.

- **Information about storage in one common storage facility :** Not required.


- **Further information about storage conditions :**

- **Recommended storage temperature :**

2-8°C

**7.3 Specific end use(s) :**

No further relevant information available.

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-215  |
|   | <b>CHROMagar™ AOLA supplement E1</b><br><b>Ref : AO882/E1 – AO883/E1</b> | Version 1.0 |
| Page 3/6  |  | 27-May-2016 |

## SECTION 8: Exposure controls / personal protection

· **Additional information about design of technical facilities :** No further data; see section 7.

### 8.1 Control parameters

· **Ingredients with limit values that require monitoring at the workplace :**

The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.

· **Additional information :**

The most current valid lists have been used as a basis for the production of this document.

### 8.2 Exposure controls

· **Personal protective equipment :**

· **General protective and hygienic measures :**

The usual precautionary measures are to be adhered to when handling chemicals.

**Respiratory protection :** In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.



**Protection of hands :** Protective gloves.  
EN374 gloves for chemicals.



**Gloves material :** Butyl rubber, BR.  
**Eye protection :** Safety glasses.



**Body protection :** Protective work clothing.



## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

· **General Information**

· **Appearance :**

**Form :** Powder

**Colour :** Yellow

· **Odour :** Characteristic

· **Odour threshold :** Not determined.


· **pH-value :** Not determined.

· **Change in condition**

**Melting point/Melting range :** Undetermined.

**Boiling point/Boiling range :** Undetermined.

· **Flash point :** Not applicable.

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-215  |
|   | <b>CHROMagar™ AOLA supplement E1</b><br><b>Ref : AO882/E1 – AO883/E1</b> | Version 1.0 |
| Page 4/6  |  | 27-May-2016 |

|  |   |
|--|---|
| · <b>Flammability (solid, gaseous) :</b>             | Not applicable.                               |
| · <b>Ignition temperature:</b>                       | Not applicable.                               |
| · <b>Decomposition temperature:</b>                  | Not determined.                               |
| · <b>Self-igniting :</b>                             | Product is not self igniting.                 |
| · <b>Danger of explosion :</b>                       | Product does not present an explosion hazard. |
| · <b>Explosion limits :</b>                          |   |
| <b>Lower :</b>                                       | Not determined.                               |
| <b>Upper :</b>                                       | Not determined.                               |
| · <b>Vapour pressure :</b>                           | Not determined.                               |
| · <b>Density at 20°C :</b>                           | 0.5 +/- 0.1 g/cm <sup>3</sup>                 |
| · <b>Relative density :</b>                          | Not determined.                               |
| · <b>Vapour density :</b>                            | Not determined.                               |
| · <b>Evaporation rate :</b>                          | Not determined.                               |
| · <b>Solubility in / Miscibility with water :</b>    | Partly soluble.                               |
| · <b>Segregation coefficient (n-octanol/water) :</b> | Not determined.                               |
| · <b>Viscosity:</b>                                  |   |
| <b>Dynamic :</b>                                     | Not determined.                               |
| <b>Kinematic :</b>                                   | Not determined.                               |
| <b>9.2 Other information :</b>                       | No further relevant information available.    |


## SECTION 10: Stability and reactivity

|   |  |
|---|--|
| <b>10.1 Reactivity :</b>                                    | No further relevant information available.                       |
| <b>10.2 Chemical stability :</b>                            |  |
| · <b>Thermal decomposition / conditions to be avoided :</b> | No decomposition if used and stored according to specifications. |
| <b>10.3 Possibility of hazardous reactions :</b>            | No dangerous reactions known.                                    |
| <b>10.4 Conditions to avoid :</b>                           | No further relevant information available.                       |
| <b>10.5 Incompatible materials :</b>                        | No further relevant information available.                       |
| <b>10.6 Hazardous decomposition products :</b>              | No dangerous decomposition products known.                       |

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

|  |   |
|--|---|
| · <b>Acute toxicity :</b>  | Based on available data, the classification criteria are not met. |
| · <b>Primary irritant effect :</b>   |   |
| · <b>On the skin :</b>   | Based on available data, the classification criteria are not met. |
| · <b>On the eyes :</b>   | Based on available data, the classification criteria are not met. |
| · <b>Respiratory or skin sensitisation :</b>   | Based on available data, the classification criteria are not met. |
| · <b>CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction) :</b> |   |
| · <b>Germ cell mutagenicity :</b>  | Based on available data, the classification criteria are not met. |
| · <b>Reproductive toxicity :</b>   | Based on available data, the classification criteria are not met. |
| · <b>STOT-single exposure :</b>  | Based on available data, the classification criteria are not met. |
| · <b>STOT-repeated exposure :</b>  | Based on available data, the classification criteria are not met. |
| · <b>Aspiration hazard :</b>   | Based on available data, the classification criteria are not met. |

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-215  |
|   | <b>CHROMagar™ AOLA supplement E1</b><br><b>Ref : AO882/E1 – AO883/E1</b> | Version 1.0 |
| Page 5/6  |  | 27-May-2016 |

## SECTION 12: Ecological information

### 12.1 Toxicity

· **Aquatic toxicity :** No further relevant information available.

**12.2 Persistence and degradability :** No further relevant information available.

**12.3 Bioaccumulative potential :** No further relevant information available.

**12.4 Mobility in soil :** No further relevant information available.

· **Additional ecological information:**

· **General notes :** Generally not hazardous for water.

### 12.5 Results of PBT and vPvB assessment

· **PBT :** Not applicable.

· **vPvB :** Not applicable.

**12.6 Other adverse effects :** No further relevant information available.

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

· **Recommendation :** Must be specially treated according to official regulations.  
Smaller quantities can be disposed of with household waste.

· **Uncleaned packaging :**

· **Recommendation :** Disposal must be made according to official regulations.

## SECTION 14: Transport information

### 14.1 UN-Number

· **ADR, ADN, IMDG, IATA :** Void

### 14.2 UN proper shipping name

· **ADR, ADN, IMDG, IATA :** Void

### 14.3 Transport hazard class(es)

· **ADR, ADN, IMDG, IATA**

· **Class :** Void

### 14.4 Packing group

· **ADR, IMDG, IATA :** Void


**14.5 Environmental hazards :** Not applicable.

**14.6 Special precautions for user :** Not applicable.

### 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code :

Not applicable.

· **UN "Model Regulation":** Void

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-215  |
|   | <b>CHROMagar™ AOLA supplement E1</b><br><b>Ref : AO882/E1 – AO883/E1</b> | Version 1.0 |
| Page 6/6  |  | 27-May-2016 |

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the mixture :

- Directive 2012/18/EU
- Named dangerous substances - ANNEX I : None of the ingredients is listed.

**15.2 Chemical safety assessment :** A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

### · Abbreviations and acronyms :

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the

International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances


ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-216  |
|   | <b>CHROMagar™ AOLA supplement E2</b><br><b>Ref : AO882/E2 – AO883/E2</b> | Version 1.0 |
| Page 1/6  |  | 27-May-2016 |

## SECTION 1: Identification of the mixture and of the company / undertaking

### 1.1 Product identifier

- Trade name: CHROMagar™ AOLA supplement E2

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

- Application of the substance / the preparation : Reagent for Laboratory Use

### 1.3 Details of the supplier of the safety data sheet

- Manufacturer/Supplier:

CHROMagar  
 4, place du 18 Juin 1940 - 75006 Paris  
 Tel: +33 1 45 48 05 05  
 Email : chromagar@chromagar.com

### 1.4 Emergency telephone number:

France (ORFILA 24h/24) - Tel : +33 (0)1 45 42 59 59

Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112

Belgique - Tel : 32 070/245 245

UK : +44 (0)-870-190-6777 (CARECHEM24)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

- Classification according to Regulation (EC) No 1272/2008 :

The product is not classified according to the CLP regulation.

### 2.2 Label elements

- Labelling according to Regulation (EC) No 1272/2008 : Not classified as hazardous for transport
- Hazard pictograms : Not classified as hazardous for transport
- Signal word : Not classified as hazardous for transport
- Hazard statements : Not classified as hazardous for transport

### 2.3 Other hazards

- Results of PBT and vPvB assessment

- PBT: Not applicable.
- vPvB: Not applicable.

## SECTION 3: Composition / information on ingredients


### 3.2 Mixtures

- Description : Mixture: consisting of the following components.
- Dangerous components : Not classified as hazardous.
- SVHC : No
- Additional information : For the wording of the listed risk phrases, refer to section 16.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

- General information : Position and transport stably in side position.  
In case of irregular breathing or respiratory arrest provide artificial respiration.
- After excessive inhalation : Provide fresh air; consult a doctor in case of medical complaints.

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-216  |
|   | <b>CHROMagar™ AOLA supplement E2</b><br><b>Ref : AO882/E2 – AO883/E2</b> | Version 1.0 |
| Page 2/6  |  | 27-May-2016 |

- **After skin contact :** In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.  
If skin irritation continues, consult a doctor.
- **After eye contact :** Rinse opened eye for several minutes under running water.
- **After swallowing :** Do not induce vomiting; call for medical help immediately.

**4.2 Most important symptoms and effects, both acute and delayed :**

No further relevant information available.

**4.3 Indication of any immediate medical attention and special treatment needed :**

No further relevant information available.

## SECTION 5: Firefighting measures

**5.1 Suitable extinguishing agents:**

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

**5.2 Special hazards arising from the substance or mixture :** No further relevant information available.

**5.3 Advice for firefighters**

- **Protective equipment :** No special measures required.

## SECTION 6: Accidental release measures

**6.1 Personal precautions, protective equipment and emergency procedures :**

No special measures required. Avoid formation of dust.

**6.2 Environmental precautions :** No special measures required.

**6.3 Methods and material for containment and cleaning up :** Pick up mechanically.

**6.4 Reference to other sections :**

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

## SECTION 7: Handling and storage

**7.1 Precautions for safe handling :**

No special measures required.

- **Information about fire - and explosion protection:**

No special measures required.

**7.2 Conditions for safe storage, including any incompatibilities :**

- **Storage :**

- **Requirements to be met by storerooms and receptacles :** Store away from light and moisture.


- **Information about storage in one common storage facility :** Not required.

- **Further information about storage conditions :** None

- **Recommended storage temperature :** 15-30°C

**7.3 Specific end use(s) :**

No further relevant information available.

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-216  |
|   | <b>CHROMagar™ AOLA supplement E2</b><br><b>Ref : AO882/E2 – AO883/E2</b> | Version 1.0 |
| Page 3/6  |  | 27-May-2016 |

## SECTION 8: Exposure controls / personal protection

· **Additional information about design of technical facilities** : No further data; see section 7.

### 8.1 Control parameters

· **Ingredients with limit values that require monitoring at the workplace** : Not required.

· **Additional information** :

The most current valid lists have been used as a basis for the production of this document.

### 8.2 Exposure controls

· **Personal protective equipment** :

· **General protective and hygienic measures** :

The usual precautionary measures are to be adhered to when handling chemicals.

**Respiratory protection** : In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.



**Protection of hands** : Protective gloves.  
EN374 gloves for chemicals.



**Gloves material** : Butyl rubber, BR.

**Eye protection** : Safety glasses.



**Body protection** : Protective work clothing.



## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

· **General Information**

· **Appearance** :

**Form** : Powder

**Colour** : White

· **Odour** : Odourless

· **Odour threshold** : Not determined.

· **pH-value** : Not applicable.

· **Change in condition**

**Melting point/Melting range** : Undetermined.


**Boiling point/Boiling range** : Undetermined.

· **Flash point** : Not applicable.

· **Flammability (solid, gaseous)** : Not determined.

· **Ignition temperature** : Not applicable.



|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-216  |
|   | <b>CHROMagar™ AOLA supplement E2</b><br><b>Ref : AO882/E2 – AO883/E2</b> | Version 1.0 |
| Page 4/6  |  | 27-May-2016 |


|  |   |
|--|---|
| · <b>Decomposition temperature:</b>                  | Not determined.                               |
| · <b>Self-igniting :</b>                             | Product is not self igniting.                 |
| · <b>Danger of explosion :</b>                       | Product does not present an explosion hazard. |
| · <b>Explosion limits :</b>                          |   |
| <b>Lower :</b>                                       | Not determined.                               |
| <b>Upper :</b>                                       | Not determined.                               |
| · <b>Vapour pressure :</b>                           | Not applicable.                               |
| · <b>Density at 20°C :</b>                           | 0.5 +/- 0.1 g/cm <sup>3</sup>                 |
| · <b>Relative density :</b>                          | Not determined.                               |
| · <b>Vapour density :</b>                            | Not applicable.                               |
| · <b>Evaporation rate :</b>                          | Not applicable.                               |
| · <b>Solubility in / Miscibility with water :</b>    | Fully soluble.                                |
| · <b>Segregation coefficient (n-octanol/water) :</b> | Not determined.                               |
| · <b>Viscosity:</b>                                  |   |
| <b>Dynamic :</b>                                     | Not applicable.                               |
| <b>Kinematic :</b>                                   | Not applicable.                               |
| <b>9.2 Other information :</b>                       | No further relevant information available.    |

## SECTION 10: Stability and reactivity

|   |  |
|---|--|
| <b>10.1 Reactivity :</b>                                    | No further relevant information available.                       |
| <b>10.2 Chemical stability :</b>                            |  |
| · <b>Thermal decomposition / conditions to be avoided :</b> | No decomposition if used and stored according to specifications. |
| <b>10.3 Possibility of hazardous reactions :</b>            | No dangerous reactions known.                                    |
| <b>10.4 Conditions to avoid :</b>                           | No further relevant information available.                       |
| <b>10.5 Incompatible materials :</b>                        | No further relevant information available.                       |
| <b>10.6 Hazardous decomposition products :</b>              | No dangerous decomposition products known.                       |

## SECTION 11: Toxicological information

|  |   |
|--|---|
| <b>11.1 Information on toxicological effects</b>                                     |   |
| · <b>Acute toxicity :</b>  | Based on available data, the classification criteria are not met. |
| · <b>Primary irritant effect :</b>   |   |
| · <b>On the skin :</b>   | Based on available data, the classification criteria are not met. |
| · <b>On the eyes :</b>   | Based on available data, the classification criteria are not met. |
| · <b>Respiratory or skin sensitisation :</b>   | Based on available data, the classification criteria are not met. |
| · <b>CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction) :</b> |   |
| · <b>Germ cell mutagenicity :</b>  | Based on available data, the classification criteria are not met. |
| · <b>Reproductive toxicity :</b>   | Based on available data, the classification criteria are not met. |
| · <b>STOT-single exposure :</b>  | Based on available data, the classification criteria are not met. |
| · <b>STOT-repeated exposure :</b>  | Based on available data, the classification criteria are not met. |
| · <b>Aspiration hazard :</b>   | Based on available data, the classification criteria are not met. |

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-216  |
|   | <b>CHROMagar™ AOLA supplement E2</b><br><b>Ref : AO882/E2 – AO883/E2</b> | Version 1.0 |
| Page 5/6  |  | 27-May-2016 |

## SECTION 12: Ecological information

### 12.1 Toxicity

· **Aquatic toxicity :** No further relevant information available.

**12.2 Persistence and degradability :** No further relevant information available.

**12.3 Bioaccumulative potential :** No further relevant information available.

**12.4 Mobility in soil :** No further relevant information available.

· **Additional ecological information:**

· **General notes :** Generally not hazardous for water.

### 12.5 Results of PBT and vPvB assessment

· **PBT :** Not applicable.

· **vPvB :** Not applicable.

**12.6 Other adverse effects :** No further relevant information available.

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

· **Recommendation :** Must be specially treated according to official regulations.  
Smaller quantities can be disposed of with household waste.

· **Uncleaned packaging :**

· **Recommendation :** Disposal must be made according to official regulations.

## SECTION 14: Transport information

### 14.1 UN-Number

· **ADR, ADN, IMDG, IATA :** Void

### 14.2 UN proper shipping name

· **ADR, ADN, IMDG, IATA :** Void

### 14.3 Transport hazard class(es)

· **ADR, ADN, IMDG, IATA**

· **Class :** Void

### 14.4 Packing group

· **ADR, IMDG, IATA :** Void

**14.5 Environmental hazards :** Not applicable.

**14.6 Special precautions for user :** Not applicable.

### 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code :


Not applicable.

· **UN "Model Regulation":** Void

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the mixture :

· **Directive 2012/18/EU**

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA supplement S</b><br><b>Ref : AO883/S</b> | AD-ASQ-273  |
|   | Page 1/8   | Version 1.0 |
|   |  | 11-Jan-2017 |

## SECTION 1: Identification of the mixture and of the company / undertaking

### 1.1 Product identifier

· Trade name : CHROMagar™ AOLA supplement S

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

· Application of the substance / the preparation : Reagent for Laboratory Use

### 1.3 Details of the supplier of the safety data sheet

· Manufacturer/Supplier :

CHROMagar  
 4, place du 18 Juin 1940 - 75006 Paris  
 Tel: +33 1 45 48 05 05  
 Email : chromagar@chromagar.com

### 1.4 Emergency telephone number :

France (ORFILA 24h/24) - Tel : +33 (0)1 45 42 59 59

Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112

Belgique - Tel : 32 070/245 245

UK : +44 (0)-870-190-6777 (CARECHEM24)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

· Classification according to Regulation (EC) No 1272/2008 :



GHS08 : Resp. Sens. 1 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.



GHS07 :

Acute Tox. 4 H302 Harmful if swallowed.

Eye Irrit. 2 H319 Causes serious eye irritation.

Skin Sens. 1 H317 May cause an allergic skin reaction.

### 2.2 Label elements

· Labelling according to Regulation (EC) No 1272/2008 :

The product is classified and labelled according to the CLP regulation

· Hazard pictograms :



GHS08



GHS07

· Signal word :


Danger

· Hazard-determining components of labelling :

Nalidixic acid sodium salt

Ceftazidime

Colistin sulfate

|   |   |             |
|---|---|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31) | AD-ASQ-273  |
|   | <b>CHROMagar™ AOLA supplement S</b><br><b>Ref : AO883/S</b>   | Version 1.0 |
| Page 2/8  |   | 11-Jan-2017 |

· **Hazard statements :**

- H302 : Harmful if swallowed.  
 H319 : Causes serious eye irritation.  
 H334 : May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
 H317 : May cause an allergic skin reaction.

· **Precautionary statements :**

- P284 : In case of inadequate ventilation wear respiratory protection.  
 P261 : Avoid breathing dust/fume/gas/mist/vapours/spray.  
 P280 : Wear protective gloves/protective clothing/eye protection/face protection.  
 P264 : Wash hands thoroughly after handling.  
 P270 : Do not eat, drink or smoke when using this product.  
 P305 + P351 + P338 : If in eyes rinse cautiously with water for several minutes.  
 Remove contact lenses, if present and easy to do. Continue rinsing.  
 P301 + P312 : If swallowed call a poison center/doctor if you feel unwell.  
 P304 + P340 : If inhaled remove person to fresh air and keep comfortable for breathing.  
 P333 + P313 : If skin irritation or rash occurs: Get medical advice/attention.  
 P337 + P313 : If eye irritation persists: Get medical advice/attention.  
 P302 + P352 : If on skin wash with plenty of water.  
 P501 : Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Labelling of packages where the contents do not exceed 125 ml**

· **Hazard pictograms**



GHS08



GHS07

· **Signal word :**

Danger

· **Hazard-determining components of labelling :**

Nalidixic acid sodium salt  
 Ceftazidime  
 Colistin sulfate

· **Hazard statements :**


- H334 : May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
 H317 : May cause an allergic skin reaction.

· **Precautionary statements :**

- P284 : In case of inadequate ventilation wear respiratory protection.  
 P261 : Avoid breathing dust/fume/gas/mist/vapours/spray.  
 P280 : Wear protective gloves/protective clothing/eye protection/face protection.  
 P264 : Wash hands thoroughly after handling.  
 P270 : Do not eat, drink or smoke when using this product.  
 P305 + P351 + P338 : If in eyes rinse cautiously with water for several minutes.  
 Remove contact lenses, if present and easy to do. Continue rinsing.  
 P301 + P312 : If swallowed call a poison center/doctor if you feel unwell.  
 P304 + P340 : If inhaled remove person to fresh air and keep comfortable for breathing.  
 P333 + P313 : If skin irritation or rash occurs: Get medical advice/attention.  
 P337 + P313 : If eye irritation persists: Get medical advice/attention.  
 P302 + P352 : If on skin wash with plenty of water.  
 P501 : Dispose of contents/container in accordance with local/regional/national/international regulations.

**2.3 Other hazards**

· **Results of PBT and vPvB assessment**





|  |   |             |
|--|---|-------------|
| <br>The Chromogenic Media Pioneer | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31) | AD-ASQ-273  |
|  | <b>CHROMagar™ AOLA supplement S</b><br><b>Ref : AO883/S</b>   | Version 1.0 |
| Page 3/8   |   | 11-Jan-2017 |

- **PBT** : Not applicable.
- **vPvB** : Not applicable.

### SECTION 3: Composition / information on ingredients

#### 3.2 Mixtures

- **Description** : Mixture: consisting of the following components.
- **Dangerous components** :

|                                     |   |        |
|-------------------------------------|---|--------|
| CAS: 78439-06-2                     |  Ceftazidime<br>Acute Tox. 4, H302; Eye Irrit. 2, H319   | 10-25% |
| CAS: 3374-05-8<br>EINECS: 222-159-7 |  Nalidixic acid sodium salt<br>Resp. Sens. 1, H334;  Acute Tox. 4, H302; Skin Sens. 1, H317 | 10-25% |
| CAS: 1264-72-8<br>EINECS: 215-034-3 |  Colistin sulfate<br>Acute Tox. 4, H302  | <5%    |

- **SVHC** : No
- **Additional information** : For the wording of the listed risk phrases, refer to section 16.

### SECTION 4: First aid measures

#### 4.1 Description of first aid measures

- **General information** : Position and transport stably in side position.  
In case of irregular breathing or respiratory arrest provide artificial respiration.
- **After excessive inhalation** : Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact** : Generally the product does not irritate the skin.  
In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.  
If skin irritation continues, consult a doctor.
- **After eye contact** : Rinse opened eye for several minutes under running water.
- **After swallowing** : Do not induce vomiting; call for medical help immediately.

#### 4.2 Most important symptoms and effects, both acute and delayed :

No further relevant information available.

#### 4.3 Indication of any immediate medical attention and special treatment needed :

No further relevant information available.

### SECTION 5: Firefighting measures

#### 5.1 Suitable extinguishing agents :


CO<sub>2</sub>, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

**5.2 Special hazards arising from the substance or mixture** : No further relevant information available.

#### 5.3 Advice for firefighters

- **Protective equipment** : No special measures required.

### SECTION 6: Accidental release measures

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA supplement S</b><br><b>Ref : AO883/S</b> | AD-ASQ-273  |
|   | Page 4/8   | Version 1.0 |
|   |  | 11-Jan-2017 |

### 6.1 Personal precautions, protective equipment and emergency procedures :

No special measures required. Avoid formation of dust.

### 6.2 Environmental precautions :

No special measures required.

### 6.3 Methods and material for containment and cleaning up : Pick up mechanically.

### 6.4 Reference to other sections :

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling :

No special measures required.

### · Information about fire - and explosion protection :

No special measures required.

### 7.2 Conditions for safe storage, including any incompatibilities :

#### · Storage :

· **Requirements to be met by storerooms and receptacles :** Store away from light and moisture.

· **Information about storage in one common storage facility :** Not required.

· **Further information about storage conditions :** None.

· **Recommended storage temperature :** 2-8°C

### 7.3 Specific end use(s) :

No further relevant information available.

## SECTION 8: Exposure controls / personal protection

· **Additional information about design of technical facilities :** No further data; see section 7.

### 8.1 Control parameters

#### · Ingredients with limit values that require monitoring at the workplace :

The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.

#### · Additional information :

The most current valid lists have been used as a basis for the production of this document.

### 8.2 Exposure controls

#### · Personal protective equipment :

#### · General protective and hygienic measures :

The usual precautionary measures are to be adhered to when handling chemicals.

#### Respiratory protection :

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.




#### Protection of hands :

Protective gloves.  
EN374 gloves for chemicals.



#### Gloves material :

Butyl rubber, BR

|   |   |             |
|---|---|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31) | AD-ASQ-273  |
|   | <b>CHROMagar™ AOLA supplement S</b><br><b>Ref : AO883/S</b>   | Version 1.0 |
| Page 5/8  |   | 11-Jan-2017 |

**Eye protection :** Safety glasses.



**Body protection :** Protective work clothing.



## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

· **General Information**

· **Appearance :**

**Form :** Powder

**Colour :** White

· **Odour :** Odourless

· **Odour threshold :** Not determined.

· **pH-value :** Not applicable.

· **Change in condition**

**Melting point/Melting range :** Undetermined.

**Boiling point/Boiling range :** Undetermined.

· **Flash point :** Not applicable.

· **Flammability (solid, gaseous) :** Not determined.

· **Ignition temperature :** Not applicable.

· **Decomposition temperature :** Not determined.

· **Self-igniting :** Product is not self igniting.

· **Danger of explosion :** Product does not present an explosion hazard.

· **Explosion limits :**

**Lower :** Not determined.

**Upper :** Not determined.

· **Vapour pressure :** Not applicable.

· **Density at 20°C :** 0.5 g/cm<sup>3</sup>

· **Relative density :** Not determined.

· **Vapour density :** Not applicable.

· **Evaporation rate :** Not applicable.

· **Solubility in / Miscibility with water :** Fully soluble.

· **Segregation coefficient (n-octanol/water) :** Not determined.

· **Viscosity :**

**Dynamic :** Not applicable.

**Kinematic :** Not applicable.

**9.2 Other information :** No further relevant information available.


## SECTION 10: Stability and reactivity

**10.1 Reactivity :** No further relevant information available.

**10.2 Chemical stability :**

· **Thermal decomposition / conditions to be avoided :**

No decomposition if used according to specifications.

|   |   |             |
|---|---|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31) | AD-ASQ-273  |
|   | <b>CHROMagar™ AOLA supplement S</b><br><b>Ref : AO883/S</b>   | Version 1.0 |
| Page 6/8  |   | 11-Jan-2017 |

- 10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- 10.4 Conditions to avoid :** No further relevant information available.
- 10.5 Incompatible materials :** No further relevant information available.
- 10.6 Hazardous decomposition products :** No dangerous decomposition products known.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

- **Acute toxicity :** Harmful if swallowed.
- **LD/LC50 values relevant for classification :**
  - 78439-06-2 Ceftazidime**  
Oral LD50 500 mg/kg (ATE)
  - 3374-05-8 Nalidixic acid sodium salt**  
Oral LD50 500 mg/kg (ATE)
  - 1264-72-8 Colistin sulfate**  
Oral LD50 452-1366 mg/kg (rat)
- **Primary irritant effect :**
- **On the skin :** Based on available data, the classification criteria are not met.
- **On the eyes :** Causes serious eye irritation.
- **Respiratory or skin sensitisation :** May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
May cause an allergic skin reaction.
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction) :**
- **Germ cell mutagenicity :** Based on available data, the classification criteria are not met.
- **Reproductive toxicity :** Based on available data, the classification criteria are not met.
- **STOT-single exposure :** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure :** Based on available data, the classification criteria are not met.
- **Aspiration hazard :** Based on available data, the classification criteria are not met.


## SECTION 12: Ecological information

### 12.1 Toxicity

- **Aquatic toxicity :** No further relevant information available.
- 12.2 Persistence and degradability :** No further relevant information available.
- 12.3 Bioaccumulative potential :** No further relevant information available.
- 12.4 Mobility in soil :** No further relevant information available.
- **Additional ecological information :**
- **General notes :** Generally not hazardous for water.
- 12.5 Results of PBT and vPvB assessment**
- **PBT :** Not applicable.
- **vPvB :** Not applicable.
- 12.6 Other adverse effects :** No further relevant information available.

## SECTION 13: Disposal considerations



|   |   |             |
|---|---|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31) | AD-ASQ-273  |
|   | <b>CHROMagar™ AOLA supplement S</b><br><b>Ref : AO883/S</b>   | Version 1.0 |
| Page 7/8  |   | 11-Jan-2017 |

### 13.1 Waste treatment methods

- **Recommendation :** Must be specially treated according to official regulations.  
Smaller quantities can be disposed of with household waste.
- **Uncleaned packaging :**
- **Recommendation :** Disposal must be made according to official regulations.

## SECTION 14: Transport information

### 14.1 UN-Number

- ADR, ADN, IMDG, IATA : Void

### 14.2 UN proper shipping name

- ADR, ADN, IMDG, IATA : Void

### 14.3 Transport hazard class(es)

- ADR, ADN, IMDG, IATA
- Class : Void

### 14.4 Packing group

- ADR, IMDG, IATA : Void

- 14.5 Environmental hazards : Not applicable.

- 14.6 Special precautions for user : Not applicable.

### 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code :

- Not applicable.
- UN "Model Regulation": Void

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the mixture :

- Directive 2012/18/EU
- Named dangerous substances - ANNEX I : None of the ingredients is listed.

- 15.2 Chemical safety assessment : A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information


This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

### · Relevant phrases :


- H302 : Harmful if swallowed.
- H317 : May cause an allergic skin reaction.
- H319 : Causes serious eye irritation.
- H334 : May cause allergy or asthma symptoms or breathing difficulties if inhaled.

### · Abbreviations and acronyms :

- ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
- IMDG: International Maritime Code for Dangerous Goods
- IATA: International Air Transport Association.
- GHS: Globally Harmonised System of Classification and Labelling of Chemicals

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)<br><b>CHROMagar™ AOLA supplement S</b><br><b>Ref : AO883/S</b> | AD-ASQ-273  |
|   |  | Version 1.0 |
| Page 8/8  |  | 11-Jan-2017 |

EINECS: European Inventory of Existing Commercial Chemical Substances  
 ELINCS: European List of Notified Chemical Substances  
 CAS: Chemical Abstracts Service (division of the American Chemical Society)  
 LC50: Lethal concentration, 50 percent  
 LD50: Lethal dose, 50 percent  
 PBT: Persistent, Bioaccumulative and Toxic  
 SVHC: Substances of Very High Concern  
 vPvB: very Persistent and very Bioaccumulative  
 Acute Tox. 4: Acute toxicity, Hazard Category 4  
 Eye Irrit. 2: Serious eye damage/eye irritation, Hazard Category 2  
 Resp. Sens. 1: Sensitisation - Respirat., Hazard Category 1  
 Skin Sens. 1: Sensitisation - Skin, Hazard Category 1

|   |  |             |
|---|--|-------------|
|  | <b>SAFETY DATA SHEET</b><br>(according 1907/2006/CE, Art. 31)            | AD-ASQ-216  |
|   | <b>CHROMagar™ AOLA supplement E2</b><br><b>Ref : AO882/E2 – AO883/E2</b> | Version 1.0 |
| Page 6/6  |  | 27-May-2016 |

· **Named dangerous substances - ANNEX I :** None of the ingredients is listed.

**15.2 Chemical safety assessment :** A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Abbreviations and acronyms :**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the

International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative