	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ C3G^R supplement Ref : CG632 / CG633	AD-ASQ-275
		Version 1.1
Page 1/7		09-Feb-2018

SECTION 1: Identification of the mixture and of the company / undertaking

1.1 Product identifier

· Trade name : CHROMagar™ C3G^R supplement

1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

· Application of the substance / the preparation : Reagent for Laboratory Use

1.3 Details of the supplier of the safety data sheet

· Manufacturer/Supplier :

CHROMagar
 4, place du 18 Juin 1940 - 75006 Paris
 Tel: +33 1 45 48 05 05
 Email : chromagar@chromagar.com

1.4 Emergency telephone number :

France (ORFILA 24h/24) - Tel : +33 (0)1 45 42 59 59

Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112

Belgique - Tel : 32 070/245 245

UK : 111 (NHS) / 0345 300 9923 (HSE)

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

· Classification according to Regulation (EC) No 1272/2008 :



GHS07

Acute Tox. 4 H302 Harmful if swallowed.

STOT SE 3 H335 May cause respiratory irritation.

2.2 Label elements

· Labelling according to Regulation (EC) No 1272/2008 :

The product is classified and labelled according to the CLP regulation

· Hazard pictograms :



GHS07

· Signal word :

Warning

· Hazard-determining components of labelling :

sodium 3-alpha,

12-alpha-dihydroxy-5-β-cholan-24-oate

· Hazard statements :

H302 : Harmful if swallowed.

H335 : May cause respiratory irritation.

· Precautionary statements :

P261 : Avoid breathing dust/fume/gas/mist/vapours/spray.

P264 : Wash hands thoroughly after handling.

P270 : Do not eat, drink or smoke when using this product.


P271 : Use only outdoors or in a well-ventilated area.

P301 + P312 : If swallowed call a poison center/doctor if you feel unwell.

P304 + P340 : If inhaled remove person to fresh air and keep comfortable for breathing.

P403 + P233 : Store in a well-ventilated place. Keep container tightly closed.

P501 : Dispose of contents/container in accordance with local/regional/national/international

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ C3G^R supplement Ref : CG632 / CG633	AD-ASQ-275
	Page 2/7	Version 1.1
		09-Feb-2018

regulations.

· **Additional information :**

EUH208 Contains antibiotic compounds. May produce an allergic reaction.

· **Labelling of packages where the contents do not exceed 125 ml**

· **Hazard pictograms**



GHS07

· **Signal word :**

Warning

· **Hazard-determining components of labelling :**

sodium 3-alpha,
12-alpha-dihydroxy-5-β-cholan-24-oate

· **Hazard statements :** Not classified as hazardous.

· **Precautionary statements :**

P301 + P312 : If swallowed call a poison center/doctor if you feel unwell.

P304 + P340 : If inhaled remove person to fresh air and keep comfortable for breathing.

P403 + P233 : Store in a well-ventilated place. Keep container tightly closed.

2.3 Other hazards

· **Results of PBT and vPvB assessment**

· **PBT :**

Not applicable.

· **vPvB :**

Not applicable.

SECTION 3: Composition / information on ingredients

3.2 Mixtures

· **Description :**

Mixture: consisting of the following components.

· **Dangerous components :**

CAS: 302-95-4

sodium 3-alpha,

50-100%

EINECS: 206-132-7

12-alpha-dihydroxy-5-β-cholan-24-oate



Acute Tox. 4, H302; STOT SE 3, H335

· **SVHC :**

No

· **Additional information :**

For the wording of the listed risk phrases, refer to section 16.

SECTION 4: First aid measures

4.1 Description of first aid measures

· **General information :**

Position and transport stably in side position.

In case of irregular breathing or respiratory arrest provide artificial respiration.

· **After excessive inhalation :**

Provide fresh air; consult a doctor in case of medical complaints.

· **After skin contact :**

Generally the product does not irritate the skin.

In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.

If skin irritation continues, consult a doctor.


· **After eye contact :**

Rinse opened eye for several minutes under running water.

· **After swallowing :**

Do not induce vomiting; call for medical help immediately.

4.2 Most important symptoms and effects, both acute and delayed :

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ C3G^R supplement Ref : CG632 / CG633	AD-ASQ-275
		Version 1.1
Page 3/7		09-Feb-2018

No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed :

No further relevant information available.

SECTION 5: Firefighting measures

5.1 Suitable extinguishing agents :

CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

5.2 Special hazards arising from the substance or mixture : No further relevant information available.

5.3 Advice for firefighters

· **Protective equipment :** No special measures required.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures :

No special measures required. Avoid formation of dust.

6.2 Environmental precautions :

Do not allow product to reach sewage system or any water course. Dilute with plenty of water.

6.3 Methods and material for containment and cleaning up : Pick up mechanically.

6.4 Reference to other sections :

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

SECTION 7: Handling and storage

7.1 Precautions for safe handling :

Keep away from heat and direct sunlight.
Thorough dedusting.

· **Information about fire - and explosion protection :**

No special measures required.

7.2 Conditions for safe storage, including any incompatibilities :

· **Storage :**

· **Requirements to be met by storerooms and receptacles :** Store away from light and moisture. No special requirements.

· **Information about storage in one common storage facility :** Not required.

· **Further information about storage conditions :**

Store in cool, dry conditions in well sealed receptacles.

· **Recommended storage temperature :**

2-8°C

7.3 Specific end use(s) :


No further relevant information available.

SECTION 8: Exposure controls / personal protection

· **Additional information about design of technical facilities :** No further data; see section 7.

8.1 Control parameters

· **Ingredients with limit values that require monitoring at the workplace :**

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ C3G^R supplement Ref : CG632 / CG633	AD-ASQ-275
		Version 1.1
Page 4/7		09-Feb-2018

The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.

· **Additional information :**

The most current valid lists have been used as a basis for the production of this document.

8.2 Exposure controls

· **Personal protective equipment :**

· **General protective and hygienic measures :**

The usual precautionary measures are to be adhered to when handling chemicals.

Respiratory protection : In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.



Protection of hands : Protective gloves.
EN374 gloves for chemicals.



Gloves material : Butyl rubber, BR

Eye protection : Safety glasses.



Body protection : Protective work clothing.



SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

· **General Information**

· **Appearance :**

Form : Powder

Colour : Yellow

· **Odour :** Odourless

· **Odour threshold :** Not determined.

· **pH-value :** Not applicable.

· **Change in condition**

Melting point/Melting range : Undetermined.

Boiling point/Boiling range : Undetermined.

· **Flash point :** Not applicable.

· **Flammability (solid, gaseous) :** Not determined.

· **Ignition temperature :** Not applicable.


· **Decomposition temperature :** Not determined.

· **Self-igniting :** Product is not self igniting.

· **Danger of explosion :** Product does not present an explosion hazard.

· **Explosion limits :**

Lower : Not determined.

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ C3G^R supplement Ref : CG632 / CG633	AD-ASQ-275
	Page 5/7	Version 1.1
		09-Feb-2018

Upper :	Not determined.
· Vapour pressure :	Not applicable.
· Density at 20°C :	0.5 +/- 0.1 g/cm ³
· Relative density :	Not determined.
· Vapour density :	Not applicable.
· Evaporation rate :	Not applicable.
· Solubility in / Miscibility with water :	Partly soluble.
· Segregation coefficient (n-octanol/water) :	Not determined.
· Viscosity :	
Dynamic :	Not applicable.
Kinematic :	Not applicable.
9.2 Other information :	No further relevant information available.

SECTION 10: Stability and reactivity


10.1 Reactivity :	No further relevant information available.
10.2 Chemical stability :	
· Thermal decomposition / conditions to be avoided :	No decomposition if used according to specifications.
10.3 Possibility of hazardous reactions :	No dangerous reactions known.
10.4 Conditions to avoid :	No further relevant information available.
10.5 Incompatible materials :	No further relevant information available.
10.6 Hazardous decomposition products :	No dangerous decomposition products known.

SECTION 11: Toxicological information

11.1 Information on toxicological effects	
· Acute toxicity :	Harmful if swallowed.
· LD/LC50 values relevant for classification:	302-95-4 sodium 3-alpha,12-alpha-dihydroxy-5-beta-cholan-24-oate
	Oral LD50 500 mg/kg (ATE)
· Primary irritant effect :	
· On the skin :	Based on available data, the classification criteria are not met.
· On the eyes :	Based on available data, the classification criteria are not met.
· Respiratory or skin sensitisation :	Based on available data, the classification criteria are not met.
· CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction) :	
· Germ cell mutagenicity :	Based on available data, the classification criteria are not met.
· Reproductive toxicity :	Based on available data, the classification criteria are not met.
· STOT-single exposure :	Based on available data, the classification criteria are not met.
· STOT-repeated exposure :	Based on available data, the classification criteria are not met.
· Aspiration hazard :	Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

12.1 Toxicity	
· Aquatic toxicity :	No further relevant information available.

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ C3G^R supplement Ref : CG632 / CG633	AD-ASQ-275
		Version 1.1
Page 6/7		09-Feb-2018

12.2 Persistence and degradability :	No further relevant information available.
12.3 Bioaccumulative potential :	No further relevant information available.
12.4 Mobility in soil :	No further relevant information available.
· Additional ecological information :	
· General notes :	Generally not hazardous for water.
12.5 Results of PBT and vPvB assessment	
· PBT :	Not applicable.
· vPvB :	Not applicable.
12.6 Other adverse effects :	No further relevant information available.

SECTION 13: Disposal considerations


13.1 Waste treatment methods	
· Recommendation :	Must be specially treated according to official regulations..
· Uncleaned packaging :	
· Recommendation :	Disposal must be made according to official regulations.

SECTION 14: Transport information

14.1 UN-Number	
· ADR, ADN, IMDG, IATA :	Void
14.2 UN proper shipping name	
· ADR, ADN, IMDG, IATA :	Void
14.3 Transport hazard class(es)	
· ADR, ADN, IMDG, IATA	
· Class :	Void
14.4 Packing group	
· ADR, IMDG, IATA :	Void
14.5 Environmental hazards :	Not applicable.
14.6 Special precautions for user :	Not applicable.
14.7 Transport in bulk according to Annex II of Marpol and the IBC Code :	Not applicable.
· UN "Model Regulation":	Void

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the mixture :	
· Directive 2012/18/EU	
· Named dangerous substances - ANNEX I :	None of the ingredients is listed.
15.2 Chemical safety assessment :	A Chemical Safety Assessment has not been carried out

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ C3G^R supplement Ref : CG632 / CG633	AD-ASQ-275
		Version 1.1
Page 7/7		09-Feb-2018

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases :**

H302: Harmful if swallowed.

H335: May cause respiratory irritation.

· **Abbreviations and acronyms :**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association.

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity, Hazard Category 4

Skin Irrit. 2: Skin corrosion/irritation, Hazard Category 2

Eye Irrit. 2: Serious eye damage/eye irritation, Hazard Category 2

Resp. Sens. 1: Sensitisation - Respirat., Hazard Category 1

Skin Sens. 1: Sensitisation - Skin, Hazard Category 1

STOT SE 3: Specific target organ toxicity - Single exposure, Hazard Category 3