	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria base</b> <b>Ref : LM852(B)</b>	AD-ASQ-280
		Version 2.0
		27-Mar-2019
Page 1/7		

## SECTION 1: Identification of the mixture and of the company / undertaking

### 1.1 Product identifier

- Trade name : CHROMagar™ Listeria base

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

- Application of the substance / the preparation : Reagent for Laboratory Use

### 1.3 Details of the supplier of the safety data sheet

- Manufacturer/Supplier :

CHROMagar  
 4, place du 18 Juin 1940 - 75006 Paris  
 Tel: +33 1 45 48 05 05  
 Email : chromagar@chromagar.com

### 1.4 Emergency telephone number :

France (ORFILA 24h/24) - Tel : +33 (0)1 45 42 59 59  
 Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
 Belgique - Tel : 32 070/245 245  
 UK : 111 (NHS) / 0345 300 9923 (HSE)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

- Classification according to Regulation (EC) No 1272/2008 :



GHS07 : Eye Irrit. 2 H319 Causes serious eye irritation.

### 2.2 Label elements

- Labelling according to Regulation (EC) No 1272/2008 :

The product is classified and labelled according to the CLP regulation.

- Hazard pictograms :



GHS07

- Signal word :

Warning

- Hazard statements :

H319 : Causes serious eye irritation

- Precautionary statements :

P280 : Wear protective gloves/protective clothing/eye protection/face protection.

P264 : Wash hands thoroughly after handling.

P305 + P351 + P338 : If in eyes rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 : If eye irritation persists: Get medical advice/attention.

- Labelling of packages where the contents do not exceed 125 ml

- Hazard pictograms




GHS07

- Signal word :

Warning

- Hazard statements :

Not classified as hazardous.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria base</b> <b>Ref : LM852(B)</b>	AD-ASQ-280
		Version 2.0
		27-Mar-2019
Page 2/7		

### 2.3 Other hazards


#### • Results of PBT and vPvB assessment

- PBT : Not applicable.
- vPvB : Not applicable.

## SECTION 3: Composition / information on ingredients

### 3.2 Mixtures

- **Description :** Mixture: consisting of the following components.
- **Dangerous components :**

CAS: 7447-41-8	Lithium chloride	10-25%
EINECS: 231-212-3	Acute Tox. 4, H302; Acute Tox. 4, H312;	
RTECS : OJ 5950000	 Eye Irrit. 2, H319; STOT SE 3, H335	
Reg.nr.: 01-2119560574-35		

- **SVHC :** No
- **Additional information :** For the wording of the listed risk phrases, refer to section 16.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

- **General information :** Position and transport stably in side position.  
In case of irregular breathing or respiratory arrest provide artificial respiration.
- **After excessive inhalation :** Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact :** Generally the product does not irritate the skin.  
In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.  
If skin irritation continues, consult a doctor.
- **After eye contact :** Rinse opened eye for several minutes under running water.
- **After swallowing :** Do not induce vomiting; call for medical help immediately.

### 4.2 Most important symptoms and effects, both acute and delayed :

No further relevant information available.

### 4.3 Indication of any immediate medical attention and special treatment needed :

No further relevant information available.

## SECTION 5: Firefighting measures

### 5.1 Suitable extinguishing agents :

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.


### 5.2 Special hazards arising from the substance or mixture : No further relevant information available.

### 5.3 Advice for firefighters

- **Protective equipment :** No special measures required.

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures :

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria base</b> <b>Ref : LM852(B)</b>	AD-ASQ-280
		Version 2.0
		27-Mar-2019
Page 3/7		

No special measures required. Avoid formation of dust.

## 6.2 Environmental precautions :

Do not allow product to reach sewage system or any water course. Dilute with plenty of water.

**6.3 Methods and material for containment and cleaning up :** Pick up mechanically.

## 6.4 Reference to other sections :

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

## SECTION 7: Handling and storage

- 7.1 Precautions for safe handling :** Keep away from heat and direct sunlight.  
Thorough dedusting.
- **Information about fire - and explosion protection :** No special measures required.
- 7.2 Conditions for safe storage, including any incompatibilities :**
- **Storage :**
- **Requirements to be met by storerooms and receptacles :** Store away from light and moisture.  
No special requirements.
- **Information about storage in one common storage facility :** Not required.
- **Further information about storage conditions :** Store in cool,  
dry conditions in well sealed receptacles.
- **Recommended storage temperature :** 2-30°C
- 7.3 Specific end use(s) :** No further relevant information available.

## SECTION 8: Exposure controls / personal protection

- **Additional information about design of technical facilities :** No further data; see section 7.

### 8.1 Control parameters

- **Ingredients with limit values that require monitoring at the workplace :**

7447-41-8 lithium chloride (10-25%)

MAK (Switzerland) Short-term value: 0.2 e mg/m³

Long-term value: 0.2 e mg/m³

SSc;als Li

- **Additional information :**

The most current valid lists have been used as a basis for the production of this document.


### 8.2 Exposure controls

- **Personal protective equipment :**
- **General protective and hygienic measures :**

The usual precautionary measures are to be adhered to when handling chemicals.

**Respiratory protection :** In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.



	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria base</b> <b>Ref : LM852(B)</b>	AD-ASQ-280
		Version 2.0
		27-Mar-2019

**Protection of hands :**

Protective gloves.  
EN374 gloves for chemicals.



**Gloves material :**

Butyl rubber, BR

**Eye protection :**

Safety glasses.



**Body protection :**

Protective work clothing.



## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

#### · General Information

#### · Appearance :

**Form :**

Powder

**Colour :**

Beige

#### · Odour :

Characteristic.

#### · Odour threshold :

Not determined.

#### · pH-value :

Not applicable.

#### · Change in condition

**Melting point/Melting range :**

Undetermined.

**Boiling point/Boiling range :**

Undetermined.

#### · Flash point :

Not applicable.

#### · Flammability (solid, gaseous) :

Not determined.

#### · Ignition temperature :

Not applicable.

#### · Decomposition temperature :

Not determined.

#### · Self-igniting :

Product is not self igniting.

#### · Danger of explosion :

Product does not present an explosion hazard.

#### · Explosion limits :

**Lower :**

Not determined.

**Upper :**

Not determined.

#### · Vapour pressure :

Not applicable.

#### · Density at 20°C :

0.5 +/- 0.1 g/cm<sup>3</sup>

#### · Relative density :

Not determined.

#### · Vapour density :

Not applicable.

#### · Evaporation rate :

Not applicable.

#### · Solubility in / Miscibility with water :

Partly soluble.

#### · Segregation coefficient (n-octanol/water) :

Not determined.

#### · Viscosity :

**Dynamic :**


Not applicable.

**Kinematic :**

Not applicable.

### 9.2 Other information :

No further relevant information available.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria base</b> <b>Ref : LM852(B)</b>	AD-ASQ-280
		Version 2.0
		27-Mar-2019
Page 5/7		

## SECTION 10: Stability and reactivity

- 10.1 Reactivity :** No further relevant information available.
- 10.2 Chemical stability :**
- **Thermal decomposition / conditions to be avoided :** No decomposition if used according to specifications.
- 10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- 10.4 Conditions to avoid :** No further relevant information available.
- 10.5 Incompatible materials :** No further relevant information available.
- 10.6 Hazardous decomposition products :** No dangerous decomposition products known.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

- **Acute toxicity :** Based on available data, the classification criteria are not met.
- **LD/LC50 values relevant for classification :**

7447-41-8 lithium chloride

Oral LD50 526 mg/kg (rat)

Dermal LD50 1100 mg/kg (ATE)

- **Primary irritant effect :**
- **On the skin :** Based on available data, the classification criteria are not met.
- **On the eyes :** Causes serious eye irritation.
- **Respiratory or skin sensitisation :** Based on available data, the classification criteria are not met.
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction) :**
- **Germ cell mutagenicity :** Based on available data, the classification criteria are not met.
- **Reproductive toxicity :** Based on available data, the classification criteria are not met.
- **STOT-single exposure :** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure :** Based on available data, the classification criteria are not met.
- **Aspiration hazard :** Based on available data, the classification criteria are not met.

## SECTION 12: Ecological information

### 12.1 Toxicity

- **Aquatic toxicity :** No further relevant information available.

**12.2 Persistence and degradability :** No further relevant information available.

**12.3 Bioaccumulative potential :** No further relevant information available.

**12.4 Mobility in soil :** No further relevant information available.


- **Additional ecological information :**

• **General notes :** Generally not hazardous for water.

### 12.5 Results of PBT and vPvB assessment

- **PBT :** Not applicable.
- **vPvB :** Not applicable.

**12.6 Other adverse effects :** No further relevant information available.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria base</b> <b>Ref : LM852(B)</b>	AD-ASQ-280
		Version 2.0
		27-Mar-2019
Page 6/7		

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

- **Recommendation :** Must be specially treated according to official regulations..
- **Uncleaned packaging :**
- **Recommendation :** Disposal must be made according to official regulations.

## SECTION 14: Transport information

### 14.1 UN-Number

- ADR, ADN, IMDG, IATA : Void

### 14.2 UN proper shipping name

- ADR, ADN, IMDG, IATA : Void

### 14.3 Transport hazard class(es)

- ADR, ADN, IMDG, IATA
- **Class :** Void

### 14.4 Packing group

- ADR, IMDG, IATA : Void

### 14.5 Environmental hazards :

Not applicable.

### 14.6 Special precautions for user :

Not applicable.

### 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code :

Not applicable.

- **UN "Model Regulation":** Void

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the mixture :

- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I :** None of the ingredients is listed.

### 15.2 Chemical safety assessment :

A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information


This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

### · **Relevant phrases :**

H302 : Harmful if swallowed.  
H312 : Harmful in contact with skin.  
H319 : Causes serious eye irritation.  
H335 : May cause respiratory irritation.

### · **Abbreviations and acronyms :**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria base</b> <b>Ref : LM852(B)</b>	AD-ASQ-280
		Version 2.0
		27-Mar-2019
Page 7/7		

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic


SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity, Hazard Category 4

Eye Irrit. 2: Serious eye damage/eye irritation, Hazard Category 2

STOT SE 3: Specific target organ toxicity - Single exposure, Hazard Category 3

	<b>SAFETY DATA SHEET</b> (according to 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria supplement</b> <b>Ref : LM852(S)</b>	AD-ASQ-234
		Version 2.0
Page 1/6		27-Mar-2019

## SECTION 1: Identification of the mixture and of the company / undertaking

### 1.1 Product identifier

- Trade name: CHROMagar™ Listeria supplement

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

- Application of the substance / the preparation : Reagent for Laboratory Use

### 1.3 Details of the supplier of the safety data sheet

- Manufacturer/Supplier:

CHROMagar  
4, place du 18 Juin 1940 - 75006 Paris  
Tel: +33 1 45 48 05 05  
Email : chromagar@chromagar.com

### 1.4 Emergency telephone number:

France (ORFILA 24h/24) - Tel : +33 (0)1 45 42 59 59  
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
Belgique - Tel : 32 070/245 245  
UK : 111 (NHS) / 0345 300 9923 (HSE)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

- Classification according to Regulation (EC) No 1272/2008 :

The product is not classified according to the CLP regulation.

### 2.2 Label elements

- Labelling according to Regulation (EC) No 1272/2008 : Not classified as hazardous
- Hazard pictograms : Not classified as hazardous
- Signal word : Not classified as hazardous
- Hazard statements : Not classified as hazardous

### 2.3 Other hazards

- Results of PBT and vPvB assessment
- PBT: Not applicable.
- vPvB: Not applicable.

## SECTION 3: Composition / information on ingredients

### 3.2 Mixtures


- Description : Mixture: consisting of the following components.
- Dangerous components : Not classified as hazardous.
- SVHC : No
- Additional information : For the wording of the listed risk phrases, refer to section 16.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

- General information : Position and transport stably in side position.  
In case of irregular breathing or respiratory arrest provide artificial respiration.



	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria supplement</b> <b>Ref : LM852(S)</b>	AD-ASQ-234
		Version 2.0
Page 2/6		27-Mar-2019

- **After excessive inhalation :** Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact :** Generally the product does not irritate the skin.  
In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.  
If skin irritation continues, consult a doctor.
- **After eye contact :** Rinse opened eye for several minutes under running water.
- **After swallowing :** Do not induce vomiting; call for medical help immediately.

#### 4.2 Most important symptoms and effects, both acute and delayed :

No further relevant information available.

#### 4.3 Indication of any immediate medical attention and special treatment needed :

No further relevant information available.

### SECTION 5: Firefighting measures

#### 5.1 Suitable extinguishing agents:

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

#### 5.2 Special hazards arising from the substance or mixture : No further relevant information available.

#### 5.3 Advice for firefighters

- **Protective equipment :** No special measures required.

### SECTION 6: Accidental release measures

#### 6.1 Personal precautions, protective equipment and emergency procedures :

No special measures required. Avoid formation of dust.

#### 6.2 Environmental precautions : Do not allow product to reach sewage system or any water course. Dilute with plenty of water.

#### 6.3 Methods and material for containment and cleaning up : Pick up mechanically.

#### 6.4 Reference to other sections :

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

### SECTION 7: Handling and storage

#### 7.1 Precautions for safe handling :

Keep away from heat and direct sunlight.  
Thorough dedusting.

#### • Information about fire - and explosion protection:

No special measures required.

#### 7.2 Conditions for safe storage, including any incompatibilities :

##### • Storage :

- **Requirements to be met by storerooms and receptacles :** Store away from light and moisture.  
No special requirements.

##### • Information about storage in one common storage facility :

Not required.

##### • Further information about storage conditions :


Store in cool, dry conditions in well sealed receptacles.

##### • Recommended storage temperature :

2-8°C

#### 7.3 Specific end use(s) :

No further relevant information available.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria supplement</b> <b>Ref : LM852(S)</b>	AD-ASQ-234
		Version 2.0
		27-Mar-2019
Page 3/6		

## SECTION 8: Exposure controls / personal protection

- **Additional information about design of technical facilities :** No further data; see section 7.

### 8.1 Control parameters

- **Ingredients with limit values that require monitoring at the workplace :**

The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.

- **Additional information :**

The most current valid lists have been used as a basis for the production of this document.

### 8.2 Exposure controls

- **Personal protective equipment :**
- **General protective and hygienic measures :**

The usual precautionary measures are to be adhered to when handling chemicals.

**Respiratory protection :** In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.



**Protection of hands :** Protective gloves.  
EN374 gloves for chemicals.



**Gloves material :** Butyl rubber, BR.  
**Eye protection :** Safety glasses.




**Body protection :** Protective work clothing.



## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

- **General Information**
- **Appearance :**
  - Form :** Powder
  - Colour :** Yellow
- **Odour :** Characteristic
- **Odour threshold :** Not determined.
- **pH-value :** Not applicable.
- **Change in condition**
  - Melting point/Melting range :** Undetermined.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria supplement</b> <b>Ref : LM852(S)</b>	AD-ASQ-234
		Version 2.0
		27-Mar-2019
Page 4/6		


<b>Boiling point/Boiling range :</b>	Undetermined.
• <b>Flash point :</b>	Not applicable.
• <b>Flammability (solid, gaseous) :</b>	Not determined.
• <b>Ignition temperature:</b>	Not applicable.
• <b>Decomposition temperature:</b>	Not determined.
• <b>Self-igniting :</b>	Product is not self igniting.
• <b>Danger of explosion :</b>	Product does not present an explosion hazard.
• <b>Explosion limits :</b>	
<b>Lower :</b>	Not determined.
<b>Upper :</b>	Not determined.
• <b>Vapour pressure :</b>	Not applicable.
• <b>Density at 20°C :</b>	0.5 +/- 0.1 g/cm³
• <b>Relative density :</b>	Not determined.
• <b>Vapour density :</b>	Not applicable.
• <b>Evaporation rate :</b>	Not applicable.
• <b>Solubility in / Miscibility with water :</b>	Partly soluble.
• <b>Segregation coefficient (n-octanol/water) :</b>	Not determined.
• <b>Viscosity:</b>	
<b>Dynamic :</b>	Not applicable.
<b>Kinematic :</b>	Not applicable.
<b>9.2 Other information :</b>	No further relevant information available.

## SECTION 10: Stability and reactivity

<b>10.1 Reactivity :</b>	No further relevant information available.
<b>10.2 Chemical stability :</b>	
• <b>Thermal decomposition / conditions to be avoided :</b>	No decomposition if used and stored according to specifications.
<b>10.3 Possibility of hazardous reactions :</b>	No dangerous reactions known.
<b>10.4 Conditions to avoid :</b>	No further relevant information available.
<b>10.5 Incompatible materials :</b>	No further relevant information available.
<b>10.6 Hazardous decomposition products :</b>	No dangerous decomposition products known.

## SECTION 11: Toxicological information

<b>11.1 Information on toxicological effects</b>	
• <b>Acute toxicity :</b>	Based on available data, the classification criteria are not met.
• <b>Primary irritant effect :</b>	
• <b>On the skin :</b>	Based on available data, the classification criteria are not met.
• <b>On the eyes :</b>	Based on available data, the classification criteria are not met.
• <b>Respiratory or skin sensitisation :</b>	Based on available data, the classification criteria are not met.
• <b>CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction) :</b>	
• <b>Germ cell mutagenicity :</b>	Based on available data, the classification criteria are not met.
• <b>Reproductive toxicity :</b>	Based on available data, the classification criteria are not met.
• <b>STOT-single exposure :</b>	Based on available data, the classification criteria are not met.
• <b>STOT-repeated exposure :</b>	Based on available data, the classification criteria are not met.
• <b>Aspiration hazard :</b>	Based on available data, the classification criteria are not met.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria supplement</b> <b>Ref : LM852(S)</b>	AD-ASQ-234
		Version 2.0
		27-Mar-2019
Page 5/6		

## SECTION 12: Ecological information

### 12.1 Toxicity

· **Aquatic toxicity :** No further relevant information available.

**12.2 Persistence and degradability :** No further relevant information available.

**12.3 Bioaccumulative potential :** No further relevant information available.

**12.4 Mobility in soil :** No further relevant information available.

· **Additional ecological information:**

· **General notes :** Generally not hazardous for water.

### 12.5 Results of PBT and vPvB assessment

· **PBT :** Not applicable.

· **vPvB :** Not applicable.

**12.6 Other adverse effects :** No further relevant information available.

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

· **Recommendation :** Must be specially treated according to official regulations.

· **Uncleaned packaging :**

· **Recommendation :** Disposal must be made according to official regulations.

## SECTION 14: Transport information

### 14.1 UN-Number

· **ADR, ADN, IMDG, IATA :** Void

### 14.2 UN proper shipping name

· **ADR, ADN, IMDG, IATA :** Void

### 14.3 Transport hazard class(es)

· **ADR, ADN, IMDG, IATA**

· **Class :** Void

### 14.4 Packing group

· **ADR, IMDG, IATA :** Void

**14.5 Environmental hazards :** Not applicable.

**14.6 Special precautions for user :** Not applicable.

**14.7 Transport in bulk according to Annex II of Marpol and the IBC Code :**

Not applicable.


· **UN "Model Regulation":** Void

## SECTION 15: Regulatory information

**15.1 Safety, health and environmental regulations/legislation specific for the mixture :**

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I :** None of the ingredients is listed.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ Listeria supplement</b> <b>Ref : LM852(S)</b>	AD-ASQ-234
		Version 2.0
Page 6/6		27-Mar-2019

**15.2 Chemical safety assessment :** A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

### • Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative