



DETECTION OF EXTENDED-SPECTRUM β -LACTAMASE PRODUCING ENTEROBACTERIACEAE



BY CHROMOGENIC ESBL SELECTIVE MEDIUM

G. Klyasova, A. Korobova, I. Frolova

National Research Center for Hematology, Moscow, Russia

EV0393


OBJECTIVES

The aim of the study was to evaluate the detection of extended-spectrum beta-lactamase producing *Enterobacteriaceae* (ESBL-E) by chromogenic ESBL selective medium and to compare results with double disk synergy test.

METHOD

Prospective study was performed from April 2013 to December 2013.

Study design was as follows:



Additional tests

- AmpC detection
Strains sensitive to Cefepime and resistant to Cefoxitin → E-test (Cefotetan and Cefotetan with Cloxacillin)
- Detection of TEM-1 and CTX-M by real-time PCR

RESULTS

We collected 1552 swabs and 1243 *Enterobacteriaceae* isolates were obtained. A total of 394 ESBL-E strains were recovered (Fig. 2), of those 123 (31%) only on chromogenic ESBL selective medium, 263 (67%) on standard medium and on CHROMagar™ESBL, 8 (2%) only on standard medium ($p < 0,0001$, Fig. 1). Characteristics of ESBL-E are summarized in Table 1. CTX-M group of ESBL were detected in the majority of ESBL-E.

Production of ESBL was confirmed for 94% (386/409) isolates obtained on chromogenic ESBL selective, other 6% (23/409) isolates were non-ESBL-E (Tab. 2). Sensitivity and specificity of CHROMagar™ESBL was 98% and 97%, respectively (Tab. 3).

CONCLUSION

Chromogenic ESBL selective medium had high sensitivity (98%) and specificity (97%) for the detection of ESBL-E and might be used in routine laboratory practice. Detection of ESBL-E from the same rectal swabs by chromogenic ESBL selective medium was significantly higher than by standard medium (98% vs. 69%, $p < 0,0001$) and data were reported to clinical department in 18-24 hours.

RESULTS

Fig. 1. ESBL-E detection




Fig. 2. Distribution of ESBL-E (N = 394)




Table 1. Characteristics of ESBL-E

ESBL-E	CTX-M, n (%)	TEM, n (%)	TEM+ CTX-M, n (%)
<i>E. coli</i>	157 (70)	107 (48)	64 (28)
<i>K. pneumoniae</i>	51 (47)	36 (33)	10 (9)
<i>Enterobacter</i> spp.	20 (80)	10 (40)	7 (28)
<i>Citrobacter</i> spp.	14 (78)	8 (44)	7 (39)
<i>K. oxytoca</i>	7 (70)	3 (30)	2 (20)
<i>P. mirabilis</i>	1 (25)	4 (100)	1 (25)
<i>Raoultella ornithinolytica</i>	2 (67)	1 (33)	1 (33)
Total	252 (64)	169 (43)	92 (23)

Table 2. False positive detection of ESBL on CHROMagar™ESBL (n=23)

Strains (N)	Species	n (%)
ampC (+) strains (N=15)	<i>Enterobacter</i> spp.	11 (73)
	<i>Citrobacter</i> spp.	2 (13)
	<i>M. morgani</i>	1 (7)
	<i>E. coli</i>	1 (7)
Strains sensitive to cephalosporins (N=8)	<i>E. coli</i>	3 (37,5)
	<i>K. oxytoca</i>	2 (25)
	<i>P. vulgaris</i>	2 (25)
	<i>Citrobacter</i> spp.	1 (12,5)

Table 3. Sensitivity and specificity of CHROMagar™ESBL medium

Species of ESBL-E	Sensitivity	Specificity
<i>E. coli</i>	99%	99%
<i>K. pneumoniae</i>	96%	100%
Other <i>Enterobacteriaceae</i>	98%	95%
Total	98%	97%

National Research Center For Hematology,
Novozykovski pr., 4, Moscow, 125167, Russia
E-mail: klyasova.g@blood.ru,
amalofeeva@yandex.ru