

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus base Ref : SA16B	AD-ASQ-260
		Version 2.2
		10-Jul-2026
Page 1/7		

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- Trade name: **CHROMagar™ Salmonella Plus base**
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
- **Sector of Use** SU3 Industrial uses: Uses of substances as such or in preparations on industrial sites
- **Application of the substance / the preparation** : Reagent for Laboratory Use
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
 CHROMagar
 29, avenue George Sand – 93210 La Plaine Saint-Denis
 Tel: +33 1 45 48 05 05
 Email : chromagar@chromagar.com
- **1.4 Emergency telephone number:**
 France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59
 Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112
 Belgium - Brüssel: +32 070/245 245
 UK : Call NHS 999 in a medical emergency. If you urgently need medical help or advice but it's not a life-threatening situation, call 111 (24/7). The NHS 111 number is currently only available in certain areas. If you're outside of these areas, you should call NHS Direct on 0845 4647
 Swiss Toxicological Information Centre (Zürich) : +41 44 251 51 51 / www.toxi.ch.

SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**
 The product is not classified, according to the CLP regulation.
- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Not classified as hazardous
- **Hazard pictograms** Not classified as hazardous
- **Signal word** Not classified as hazardous
- **Hazard statements** Not classified as hazardous
- **Additional information:**
 EUH210 Safety data sheet available on request.
- **2.3 Other hazards**
 The mixture shall not contain substances listed in accordance with Article 59(1) of REACH as having endocrine disrupting properties, or be recognised as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 in a concentration equal to or greater than 0.1%.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

- **3.2 Mixtures**

 The Chromogenic Media Pioneer	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus base Ref : SA16B	AD-ASQ-260
		Version 2.2
		10-Jul-2026
Page 2/7		

- **Description:** Preparation: composed of the substances listed below.
- **Dangerous components:**
 - CAS: 302-95-4 sodium 3- α ,12- α -dihydroxy-5- β -cholan-24-oate < 5%
 - EINECS: 206-132-7  Acute Tox. 4, H302; STOT SE 3, H335
- **SVHC** -
- **Additional information:** For the wording of the listed risk phrases, refer to section 16.

SECTION 4: First aid measures

- **4.1 Description of first aid measures**
- **General information:**
 - Position and transport stably in side position.
 - In case of irregular breathing or respiratory arrest provide artificial respiration.
- **After inhalation:** Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact:**
 - Generally the product does not irritate the skin.
 - In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.
 - If skin irritation continues, consult a doctor.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** Do not induce vomiting; call for medical help immediately.
- **4.2 Most important symptoms and effects, both acute and delayed :**
 - No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed :**
 - No further relevant information available.

SECTION 5: Firefighting measures

- **5.1 Suitable extinguishing agents:**
 - CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.
- **5.3 Advice for firefighters**
- **Protective equipment:** No special measures required.

SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures :**
 - No special measures required.
 - Avoid formation of dust.
- **6.2 Environmental precautions:**
 - Do not allow product to reach sewage system or any water course.
 - Dilute with plenty of water.
 - Do not allow the product to enter sewers/ surface water or ground water.
- **6.3 Methods and material for containment and cleaning up:** Pick up mechanically.
- **6.4 Reference to other sections :**
 - See Section 7 for information on safe handling.
 - See Section 8 for information on personal protective equipment.
 - See Section 13 for disposal information.

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus base Ref : SA16B	AD-ASQ-260
Page 3/7		Version 2.2
		10-Jul-2026

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :**
Keep away from heat and direct sunlight.
Thorough dedusting.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:**
Store away from light and moisture.
No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Store in cool, dry conditions in well sealed receptacles.
- **Recommended storage temperature:** 15-30 °C
- **7.3 Specific end use(s) :** No further relevant information available.

SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**
- **Ingredients with limit values that require monitoring at the workplace:**
The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:**
The most current valid lists have been used as a basis for the production of this document.
- **8.2 Exposure controls**
- **Appropriate engineering controls** No further data; see section 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:**



In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

Filter A/P2

- **Hand protection**



Protective gloves

Chemical resistant gloves EN 374

- **Gloves material :** Butyl rubber, BR
- **Eye/face protection**



Safety glasses

- **Body protection:**



Protective work clothing

 The Chromogenic Media Pioneer	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus base Ref : SA16B	AD-ASQ-260
		Version 2.2
		10-Jul-2026
Page 4/7		

SECTION 9: Physical and chemical properties

· 9.1 Information on basic physical and chemical properties

· General Information

· Physical state	Solid
· Form:	Powder
· Colour:	Beige
· Odour:	Characteristic
· Odour threshold:	Not determined.
· Melting point/freezing point:	Undetermined.
· Boiling point or initial boiling point and boiling range	Undetermined.
· Flammability	Not determined.
· Lower and upper explosion limit	
· Lower:	Not determined.
· Upper:	Not determined.
· Flash point:	Not applicable.
· Auto-ignition temperature:	Not applicable.
· Decomposition temperature:	Not determined.
· pH	Not applicable.
· Viscosity:	
· Kinematic viscosity	Not applicable.
· Kinematic:	Not applicable.
· Dynamic:	Not applicable.
· Solubility	
· water:	Partly soluble.
· Partition coefficient n-octanol/water (log value)	Not determined.
· Vapour pressure:	Not applicable.
· Density and/or relative density	
· Density at 20 °C:	0.5+/-0.1 g/cm ³
· Relative density :	Not determined.
· Vapour density :	Not applicable.
· Particle characteristics	
See section 3.	

· 9.2 Other information :

· Important information on protection of health and environment, and on safety.

· Ignition temperature:	Product is not self igniting.
· Explosive properties:	Product does not present an explosion hazard.
· Change in condition	
· Evaporation rate :	Not applicable.

· Information with regard to physical hazard classes

· Explosives	Not classified as hazardous
· Flammable gases	Not classified as hazardous
· Aerosols	Not classified as hazardous
· Oxidising gases	Not classified as hazardous
· Gases under pressure	Not classified as hazardous
· Flammable liquids	Not classified as hazardous
· Flammable solids	Not classified as hazardous

 The Chromogenic Media Pioneer	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus base Ref : SA16B	AD-ASQ-260
		Version 2.2
		10-Jul-2026
Page 5/7		

- | | |
|--|-----------------------------|
| · Self-reactive substances and mixtures | Not classified as hazardous |
| · Pyrophoric liquids | Not classified as hazardous |
| · Pyrophoric solids | Not classified as hazardous |
| · Self-heating substances and mixtures | Not classified as hazardous |
| · Substances and mixtures, which emit flammable gases in contact with water | Not classified as hazardous |
| · Oxidising liquids | Not classified as hazardous |
| · Oxidising solids | Not classified as hazardous |
| · Organic peroxides | Not classified as hazardous |
| · Corrosive to metals | Not classified as hazardous |
| · Desensitised explosives | Not classified as hazardous |

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used and stored according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **On the skin:** Based on available data, the classification criteria are not met.
- **On the eyes:** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **11.2 Information on other hazards**
- **Endocrine disrupting properties:**
None of the ingredients is listed.

SECTION 12: Ecological information

- **12.1 Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability :** No further relevant information available.
- **12.3 Bioaccumulative potential :** No further relevant information available.
- **12.4 Mobility in soil :** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus base Ref : SA16B	AD-ASQ-260
		Version 2.2
		10-Jul-2026
Page 6/7		

- **vPvB:** Not applicable.
- **12.6 Endocrine disrupting properties** For information on endocrine disrupting properties see section 11.
- **12.7 Other adverse effects :** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
 Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.
 Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation :** Must be specially treated according to official regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

SECTION 14: Transport information

- **14.1 UN number or ID number**
- **DOT, ADR, ADN, IMDG, IATA** none
- **14.2 UN proper shipping name**
- **ADR, ADN, IMDG, IATA** none
- **14.3 Transport hazard class(es)**
- **ADR, ADN, IMDG, IATA**
- **Class** none
- **14.4 Packing group**
- **ADR, IMDG, IATA** none
- **14.5 Environmental hazards:** Not applicable.
- **14.6 Special precautions for user** Not applicable.
- **14.7 Maritime transport in bulk according to IMO instruments** Not applicable.
- **UN "Model Regulation":** none

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II:**
 None of the ingredients is listed.
- **REGULATION (EU) 2019/1148**
- **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3)):**
 None of the ingredients is listed.
- **Annex II - REPORTABLE EXPLOSIVES PRECURSORS:**
 None of the ingredients is listed.

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus base Ref : SA16B	AD-ASQ-260
		Version 2.2
		10-Jul-2026
Page 7/7		

· **Regulation (EC) No 273/2004 on drug precursors:**

None of the ingredients is listed.

· **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors:**

None of the ingredients is listed.

· **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Classification in accordance with Delegated Regulation (EU) 2022/692 of 16 February 2022 adapting to technical and scientific progress GB CLP for the 18th time and amending table 3 of Annex VI part 3 of Regulation (EC) 1272/2008.

This Safety Data Sheet is in compliance with Regulation (EC) No 1907/2006, Article 31 as amended by Regulation (EU) 2020/878.

· **Relevant phrases**

H302 Harmful if swallowed.

H335 May cause respiratory irritation.

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity – Category 4

STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
Page 1/10		Version 1.1
		10-Jul-2026

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

- Trade name: **CHROMagar™ Salmonella Plus supplement**
- UFI: PX00-W06S-T00H-EUDV

1.2 Relevant identified uses of the substance or mixture and uses advised against

- Sector of Use** SU3 Industrial uses: Uses of substances as such or in preparations on industrial sites
- Application of the substance / the preparation** : Reagent for Laboratory Use

1.3 Details of the supplier of the safety data sheet

Manufacturer/Supplier:

CHROMagar
29 Avenue George Sand, 93210 La Plaine St-Denis
Tel: +33 1 45 48 05 05
Email : chromagar@chromagar.com

1.4 Emergency telephone number:

France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112
Belgium - Brüssel: +32 070/245 245
UK : Call NHS 999 in a medical emergency. If you urgently need medical help or advice but it's not a life-threatening situation, call 111 (24/7). The NHS 111 number is currently only available in certain areas. If you're outside of these areas, you should call NHS Direct on 0845 4647
Swiss Toxicological Information Centre (Zürich) : +41 44 251 51 51 / www.toxi.ch.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008 :



GHS05

Eye Dam. 1 H318 Causes serious eye damage.

Aquatic Chronic 3 H412 Harmful to aquatic life with long lasting effects.

ED ENV 1 EUH430 May cause endocrine disruption in the environment.

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008

The product is classified and labelled according to the C.L.P. regulation.

Hazard pictograms



GHS05

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
		Version 1.1
		10-Jul-2026

· **Signal word** Danger

· **Hazard-determining components of labelling:**

Surfactants

· **Hazard statements**

H318 Causes serious eye damage.

H412 Harmful to aquatic life with long lasting effects.

EUH430 May cause endocrine disruption in the environment.

· **Precautionary statements**

P273 Avoid release to the environment.

P280 Wear protective gloves / eye protection / face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER/doctor.

P405 Store locked up.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Labelling of packages where the contents do not exceed 125 ml**

· **Hazard pictograms**



GHS05

· **Signal word** Danger

· **Hazard-determining components of labelling:**

Surfactants

· **Hazard statements**

H318 Causes serious eye damage.

H412 Harmful to aquatic life with long lasting effects.

· **Precautionary statements**

· P273 Avoid release to the environment.

· P280 Wear eye protection / face protection.

· P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

· P310 Immediately call a POISON CENTER/doctor.

· **2.3 Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

· **Determination of endocrine-disrupting properties**

Nonionic surfactant, branched: List I


SECTION 3: Composition/information on ingredients

· **3.2 Mixtures**

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
		Version 1.1
		10-Jul-2026

· **Description:** Preparation: composed of the substances listed below.

· **Dangerous components:**

CAS:	Nonionic surfactant	10-25%
NLP:		

· **SVHC** Not applicable

· **Additional information:** For the wording of the listed risk phrases, refer to section 16.

SECTION 4: First aid measures

· **4.1 Description of first aid measures**

· **General information:** Position and transport stably in side position.

· **After inhalation:** Provide fresh air; consult a doctor in case of medical complaints.

· **After skin contact:**

In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.

If skin irritation continues, consult a doctor.

· **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.

· **After swallowing:** Do not induce vomiting; call for medical help immediately.

· **4.2 Most important symptoms and effects, both acute and delayed :**

No further relevant information is available.

· **4.3 Indication of any immediate medical attention and special treatment needed :**

No further relevant information is available.

SECTION 5: Firefighting measures

· **5.1 Suitable extinguishing agents:**

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

· **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.

· **5.3 Advice for firefighters**

· **Protective equipment:** No special measures required.

SECTION 6: Accidental release measures

· **6.1 Personal precautions, protective equipment and emergency procedures :**

No special measures required.

· **6.2 Environmental precautions:**

Inform respective authorities in case of seepage into water course or sewage system.

Do not allow the product to enter sewers/ surface water or ground water.

· **6.3 Methods and material for containment and cleaning up:**

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
Page 4/10		Version 1.1
		10-Jul-2026

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust). Dispose contaminated material as waste according to section 13.

Ensure adequate ventilation.

· **6.4 Reference to other sections :**

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

SECTION 7: Handling and storage

· **7.1 Precautions for safe handling :**

Ensure good ventilation/exhaustion at the workplace. Prevent formation of aerosols.

· **Information about fire - and explosion protection:** No special measures required.

· **7.2 Conditions for safe storage, including any incompatibilities**

· **Storage:**

· **Requirements to be met by storerooms and receptacles:** Store away from light and moisture.

· **Information about storage in one common storage facility:** Not required.

· **Further information about storage conditions:** Keep container tightly sealed.

· **Recommended storage temperature:** 15-30°C.

· **7.3 Specific end use(s) :** No further relevant information available.

SECTION 8: Exposure controls/personal protection

· **8.1 Control parameters**

· **Ingredients with limit values that require monitoring at the workplace:**

The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.

· **Additional information:**

The most current valid lists have been used as a basis for the production of this document.

· **8.2 Exposure controls**

· **Appropriate engineering controls** No further data; see section 7.

· **Individual protection measures, such as personal protective equipment**

· **General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

· **Respiratory protection:**



In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

Filter A/P2

· **Hand protection**



Protective gloves

Chemical resistant gloves EN 374

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
Page 5/10		Version 1.1
		10-Jul-2026

- **Gloves material** : Butyl rubber, BR
- **Eye/face protection**



Safety glasses

- **Body protection:**



Protective work clothing

SECTION 9: Physical and chemical properties

· 9.1 Information on basic physical and chemical properties

· General Information

· Physical state	Liquid
· Form:	Liquid
· Colour:	Colourless
· Odour:	Characteristic
· Odour threshold:	Not determined.
· Melting point/freezing point:	Not determined.
· Boiling point or initial boiling point and boiling range	Not determined.
· Flammability	Not determined.
· Lower and upper explosion limit	
· Lower:	Not determined.
· Upper:	Not determined.
· Flash point:	Not applicable.
· Auto-ignition temperature:	Not applicable.
· Decomposition temperature:	Not determined.
· pH	Not applicable.
· Viscosity:	
· Kinematic viscosity	Not applicable.
· Dynamic:	Not applicable.
· Solubility	
· water:	Fully miscible.
· Partition coefficient n-octanol/water (log value)	Not determined.
· Vapour pressure:	Not applicable.
· Density and/or relative density	
· Density at 20 °C:	1 +/-10% g/cm ³
· Relative density :	Not determined.
· Vapour density :	Not applicable.

· 9.2 Other information :

· Important information on protection of health and environment, and on safety.

- **Ignition temperature:** Product is not self igniting.
- **Explosive properties:** Product does not present an explosion hazard.

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
Page 6/10		Version 1.1
		10-Jul-2026

- **Solvent content:**
- **COV (CH) :** 0.00 %
- **Change in condition**
- **Evaporation rate :** Not applicable.
- **Information with regard to physical hazard classes**
 - **Explosives** Not classified as hazardous
 - **Flammable gases** Not classified as hazardous
 - **Aerosols** Not classified as hazardous
 - **Oxidising gases** Not classified as hazardous
 - **Gases under pressure** Not classified as hazardous
 - **Flammable liquids** Not classified as hazardous
 - **Flammable solids** Not classified as hazardous
 - **Self-reactive substances and mixtures** Not classified as hazardous
 - **Pyrophoric liquids** Not classified as hazardous
 - **Pyrophoric solids** Not classified as hazardous
 - **Self-heating substances and mixtures** Not classified as hazardous
 - **Substances and mixtures, which emit flammable gases in contact with water** Not classified as hazardous
 - **Oxidising liquids** Not classified as hazardous
 - **Oxidising solids** Not classified as hazardous
 - **Organic peroxides** Not classified as hazardous
 - **Corrosive to metals** Not classified as hazardous
 - **Desensitised explosives** Not classified as hazardous

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used and stored according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **Primary irritant effect:**
- **On the skin:**
Causes skin irritation.
- **On the eyes:** Causes serious eye damage.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
Page 7/10		Version 1.1
		10-Jul-2026

· **Aspiration hazard** Based on available data, the classification criteria are not met.

· **11.2 Information on other hazards**

· **Endocrine disrupting properties**

Nonionic surfactant, branched: List I

SECTION 12: Ecological information

· **12.1 Toxicity**

· **Aquatic toxicity:** No further relevant information available.

· **12.2 Persistence and degradability :** No further relevant information available.

· **12.3 Bioaccumulative potential :** No further relevant information available.

· **12.4 Mobility in soil :** No further relevant information available.

· **12.5 Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

· **12.6 Endocrine disrupting properties** For information on endocrine disrupting properties see section 11.

· **12.7 Other adverse effects :** No further relevant information available.

· **Remark:** Harmful to fish

· **Additional ecological information:**

· **General notes:**

Water hazard class 3 (German Regulation) (Self-assessment): extremely hazardous for water

Do not allow product to reach ground water, water course or sewage system, even in small quantities.

Danger to drinking water if even extremely small quantities leak into the ground.

Harmful to aquatic organisms

SECTION 13: Disposal considerations

· **13.1 Waste treatment methods**

· **Recommendation :** Must be specially treated according to official regulations.

· **European waste catalogue**

HP4 Irritant - skin irritation and eye damage

HP14 Ecotoxic

· **Uncleaned packaging:**

· **Recommendation:** Disposal must be made according to official regulations.

SECTION 14: Transport information

· **14.1 UN number or ID number**

· **DOT, ADR, IMDG, IATA**

UN3082

· **14.2 UN proper shipping name**

ADR

UN3082 ENVIRONMENTALLY HAZARDOUS
SUBSTANCE, LIQUID, N.O.S.

IMDG, IATA

ENVIRONMENTALLY HAZARDOUS SUBSTANCE,
LIQUID, N.O.S.

· **14.3 Transport hazard class(es)**

· **ADR**



	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
Page 8/10		Version 1.1
		10-Jul-2026

Class	9 (M6) Miscellaneous dangerous substances and articles.
· Label	9
· IMDG	
	
· Class	9 Miscellaneous dangerous substances and articles.
· Label	9
· IATA	
	
· Class	9 Miscellaneous dangerous substances and articles.
· Label	9
· 14.4 Packing group	
· ADR, IMDG, IATA	III
· 14.5 Environmental hazards:	
· Special marking (ADR):	Symbol (fish and tree)
· Special marking (IATA):	Symbol (fish and tree)
· 14.6 Special precautions for user	Warning: Miscellaneous dangerous substances and articles.
· Hazard identification number (Kemler code):	90
· EMS Number:	F-A,S-F
· Stowage Category	A
· 14.7 Maritime transport in bulk according to IMO instruments	Not applicable.
· Transport/Additional information:	
· ADR	
· Limited quantities (LQ)	5L
· Excepted quantities (EQ)	Code: E1
Maximum net quantity per inner packaging: 30 ml Maximum net quantity per outer packaging: 1000 ml	
· Transport category	3
· Tunnel restriction code	(-)
· IMDG	
· Limited quantities (LQ)	5L
· Excepted quantities (EQ)	Code: E1
Maximum net quantity per inner packaging: 30 ml	
Maximum net quantity per outer packaging: 1000 ml	
UN "Model Regulation":	UN 3082 ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. 9, III

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
		Version 1.1
		10-Jul-2026
Page 9/10		

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Poisons Act**
- **Regulated explosives precursors**
None of the ingredients is listed.
- **Regulated poisons**
None of the ingredients is listed.
- **Reportable explosives precursors**
None of the ingredients is listed.
- **Reportable poisons**
None of the ingredients is listed.
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **LIST OF SUBSTANCES SUBJECT TO AUTHORISATION (ANNEX XIV)**
Nonionic surfactant, branched: Sunset date: 2021-01-04
- **LIST OF SUBSTANCES SUBJECT TO AUTHORISATION (UK ANNEX XIV)**
Nonionic surfactant, branched: Sunset date: 2021-01-04
- **REGULATION (EC) No 1907/2006 ANNEX XVII** Conditions of restriction: 3
- **Regulation (EU) No 649/2012**
Nonionic surfactant, branched: Annex I Part 1 Annex I Part 2
- **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II**
None of the ingredients is listed.
- **REGULATION (EU) 2019/1148**
- **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3))**
None of the ingredients is listed.
- **Annex II - REPORTABLE EXPLOSIVES PRECURSORS**
None of the ingredients is listed.
- **Regulation (EC) No 273/2004 on drug precursors**
None of the ingredients is listed.
- **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors**
None of the ingredients is listed.
- **National regulations:**
- **Substances of very high concern (SVHC) according to UK REACH**
Nonionic surfactant, branched: List I
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any

	SAFETY DATA SHEET (according 1907/2006/CE, Art. 31) CHROMagar™ Salmonella Plus supplement Ref : SA16SN	AD-ASQ-401
		Version 1.1
		10-Jul-2026
Page 10/10		

specific product features and shall not establish a legally valid contractual relationship.

Classification in accordance with Delegated Regulation (EU) 2022/692 of 16 February 2022 adapting to technical and scientific progress C.L.P for the 18th time and amending table 3 of Annex VI part 3 of Regulation (EC) 1272/2008.

Classification in accordance with Delegated Regulation (EU) 2023/1434 of 25 April 2023 adapting to technical and scientific progress for the 19th time C.L.P and amending Annex VI part 1, section 1.1.3 of Regulation (EC) 1272/2008.

Classification in accordance with Delegated Regulation (EU) 2023/1435 of 2 May 2023 adapting to technical and scientific progress for the 20th time C.L.P and amending Annex VI part 3 of Regulation (EC) 1272/2008.

Classification in accordance with Delegated Regulation (EU) 2024/197 of 19 October 2023 adapting to technical and scientific progress for the 21st time C.L.P and amending Annex VI part 3 of Regulation (EC) 1272/2008.

This Safety Data Sheets is in compliance with Regulation (EC) No 1907/2006, Article 31 as amended by Regulation (EU) 2023/707.

· **Relevant phrases**

EUH430 May cause endocrine disruption in the environment. H302 Harmful if swallowed.

H302 Harmful if swallowed.

H318 Causes serious eye damage.

H332 Harmful if inhaled.

H400 Very toxic to aquatic life.

H411 Toxic to aquatic life with long lasting effects.

· **Classification according to Regulation (EC) No 1272/2008**

Skin corrosion/irritation Serious eye damage/irritation

Hazardous to the aquatic environment - long-term (chronic) aquatic hazard

Endocrine disruptor for the environment

· **Date of previous version:** 26.02.2025

· **Abbreviations and acronyms:**

The classification of the mixture is generally based on the calculation method using substance data according to Regulation (EC) No 1272/2008.

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

VOCV: Lenkungsabgabe auf flüchtigen organischen Verbindungen, Schweiz (Swiss Ordinance on volatile organic compounds)

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity – Category 4

Eye Dam. 1: Serious eye damage/eye irritation – Category 1

Aquatic Acute 1: Hazardous to the aquatic environment - acute aquatic hazard – Category

Aquatic Chronic 2: Hazardous to the aquatic environment - long-term aquatic hazard – Category 2

Aquatic Chronic 3: Hazardous to the aquatic environment - long-term aquatic hazard – Category 3