
	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens base</b> <b>Ref : PF652(B)</b>	AD-ASQ-305
		Version 2.0
		25-Aug-2023
Page 1/7		

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name:** CHROMagar™ C.perfringens base
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
- **Sector of Use** SU3 Industrial uses: Uses of substances as such or in preparations on industrial sites
- **Application of the substance / the preparation :** Reagent for Laboratory Use
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
 CHROMagar  
 4, place du 18 Juin 1940 - 75006 Paris  
 Tel: +33 1 45 48 05 05  
 Email : chromagar@chromagar.com
- **1.4 Emergency telephone number:**  
 France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59  
 Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
 Belgium - Brüssel: +32 070/245 245  
 UK : Call NHS 999 in a medical emergency. If you urgently need medical help or advice but it's not a life-threatening situation, call 111 (24/7). The NHS 111 number is currently only available in certain areas. If you're outside of these areas, you should call NHS Direct on 0845 4647  
 Swiss Toxicological Information Centre (Zürich) : +41 44 251 51 51 / www.toxi.ch.

## SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**  
 The product is not classified, according to the CLP regulation.
- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Not classified as hazardous
- **Hazard pictograms** Not classified as hazardous
- **Signal word** Not classified as hazardous
- **Hazard statements** Not classified as hazardous
- **Additional information:**  
 EUH210 Safety data sheet available on request.
- **2.3 Other hazards**  
 The mixture shall not contain substances listed in accordance with Article 59(1) of REACH as having endocrine disrupting properties, or be recognised as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 in a concentration equal to or greater than 0.1%.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens base</b> Ref : PF652(B)	AD-ASQ-305
		Version 2.0
Page 2/7		25-Aug-2023

### SECTION 3: Composition/information on ingredients

#### · 3.2 Mixtures

· **Description:** Preparation: composed of the substances listed below.

#### · Dangerous components:

CAS: 497-19-8                      sodium carbonate                      ≤2.5%

EINECS: 207-838-8               Eye Irrit. 2, H319

Index number: 011-005-00-2

Reg.nr.: 01-2119485498-19

CAS: 7757-79-1                      potassium nitrate                      ≤2.5%

EINECS: 231-818-8               Ox. Sol. 2, H272

RTECS: TT 3700000

· **SVHC** -

· **Additional information:** For the wording of the listed risk phrases, refer to section 16.

### SECTION 4: First aid measures

#### · 4.1 Description of first aid measures

#### · **General information:**

Position and transport stably in side position.

In case of irregular breathing or respiratory arrest provide artificial respiration.

· **After inhalation:** Provide fresh air; consult a doctor in case of medical complaints.

#### · **After skin contact:**

Generally the product does not irritate the skin.

In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.

If skin irritation continues, consult a doctor.

· **After eye contact:** Rinse opened eye for several minutes under running water.

· **After swallowing:** Do not induce vomiting; call for medical help immediately.

#### · **4.2 Most important symptoms and effects, both acute and delayed :**

No further relevant information available.

#### · **4.3 Indication of any immediate medical attention and special treatment needed :**

No further relevant information available.

### SECTION 5: Firefighting measures

#### · **5.1 Suitable extinguishing agents:**

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

· **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.

#### · **5.3 Advice for firefighters**

· **Protective equipment:** No special measures required.


### SECTION 6: Accidental release measures

#### · **6.1 Personal precautions, protective equipment and emergency procedures :**

No special measures required.

Avoid formation of dust.

#### · **6.2 Environmental precautions:**

 The Chromogenic Media Pioneer	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens base</b> Ref : PF652(B)	AD-ASQ-305
		Version 2.0
Page 3/7		25-Aug-2023

Do not allow product to reach sewage system or any water course.

Inform respective authorities in case of seepage into water course or sewage system.

Do not allow the product to enter sewers/ surface water or ground water.

- **6.3 Methods and material for containment and cleaning up:** Pick up mechanically.
- **6.4 Reference to other sections :**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protective equipment.  
See Section 13 for disposal information.

## SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :**  
Keep away from heat and direct sunlight.  
Thorough dedusting.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:**  
Store away from light and moisture.  
No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Store in cool, dry conditions in well sealed receptacles.
- **Recommended storage temperature:** 15-30°C
- **7.3 Specific end use(s) :** No further relevant information available.

## SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**
- **Ingredients with limit values that require monitoring at the workplace:**  
The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:**  
The most current valid lists have been used as a basis for the production of this document.
- **8.2 Exposure controls**
- **Appropriate engineering controls** No further data; see section 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**  
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:**



Filter A/P2

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.


- **Hand protection**



Protective gloves

Chemical resistant gloves EN 374

- **Gloves material :** Butyl rubber, BR

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens base</b> Ref : PF652(B)	AD-ASQ-305
		Version 2.0
		25-Aug-2023

· **Eye/face protection**



Safety glasses

· **Body protection:**



Protective work clothing

## SECTION 9: Physical and chemical properties

· **9.1 Information on basic physical and chemical properties**


· **General Information**

· <b>Physical state</b>	Solid
· <b>Form:</b>	Powder
· <b>Colour:</b>	Beige
· <b>Odour:</b>	Characteristic
· <b>Odour threshold:</b>	Not determined.
· <b>Melting point/freezing point:</b>	Undetermined.
· <b>Boiling point or initial boiling point and boiling range</b>	Undetermined.
· <b>Flammability</b>	Not determined.
· <b>Lower and upper explosion limit</b>	
· <b>Lower:</b>	Not determined.
· <b>Upper:</b>	Not determined.
· <b>Flash point:</b>	Not applicable.
· <b>Auto-ignition temperature:</b>	Not applicable.
· <b>Decomposition temperature:</b>	Not determined.
· <b>pH</b>	Not applicable.
· <b>Viscosity:</b>	
· <b>Kinematic viscosity</b>	Not applicable.
· <b>Kinematic:</b>	Not applicable.
· <b>Dynamic:</b>	Not applicable.
· <b>Solubility</b>	
· <b>water:</b>	Partly soluble.
· <b>Partition coefficient n-octanol/water (log value)</b>	Not determined.
· <b>Vapour pressure:</b>	Not applicable.
· <b>Density and/or relative density</b>	
· <b>Density at 20 °C:</b>	0.5+/-0.1 g/cm <sup>3</sup>
· <b>Relative density :</b>	Not determined.
· <b>Vapour density :</b>	Not applicable.
· <b>Particle characteristics</b>	
See section 3.	

· **9.2 Other information :**

· **Important information on protection of health and environment, and on safety.**

· <b>Ignition temperature:</b>	Product is not self igniting.
· <b>Explosive properties:</b>	Product does not present an explosion hazard.
· <b>Change in condition</b>	
· <b>Evaporation rate :</b>	Not applicable.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens base</b> <b>Ref : PF652(B)</b>	AD-ASQ-305
		Version 2.0
Page 5/7		25-Aug-2023

· **Information with regard to physical hazard classes**


· <b>Explosives</b>	Not classified as hazardous
· <b>Flammable gases</b>	Not classified as hazardous
· <b>Aerosols</b>	Not classified as hazardous
· <b>Oxidising gases</b>	Not classified as hazardous
· <b>Gases under pressure</b>	Not classified as hazardous
· <b>Flammable liquids</b>	Not classified as hazardous
· <b>Flammable solids</b>	Not classified as hazardous
· <b>Self-reactive substances and mixtures</b>	Not classified as hazardous
· <b>Pyrophoric liquids</b>	Not classified as hazardous
· <b>Pyrophoric solids</b>	Not classified as hazardous
· <b>Self-heating substances and mixtures</b>	Not classified as hazardous
· <b>Substances and mixtures, which emit flammable gases in contact with water</b>	Not classified as hazardous
· <b>Oxidising liquids</b>	Not classified as hazardous
· <b>Oxidising solids</b>	Not classified as hazardous
· <b>Organic peroxides</b>	Not classified as hazardous
· <b>Corrosive to metals</b>	Not classified as hazardous
· <b>Desensitised explosives</b>	Not classified as hazardous

## SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**  
No decomposition if used and stored according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

## SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **LD/LC50 values relevant for classification :**  
**CAS: 497-19-8 sodium carbonate**  
 Oral LD50 2,800 mg/kg (rat)  
 Dermal LD50 >2,000 mg/kg (rabbit)  
 Inhalative LC50/4 h 2.3 mg/l (rat)  
**CAS: 7757-79-1 potassium nitrate**  
 Oral LD50 3,750 mg/kg (rat)
- **On the skin:** Based on available data, the classification criteria are not met.
- **On the eyes:** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens base</b> <b>Ref : PF652(B)</b>	AD-ASQ-305
		Version 2.0
Page 6/7		25-Aug-2023

· **STOT-repeated exposure** Based on available data, the classification criteria are not met.

· **Aspiration hazard** Based on available data, the classification criteria are not met.

· **11.2 Information on other hazards**

· **Endocrine disrupting properties :**

None of the ingredients is listed.

## SECTION 12: Ecological information

· **12.1 Toxicity**

· **Aquatic toxicity :**

**CAS: 497-19-8 sodium carbonate**

EC50 48h 220 mg/l (daphnia)

EC50 96h 300 mg/l (fish)

· **12.2 Persistence and degradability :** No further relevant information available.

· **12.3 Bioaccumulative potential :** No further relevant information available.

· **12.4 Mobility in soil :** No further relevant information available.

· **12.5 Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

· **12.6 Endocrine disrupting properties** For information on endocrine disrupting properties see section 11.

· **12.7 Other adverse effects :** No further relevant information available.

· **Additional ecological information:**

· **General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

## SECTION 13: Disposal considerations

· **13.1 Waste treatment methods**

· **Recommendation :** Must be specially treated according to official regulations.

· **Uncleaned packaging:**

· **Recommendation:** Disposal must be made according to official regulations.

## SECTION 14: Transport information

· **14.1 UN number or ID number**

· **DOT, ADR, ADN, IMDG, IATA** none

· **14.2 UN proper shipping name**

· **ADR, ADN, IMDG, IATA** none

· **14.3 Transport hazard class(es)**

· **ADR, ADN, IMDG, IATA**

· **Class** none

· **14.4 Packing group**


· **ADR, IMDG, IATA** none

· **14.5 Environmental hazards:** Not applicable.

· **14.6 Special precautions for user** Not applicable.

· **14.7 Maritime transport in bulk according to IMO instruments** Not applicable.

· **UN "Model Regulation":** none

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens base</b> Ref : PF652(B)	AD-ASQ-305
		Version 2.0
		25-Aug-2023
Page 7/7		

## SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II :**  
None of the ingredients is listed.
- **REGULATION (EU) 2019/1148**
- **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3)) :**  
None of the ingredients is listed.
- **Annex II - REPORTABLE EXPLOSIVES PRECURSORS :**  
CAS: 7757-79-1 potassium nitrate
- **Regulation (EC) No 273/2004 on drug precursors :**  
None of the ingredients is listed.
- **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors**  
None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Classification in accordance with Delegated Regulation (EU) 2022/692 of 16 February 2022 adapting to technical and scientific progress GB CLP for the 18th time and amending table 3 of Annex VI part 3 of Regulation (EC) 1272/2008.

This Safety Data Sheet is in compliance with Regulation (EC) No 1907/2006, Article 31 as amended by Regulation (EU) 2020/878.

- **Relevant phrases**  
H272 May intensify fire; oxidiser.  
H319 Causes serious eye irritation.
- **Abbreviations and acronyms:**  
ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
DOT: US Department of Transportation  
IATA: International Air Transport Association  
GHS: Globally Harmonised System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service (division of the American Chemical Society)  
LC50: Lethal concentration, 50 percent  
LD50: Lethal dose, 50 percent  
PBT: Persistent, Bioaccumulative and Toxic  
SVHC: Substances of Very High Concern  
vPvB: very Persistent and very Bioaccumulative  
Ox. Sol. 2: Oxidizing solids – Category 2  
Eye Irrit. 2: Serious eye damage/eye irritation – Category 2



	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S1</b> <b>Ref : PF652(S1)</b></p>	AD-ASQ-306
Page 1/8		Version 2.0
		25-Aug-2023

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name:** CHROMagar™ C.perfringens supplement S1
- **UFI:** YE00-V02D-Y00J-FSWG
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
- **Sector of Use** SU3 Industrial uses: Uses of substances as such or in preparations on industrial sites
- **Application of the substance / the preparation :** Reagent for Laboratory Use
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
CHROMagar  
4, place du 18 Juin 1940 - 75006 Paris  
Tel: +33 1 45 48 05 05  
Email : chromagar@chromagar.com
- **1.4 Emergency telephone number:**  
France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59  
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
Belgium - Brüssel: +32 070/245 245  
UK : Call NHS 999 in a medical emergency. If you urgently need medical help or advice but it's not a life-threatening situation, call 111 (24/7). The NHS 111 number is currently only available in certain areas. If you're outside of these areas, you should call NHS Direct on 0845 4647  
Swiss Toxicological Information Centre (Zürich) : +41 44 251 51 51 / www.toxi.ch.

## SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**



GHS07

Eye Irrit. 2 H319 Causes serious eye irritation.

- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008**  
The product is classified and labelled according to the CLP regulation.
- **Hazard pictograms**



GHS07

- **Signal word** Warning
- **Hazard statements**  
H319 Causes serious eye irritation.
- **Precautionary statements**  
P264 Wash thoroughly after handling.  
P280 Wear eye protection / face protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.



	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S1</b> <b>Ref : PF652(S1)</b></p>	AD-ASQ-306
Page 2/8		Version 2.0
		25-Aug-2023

- Labelling of packages where the contents do not exceed 125 ml
- Hazard pictograms



GHS07

- **Signal word** Warning
- **Hazard statements** Not classified as hazardous
- **2.3 Other hazards**

The mixture shall not contain substances listed in accordance with Article 59(1) of REACH as having endocrine disrupting properties, or be recognised as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 in a concentration equal to or greater than 0.1%.

- **Results of PBT and vPvB assessment**

- **PBT:** Not applicable.
- **vPvB:** Not applicable.

### SECTION 3: Composition/information on ingredients

- **3.2 Mixtures**

- **Description:** Preparation: composed of the substances listed below.

- **Dangerous components :**

CAS: 497-19-8                      sodium carbonate    25-<50%

EINECS: 207-838-8               Eye Irrit. 2, H319

Index number: 011-005-00-2

Reg.nr.: 01-2119485498-19

- **SVHC -**

- **Additional information:** For the wording of the listed risk phrases, refer to section 16.

### SECTION 4: First aid measures

- **4.1 Description of first aid measures**

- **General information:**

Position and transport stably in side position.

In case of irregular breathing or respiratory arrest provide artificial respiration.

- **After inhalation:** Provide fresh air; consult a doctor in case of medical complaints.

- **After skin contact:**

Generally the product does not irritate the skin.

In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.

If skin irritation continues, consult a doctor.

- **After eye contact:** Rinse opened eye for several minutes under running water.

- **After swallowing:** Do not induce vomiting; call for medical help immediately.

- **4.2 Most important symptoms and effects, both acute and delayed :**

No further relevant information available.

- **4.3 Indication of any immediate medical attention and special treatment needed :**

No further relevant information available.

### SECTION 5: Firefighting measures

	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S1</b> <b>Ref : PF652(S1)</b></p>	AD-ASQ-306
Page 3/8		Version 2.0
		25-Aug-2023

- **5.1 Suitable extinguishing agents:**  
CO<sub>2</sub>, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.
- **5.3 Advice for firefighters**
- **Protective equipment:** No special measures required.

## SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures :**  
No special measures required.  
Avoid formation of dust.
- **6.2 Environmental precautions:**  
Do not allow product to reach sewage system or any water course.  
Inform respective authorities in case of seepage into water course or sewage system.  
Do not allow the product to enter sewers/ surface water or ground water.
- **6.3 Methods and material for containment and cleaning up:** Pick up mechanically.
- **6.4 Reference to other sections :**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protective equipment.  
See Section 13 for disposal information.

## SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :**  
Keep away from heat and direct sunlight.  
Thorough dedusting.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:**  
Store away from light and moisture.  
No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Store in cool, dry conditions in well sealed receptacles.
- **Recommended storage temperature:** 2-8 °C
- **7.3 Specific end use(s) :** No further relevant information available.

## SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**
- **Ingredients with limit values that require monitoring at the workplace :**  
**CAS: 1185-57-5 ammonium iron(III) citrate (50-100%)**  
WEL (Great Britain) Short-term value: 2 mg/m<sup>3</sup>  
Long-term value: 1 mg/m<sup>3</sup>  
as Fe  
MAK (Switzerland) Long-term value: 1e mg/m<sup>3</sup>  
als Fe berechnet
- **Additional information:**  
The most current valid lists have been used as a basis for the production of this document.
- **8.2 Exposure controls**
- **Appropriate engineering controls** No further data; see section 7.

	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S1</b> <b>Ref : PF652(S1)</b></p>	AD-ASQ-306
Page 4/8		Version 2.0
		25-Aug-2023

- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**  
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:**



In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

Filter A/P2

- **Hand protection**



Protective gloves

Chemical resistant gloves EN 374

- **Gloves material :** Butyl rubber, BR
- **Eye/face protection**



Safety glasses

- **Body protection:**



Protective work clothing

## SECTION 9: Physical and chemical properties

- **9.1 Information on basic physical and chemical properties**
- **General Information**
- **Physical state** Solid
- **Form:** Powder
- **Colour:** Brown
- **Odour:** Characteristic
- **Odour threshold:** Not determined.
- **Melting point/freezing point:** Undetermined.
- **Boiling point or initial boiling point and boiling range** Undetermined.
- **Flammability** Not determined.
- **Lower and upper explosion limit**
- **Lower:** Not determined.
- **Upper:** Not determined.
- **Flash point:** Not applicable.
- **Auto-ignition temperature:** Not applicable.
- **Decomposition temperature:** Not determined.
- **pH** Not applicable.
- **Viscosity:**
- **Kinematic viscosity** Not applicable.
- **Kinematic:** Not applicable.
- **Dynamic:** Not applicable.
- **Solubility**
- **water:** Soluble.

	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S1</b> <b>Ref : PF652(S1)</b></p>	AD-ASQ-306
Page 5/8		Version 2.0
		25-Aug-2023

- **Partition coefficient n-octanol/water (log value)** Not determined.
- **Vapour pressure:** Not applicable.
- **Density and/or relative density**
- **Density at 20 °C:** 0.5+/-0.1 g/cm<sup>3</sup>
- **Relative density :** Not determined.
- **Vapour density :** Not applicable.
- **Particle characteristics**  
See section 3.
- **9.2 Other information :**
- **Important information on protection of health and environment, and on safety.**
- **Ignition temperature:** Product is not self igniting.
- **Explosive properties:** Product does not present an explosion hazard.
- **Change in condition**
- **Evaporation rate :** Not applicable.
- **Information with regard to physical hazard classes**
- **Explosives** Not classified as hazardous
- **Flammable gases** Not classified as hazardous
- **Aerosols** Not classified as hazardous
- **Oxidising gases** Not classified as hazardous
- **Gases under pressure** Not classified as hazardous
- **Flammable liquids** Not classified as hazardous
- **Flammable solids** Not classified as hazardous
- **Self-reactive substances and mixtures** Not classified as hazardous
- **Pyrophoric liquids** Not classified as hazardous
- **Pyrophoric solids** Not classified as hazardous
- **Self-heating substances and mixtures** Not classified as hazardous
- **Substances and mixtures, which emit flammable gases in contact with water** Not classified as hazardous
- **Oxidising liquids** Not classified as hazardous
- **Oxidising solids** Not classified as hazardous
- **Organic peroxides** Not classified as hazardous
- **Corrosive to metals** Not classified as hazardous
- **Desensitised explosives** Not classified as hazardous

## SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**  
No decomposition if used and stored according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

 The Chromogenic Media Pioneer	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S1</b> <b>Ref : PF652(S1)</b>	AD-ASQ-306
Page 6/8		Version 2.0
		25-Aug-2023

## SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **LD/LC50 values relevant for classification :**  
**CAS: 497-19-8 sodium carbonate**  
Oral LD50 2,800 mg/kg (rat)  
Dermal LD50 >2,000 mg/kg (rabbit)  
Inhalative LC50/4 h 2.3 mg/l (rat)
- **On the skin:** Based on available data, the classification criteria are not met.
- **On the eyes:** Causes serious eye irritation.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **11.2 Information on other hazards**
- **Endocrine disrupting properties :**  
None of the ingredients is listed.

## SECTION 12: Ecological information

- **12.1 Toxicity**
- **Aquatic toxicity :**  
**CAS: 497-19-8 sodium carbonate**  
EC50 48h 220 mg/l (daphnia)  
EC50 96h 300 mg/l (fish)
- **12.2 Persistence and degradability :** No further relevant information available.
- **12.3 Bioaccumulative potential :** No further relevant information available.
- **12.4 Mobility in soil :** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Endocrine disrupting properties** For information on endocrine disrupting properties see section 11.
- **12.7 Other adverse effects :** No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

## SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation :** Must be specially treated according to official regulations.
- **European waste catalogue :**  
HP4 Irritant - skin irritation and eye damage

 The Chromogenic Media Pioneer	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S1</b> <b>Ref : PF652(S1)</b>	AD-ASQ-306
Page 7/8		Version 2.0
		25-Aug-2023

- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

## SECTION 14: Transport information

- **14.1 UN number or ID number**
- **DOT, ADR, ADN, IMDG, IATA** none
- **14.2 UN proper shipping name**
- **ADR, ADN, IMDG, IATA** none
- **14.3 Transport hazard class(es)**
- **ADR, ADN, IMDG, IATA**
- **Class** none
- **14.4 Packing group**
- **ADR, IMDG, IATA** none
- **14.5 Environmental hazards:** Not applicable.
- **14.6 Special precautions for user** Not applicable.
- **14.7 Maritime transport in bulk according to IMO instruments** Not applicable.
- **UN "Model Regulation":** none

## SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II :**  
None of the ingredients is listed.
- **REGULATION (EU) 2019/1148**
- **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3)) :**  
None of the ingredients is listed.
- **Annex II - REPORTABLE EXPLOSIVES PRECURSORS :**  
None of the ingredients is listed.
- **Regulation (EC) No 273/2004 on drug precursors :**  
None of the ingredients is listed.
- **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors :**  
None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Classification in accordance with Delegated Regulation (EU) 2022/692 of 16 February 2022 adapting to technical and scientific progress GB CLP for the 18th time and amending table 3 of Annex VI part 3 of Regulation (EC) 1272/2008.

This Safety Data Sheets is in compliance with Regulation (EC) No 1907/2006, Article 31 as amended by Regulation (EU) 2020/878.

	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S1</b> <b>Ref : PF652(S1)</b></p>	AD-ASQ-306
Page 8/8		Version 2.0
		25-Aug-2023

· **Relevant phrases**

H319 Causes serious eye irritation.

· **Classification according to Regulation (EC) No 1272/2008 :**

Serious eye damage/irritation The classification of the mixture is generally based on the calculation method using substance data according to Regulation (EC) No 1272/2008.

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Eye Irrit. 2: Serious eye damage/eye irritation – Category 2



	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S2</b> <b>Ref : PF652(S2)</b>	AD-ASQ-307
		Version 2.0
		25-Aug-2023
Page 1/7		

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name:** CHROMagar™ C.perfringens supplement S2
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
- **Sector of Use** SU3 Industrial uses: Uses of substances as such or in preparations on industrial sites
- **Application of the substance / the preparation :** Reagent for Laboratory Use
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
 CHROMagar  
 4, place du 18 Juin 1940 - 75006 Paris  
 Tel: +33 1 45 48 05 05  
 Email : chromagar@chromagar.com
- **1.4 Emergency telephone number:**  
 France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59  
 Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
 Belgium - Brüssel: +32 070/245 245  
 UK : Call NHS 999 in a medical emergency. If you urgently need medical help or advice but it's not a life-threatening situation, call 111 (24/7). The NHS 111 number is currently only available in certain areas. If you're outside of these areas, you should call NHS Direct on 0845 4647  
 Swiss Toxicological Information Centre (Zürich) : +41 44 251 51 51 / www.toxi.ch.

## SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**  
 The product is not classified, according to the CLP regulation.
- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Not classified as hazardous
- **Hazard pictograms** Not classified as hazardous
- **Signal word** Not classified as hazardous
- **Hazard statements** Not classified as hazardous
- **Additional information:**  
 EUH210 Safety data sheet available on request.
- **2.3 Other hazards**  
 The mixture shall not contain substances listed in accordance with Article 59(1) of REACH as having endocrine disrupting properties, or be recognised as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 in a concentration equal to or greater than 0.1%.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S2</b> <b>Ref : PF652(S2)</b>	AD-ASQ-307
		Version 2.0
		25-Aug-2023
Page 2/7		

### SECTION 3: Composition/information on ingredients

#### · 3.2 Mixtures

· **Description:** Preparation: composed of the substances listed below.

#### · Dangerous components :

CAS: 1405-20-5 Polymyxin B, sulfate (salt) 10-25%

⚠ Acute Tox. 4, H302

· **SVHC -**

· **Additional information:** For the wording of the listed risk phrases, refer to section 16.

### SECTION 4: First aid measures

#### · 4.1 Description of first aid measures

##### · **General information:**

Position and transport stably in side position.

In case of irregular breathing or respiratory arrest provide artificial respiration.

· **After inhalation:** Provide fresh air; consult a doctor in case of medical complaints.

##### · **After skin contact:**

Generally the product does not irritate the skin.

In the event of irritation, to withdraw soiled clothing, to wash the parts reached with soapy water, to rinse with large water.

If skin irritation continues, consult a doctor.

· **After eye contact:** Rinse opened eye for several minutes under running water.

· **After swallowing:** Do not induce vomiting; call for medical help immediately.

#### · **4.2 Most important symptoms and effects, both acute and delayed :**

No further relevant information available.

#### · **4.3 Indication of any immediate medical attention and special treatment needed :**

No further relevant information available.

### SECTION 5: Firefighting measures

#### · **5.1 Suitable extinguishing agents:**

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

· **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.

#### · **5.3 Advice for firefighters**

· **Protective equipment:** No special measures required.

### SECTION 6: Accidental release measures

#### · **6.1 Personal precautions, protective equipment and emergency procedures :**

No special measures required.

Avoid formation of dust.

#### · **6.2 Environmental precautions:**

Do not allow product to reach sewage system or any water course.

Inform respective authorities in case of seepage into water course or sewage system.

Do not allow the product to enter sewers/ surface water or ground water.

· **6.3 Methods and material for containment and cleaning up:** Pick up mechanically.

#### · **6.4 Reference to other sections :**

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S2</b> Ref : PF652(S2)</p>	AD-ASQ-307
		Version 2.0
Page 3/7		25-Aug-2023

## SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :**  
Keep away from heat and direct sunlight.  
Thorough dedusting.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:**  
Store away from light and moisture.  
No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Store in cool, dry conditions in well sealed receptacles.
- **Recommended storage temperature:** 2-8°C
- **7.3 Specific end use(s) :** No further relevant information available.

## SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**
- **Ingredients with limit values that require monitoring at the workplace:**  
The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:**  
The most current valid lists have been used as a basis for the production of this document.
- **8.2 Exposure controls**
- **Appropriate engineering controls** No further data; see section 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**  
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:**



In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

Filter A/P2

- **Hand protection**



Protective gloves

Chemical resistant gloves EN 374

- **Gloves material :** Butyl rubber, BR
- **Eye/face protection**



Safety glasses

	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S2</b> <b>Ref : PF652(S2)</b></p>	AD-ASQ-307
Page 4/7		Version 2.0
		25-Aug-2023

· **Body protection:**



Protective work clothing

## SECTION 9: Physical and chemical properties

· **9.1 Information on basic physical and chemical properties**

· **General Information**

· <b>Physical state</b>	Solid
· <b>Form:</b>	Powder
· <b>Colour:</b>	White
· <b>Odour:</b>	Characteristic
· <b>Odour threshold:</b>	Not determined.
· <b>Melting point/freezing point:</b>	Undetermined.
· <b>Boiling point or initial boiling point and boiling range</b>	Undetermined.
· <b>Flammability</b>	Not determined.
· <b>Lower and upper explosion limit</b>	
· <b>Lower:</b>	Not determined.
· <b>Upper:</b>	Not determined.
· <b>Flash point:</b>	Not applicable.
· <b>Auto-ignition temperature:</b>	Not applicable.
· <b>Decomposition temperature:</b>	Not determined.
· <b>pH</b>	Not applicable.
· <b>Viscosity:</b>	
· <b>Kinematic viscosity</b>	Not applicable.
· <b>Kinematic:</b>	Not applicable.
· <b>Dynamic:</b>	Not applicable.
· <b>Solubility</b>	
· <b>water:</b>	Soluble.
· <b>Partition coefficient n-octanol/water (log value)</b>	Not determined.
· <b>Vapour pressure:</b>	Not applicable.
· <b>Density and/or relative density</b>	
· <b>Density at 20 °C:</b>	0.5+/-0.1 g/cm <sup>3</sup>
· <b>Relative density :</b>	Not determined.
· <b>Vapour density :</b>	Not applicable.
· <b>Particle characteristics</b>	
See section 3.	

· **9.2 Other information :**

· **Important information on protection of health and environment, and on safety.**

· <b>Ignition temperature:</b>	Product is not self igniting.
· <b>Explosive properties:</b>	Product does not present an explosion hazard.
· <b>Change in condition</b>	
· <b>Evaporation rate :</b>	Not applicable.

· **Information with regard to physical hazard classes**

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S2</b> <b>Ref : PF652(S2)</b>	AD-ASQ-307
		Version 2.0
		25-Aug-2023

· <b>Explosives</b>	Not classified as hazardous
· <b>Flammable gases</b>	Not classified as hazardous
· <b>Aerosols</b>	Not classified as hazardous
· <b>Oxidising gases</b>	Not classified as hazardous
· <b>Gases under pressure</b>	Not classified as hazardous
· <b>Flammable liquids</b>	Not classified as hazardous
· <b>Flammable solids</b>	Not classified as hazardous
· <b>Self-reactive substances and mixtures</b>	Not classified as hazardous
· <b>Pyrophoric liquids</b>	Not classified as hazardous
· <b>Pyrophoric solids</b>	Not classified as hazardous
· <b>Self-heating substances and mixtures</b>	Not classified as hazardous
· <b>Substances and mixtures, which emit flammable gases in contact with water</b>	Not classified as hazardous
· <b>Oxidising liquids</b>	Not classified as hazardous
· <b>Oxidising solids</b>	Not classified as hazardous
· <b>Organic peroxides</b>	Not classified as hazardous
· <b>Corrosive to metals</b>	Not classified as hazardous
· <b>Desensitised explosives</b>	Not classified as hazardous

## SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**  
No decomposition if used and stored according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

## SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **On the skin:** Based on available data, the classification criteria are not met.
- **On the eyes:** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **11.2 Information on other hazards**
- **Endocrine disrupting properties :**  
one of the ingredients is listed.

## SECTION 12: Ecological information

- **12.1 Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability :** No further relevant information available.

 The Chromogenic Media Pioneer	<p align="center"><b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S2</b> <b>Ref : PF652(S2)</b></p>	AD-ASQ-307
Page 6/7		Version 2.0
		25-Aug-2023

- **12.3 Bioaccumulative potential** : No further relevant information available.
- **12.4 Mobility in soil** : No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
- **PBT**: Not applicable.
- **vPvB**: Not applicable.
- **12.6 Endocrine disrupting properties** For information on endocrine disrupting properties see section 11.
- **12.7 Other adverse effects** : No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

### SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation** : Must be specially treated according to official regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

### SECTION 14: Transport information

- **14.1 UN number or ID number**
- **DOT, ADR, ADN, IMDG, IATA** none
- **14.2 UN proper shipping name**
- **ADR, ADN, IMDG, IATA** none
- **14.3 Transport hazard class(es)**
- **ADR, ADN, IMDG, IATA**
- **Class** none
- **14.4 Packing group**
- **ADR, IMDG, IATA** none
- **14.5 Environmental hazards:** Not applicable.
- **14.6 Special precautions for user** Not applicable.
- **14.7 Maritime transport in bulk according to IMO instruments** Not applicable.
- **UN "Model Regulation":** none

### SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II :**  
None of the ingredients is listed.
- **REGULATION (EU) 2019/1148**
- **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3)) :**  
None of the ingredients is listed.

	<b>SAFETY DATA SHEET</b> (according 1907/2006/CE, Art. 31) <b>CHROMagar™ C.perfringens supplement S2</b> <b>Ref : PF652(S2)</b>	AD-ASQ-307
		Version 2.0
		25-Aug-2023
Page 7/7		

· **Annex II - REPORTABLE EXPLOSIVES PRECURSORS :**

None of the ingredients is listed.

· **Regulation (EC) No 273/2004 on drug precursors :**

None of the ingredients is listed.

· **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors :**

None of the ingredients is listed.

· **15.2 Chemical safety assessment :** A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Classification in accordance with Delegated Regulation (EU) 2022/692 of 16 February 2022 adapting to technical and scientific progress GB CLP for the 18th time and amending table 3 of Annex VI part 3 of Regulation (EC) 1272/2008.

This Safety Data Sheets is in compliance with Regulation (EC) No 1907/2006, Article 31 as amended by Regulation (EU) 2020/878.

· **Relevant phrases**

H302 Harmful if swallowed.

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity – Category 4