



Usefulness of Chromogenic Media with Fluconazole Supplementation for Presumptive Identification of *Candida auris*

Alba Ruiz-Gaitán ^{1,2,*,†}[®], Ignacio Sigona-Giangreco ^{2,†}[®], José Manuel Pérez-Royo ²[®], Victor Garcia-Bustos ^{2,3,*}, Marta García-Hita ²[®], Eulogio Valentín-Gómez ⁴[®], Salvador Giner Almaraz ¹, Piet W. J. de Groot ⁵[®] and Javier Pemán ^{1,2}

- ¹ Department of Clinical Microbiology, La Fe University and Polytechnic Hospital, 46026 Valencia, Spain
- ² Severe Infection Research Group, Medical Research Institute La Fe, 46026 Valencia, Spain
- ³ Department of Internal Medicine, La Fe University and Polytechnic Hospital, 46026 Valencia, Spain
- ⁴ Department of Microbiology and Ecology, University of Valencia, 46010 Valencia, Spain
- ⁵ Regional Center for Biomedical Research, Castilla-La Mancha Science & Technology Park, University of Castilla-La Mancha, 02008 Albacete, Spain
- * Correspondence: albacruiz@gmail.com or ruiz_alb@gva.es (A.R.-G.); victorgarciabustos@gmail.com or garcia_vicbus@gva.es (V.G.-B.); Tel.: +34-961246690 (A.R.-G.)
- + These authors contributed equally to this work.

Abstract: Introduction: Candida auris is a major threat to public health. Rapid detection is essential for early treatment and transmission control. The use of chromogenic media allows the presumptive identification of this new species. The aim of this study is to describe the morphological characteristics of C. auris colonies on three commercial chromogenic media. Methods: Nineteen C. auris isolates from different countries/clades and 18 isolates of other species were cultivated in CHROMagarTM Candida Plus, HiCrome[™] Candida, CHROMagar-Candida, and fluconazole-supplemented (32 mg/L) CHROMagar-Candida media. *Results*: On CHROMagarTM Candida Plus and HiCromeTM Candida, C. auris isolates from Colombia, Venezuela, India, Korea, and Japan displayed blue-shaded colonies, while isolates from Spain and Germany exhibited light pink shades with a bluish halo. All isolates showed white to pink colonies on CHROMagar-Candida. On CHROMagar Candida supplemented with fluconazole, whilst C. auris, C. glabrata, or C. krusei showed a similar pink color at 48 h incubation, phenotypic differentiation was possible by the rough, paraffin-like texture or the intense purple color acquired by C. krusei and C. glabrata, respectively. Moreover, in this medium, the presence of C. auris in combination with other species of similar color was not limiting for its early identification, due to this medium selecting only strains resistant to this antifungal. Conclusions: The use of chromogenic media such as CHROMagarTM Candida Plus facilitates a presumptive identification of *C. auris*. However, this identification can be difficult in the presence of mixed cultures. In these cases, the use of CHROMagarTM Candida medium with 32 mg/L fluconazole offers better performance for the identification of C. auris by inhibiting fluconazole-susceptible strains and selecting rare or high fluconazole MIC (>32 mg/L) isolates.

Keywords: *Candida auris*; chromogenic media; CHROMagar[™] Candida Plus; HiCromeTM Candida; fluconazole supplementation; yeast identification

1. Introduction

Since its first description in 2009 [1], *Candida auris* has been responsible for nosocomial outbreaks in several countries around the world [2]. Distinguishing features of this species include its multi-resistance to antifungal drugs, its potential to cause deep-seated infections with elevated mortality rates [3], its high transmissibility and persistence in the hospital environment [3,4], and the difficulty of identification by conventional biochemical and



Citation: Ruiz-Gaitán, A.; Sigona-Giangreco, I.; Pérez-Royo, J.M.; Garcia-Bustos, V.; García-Hita, M.; Valentín-Gómez, E.; Almaraz, S.G.; de Groot, P.W.J.; Pemán, J. Usefulness of Chromogenic Media with Fluconazole Supplementation for Presumptive Identification of *Candida auris. Diagnostics* **2023**, *13*, 231. https://doi.org/10.3390/ diagnostics13020231

Academic Editors: Andor W.J.M. Glaudemans, Mike M. Sathekge and Alfred O. Ankrah

Received: 31 October 2022 Revised: 23 December 2022 Accepted: 5 January 2023 Published: 8 January 2023



Copyright: © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). phenotypic techniques. These characteristics have led the CDC to classify *C. auris* as a global public health threat [5,6].

Phenotypically, *C. auris* is difficult to distinguish from other *Candida* species. On Sabouraud agar, *C. auris* forms white to cream colonies, and on chromogenic media such as CHROMagar Candida, it forms light pink colonies [7,8].

Available commercial identification systems generally misidentify *C. auris* as many other species with greater or lesser phylogenetic closeness, such as *Candida haemulonii*, *Rhodotorula glutinis*, *Saccharomyces cerevisiae*, or, less frequently, *Candida famata*, *Candida dobushaemulonii*, *Candida sake*, *Candida lusitaniae*, *Candida albicans*, *Candida guilliermondii*, or *Candida parapsilosis* [5].

Since the characterization of this yeast as an emerging global threat to public health and the acknowledgment of the difficulties for its identification using conventional biochemical techniques, such as API ID20C (bioMérieux) and AuxaColor™ 2 (BioRad-Laboratories, Marnes-la-Coquette, France), efforts have been driven to develop rapid detection methods with advanced techniques capable of avoiding the limitations of phenotypic tests, such as mass spectrometry (MALDI-TOF MS), or molecular techniques such as internal transcribed spacer sequencing (ITS). Despite their accuracy, these techniques are expensive and require qualified personnel, hampering their use at non-specialized centers or especially resourcelimited settings. This fact limits not only the early initiation of appropriate treatment, but also the early detection of cases and contacts [9,10]. Furthermore, all suspected cases of *C. auris* colonization or infection should be immediately isolated owing to its high transmissibility in the hospital environment and contact cases must be promptly identified and screened for *C. auris*. In this matter, it is essential to avoid the delay in the microbiological identification of C. auris, limit or control the generation of hospital outbreaks, and diminish patient morbimortality as well as the healthcare-associated burden [11]. The still-present lack of knowledge on easily reproducible, low-cost but sensitive techniques for prompt C. auris identification widens the gap in equity between developing countries [12].

Most microbiology laboratories perform presumptive identification of *Candida* spp. using chromogenic media [13,14]. Previous reports on CHROMagarTM *Candida* medium have shown that it is not useful for differentially identifying *C. auris* and *C. parapsilosis* colonies [15]. Consequently, other selective media such as CHROMagarTM Candida Plus medium (CHROMagar, Paris, France) and HiCromeTM Candida have been developed to further discriminate *C. auris* colonies from other species [16–18]. Despite recent progress in this area, microbiological data with these modified media for presumptive identification of *C. auris* are still limited to studies with a relatively low number of species and strains tested, discrepancies between morphotypes, or lack of information on mixed cultures.

Therefore, the aim of this study was to describe the phenotypic characteristic of *C. auris* colonies on four commercially available chromogenic media, as well as to evaluate their performance for presumptive identification of *C. auris* compared to other non-*C. auris* species and mixed cultures, including a modified fluconazole-supplemented CHROMagarTM Candida.

2. Materials and Methods

A total of 35 suspensions of pathogenic yeast species were inoculated in CHROMagarTM Candida Plus (CHROMagar, Paris, France), HiCromeTM Candida (HIMEDIA LBS, Mumbai, India), CHROMagarTM- Candida (Biolife, Milan, Italy), and CHROMagar Candida supplemented with 32 mg/L of fluconazole (CHROMagar Candida plus fluconazole) (MAIM, Barcelona, Spain). Of these, 19 *C. auris* strains from different countries were selected to represent the four main global clonal lineages. Thirteen of these strains were from the Westerdijk Fungal Biodiversity Institute (CBS-KNAW): CBS-12766, CBS-12882, CBS-12766 (clade I; India); CBS-12372, CBS-12373 (clade II; Korea); CBS-10913, JM-15448 (clade II: Japan); CBS-15603, CBS-15604, CBS-15605, CBS-15606, CBS-15607, CBS-15610 (clade III: Spain); COL-01, COL-02, COL-03 (clade IV: Colombia), 10-05-15-23, 10-05-15-19, 10-05-15-22 (clade IV: Venezuela). Also, 5 reference *Candida* strains were analyzed (*C. albicans* ATCC-90028,

Candida glabrata ATCC-90030, *Candida tropicalis* ATCC-200956, *C. lusitaniae* ATCC-135, and *Candida krusei* ATCC-6258), as well as one each of *Candida utilis*, *Candida rugosa*, *Candida orthopsilosis*, *Candida metapsilosis*, *Candida boidinii*, *Candida nivariensis*, *C. guillermondii*, *S. cerevisiae*, *Geotrichum capitatum*, *Trichosporon asahii*, *Kodamaea ohmeri*, and *Rodothorula glutinis* from clinical samples. Identification of all strains were confirmed by MALDI-TOF VITEK-MS[®] (bioMérieux, Marcy-l'Étoile, France) and sequencing of ITS 1-2 [19].

For monocultures, inoculum suspensions were prepared at 0.5 McFarland, adjusted to a concentration of 10^4 colony-forming units (CFU), and inoculated on three chromogenic media using a calibrated 5 µL pipette. All strains in monoculture were incubated at 37 °C for 48 h. Four mixed cultures containing a *C. auris* strain with 3 other non-*C. auris* species were performed to assess their ability to differentiate among *Candida* species. For all mixed cultures, inoculum suspensions were prepared at 0.5 McFarland and adjusted to 10^4 CFU. Then, the inoculum was adjusted to 10^3 CFU, and suspensions were inoculated on four chromogenic media using a calibrated 100 µL pipette. Mixed cultures were incubated at 37 °C and read at 24 h and 48 h by a clinical mycologist. Recording of color, size, and texture of colonies in each medium was performed both visually and using a magnifying lens (Nikon SMZ1500).

3. Results

In the four chromogenic media analyzed, the most specific morphology was observed at 48 h in all strains (Figure 1). Recovered colonies were presumptively identified according to the colony color described by the manufacturer.



Figure 1. Appearance of *Candida auris* colonies of the four clades described on CHROMagar[™] Candida Plus medium, CHROMagar[™] Candida medium, and HiCrome[™] Candida after 48 h incubation at 35 °C.

C. auris strains developed similar colonies in size, color, and texture in each medium. In CHROMagarTM Candida Plus and HiCromeTM Candida, *C. auris* strains from Colombia, Venezuela, India, Korea, and Japan exhibited blue-shaded colonies, while strains from Spain exhibited light pink with a bluish halo that diffuses into the surrounding agar. By analysis using HiCromeTM Candida, *C. auris* showed pale pink shades without a halo (see Figure 1).

On CHROMagarTM-Candida, all *C. auris* strains showed white and pink colonies at 24 and 48 h, respectively. Colonies of *C. auris*, *C. glabrata*, or *C. krusei* on CHROMagar Candida plus fluconazole showed a similar pink color after 48 h incubation; however, phenotypic differentiation was possible by the rough, paraffin-like texture or the intense purple color acquired by *C. krusei* and *C. glabrata*, respectively (Figures 1 and 2).



Figure 2. Appearance of colonies of most pathogenic *Candida* species on CHROMagar[™] Candida Plus medium, CHROMagar[™] Candida medium, and HiCrome[™] Candida after 48 h incubation at 35 °C.

After 24 h incubation on CHROMagar Candida plus fluconazole, *C. auris* showed smooth, convex, shiny, and light pink colonies (colony size 0.7–3 mm). After 48 h, they acquired an egg yolk morphology (white periphery and purple center). Although *C. auris*, *C. glabrata*, or *C. krusei* showed a similar pink color at 48 h incubation, phenotypic differentiation was possible due to the same characteristics observed for *C. krusei* and *C. glabrata* on CHROMagarTM-Candida. All *C. auris* strains but COL-01 and CBS-10913 were able to grow (these two isolates were susceptible to fluconazole), while *S. cerevisiae*, *C. albicans*, and *C. famata* showed inhibition of growth in this fluconazole-supplemented medium (Figures 2 and 3).

Regarding the colony morphology, the strains grown on CHROMagarTM Candida Plus developed individualized colonies, in contrast to the other three media where the growth is agglomerated at 48 h (Figures 2 and 3). On CHROMagarTM Candida and HiCromeTM Candida, after this incubation time, differentiation between the mauve/light pink species (*C. auris, C. ortopsilosis, C. metapsilosis, C. nivariensis, C. guillermondii*, and *C. lusitaniae*) was difficult. In CHROMagarTM Candida Plus, some species (*C. albicans, C. tropicalis, C. rugose*, and *K. ohmeri*) developed a similar metallic blue color, the distinction of which was difficult without significant changes with incubations over 48 h (Figure 2). Similar results were reported for species such as *C. auris, C. albicans, C. parapsilosis/orthopsilosis*, or *C. tropicalis* that presented a bluish color. However, differentiation was possible through the variety of shades or halos with different colors [17].

Convex colonies; shiny, metallic blue without halo 1 to 2 mm in size. pink without halo. 1 to 2 mm in size. shiny, no halo pre-2 to 3 mm in size. HROMaga Candida Pie Lenticular colonies; shiny, blue with pale margins. No halo present. 1 to 3 mm in size. Lenticular colonies; bright pink translucent appearan 1 to 2 mm in size. Lenticular colonies; light p translucent appearance. 1 to 2 mm in size. hromAga Candida hromAgar Candida Lenticular colonies; blue lighter edges. presence o bluish halo. 1 to 2 mm in size. Lenticular colonies; bright pink translucent and clear Lenticular colonies; light p translucent appearance. 1 to 2 mm in size HiCrome[†] Candida HiCrome Candida edges. 1 to 2 mm in size. Flat colonies; slightly bluish white, opaque and smooth, without halo. Convex colonies; shiny, ligh greenish blue. No halo Convex colonies; shiny, pink No halo present. 1 to 2 mm in size. HROMagar[®] Candida Plus without halo. 1 to 2 mm in size. present. 1 to 2 mm in size. Circular colonies; white waxy, rough-granular with jagged edges. 1 to 2 mm in size. Convex colonies; shiny pink, mauve with light edges, no halo present. 1 to 2 mm in size. Lenticular colonies; White translucent appearance. 1 to 2 mm in size. hromAga hromAga Condid Lenticular colonies; white translucent appearance. 1 to 2 mm in size. pink, translucent appearance. 1 to 2 mm in size. Lenticular colonies: White HICrome⁺ Candida HICrome¹ Candida HICrome¹ Candida translucent appearance. 1 to 2 mm in size. oitatum Shiny cobalt blue colonies; with a wide margin containing deep transverse Convex colonies; Shiny, vith dark centre ("Bull's eye"), no halo present. I to 3 mm in size. Off-white colonies; dull and smooth, with smooth edges 1 to 2 mm in size. CHROMagai Candida Plu fissures. 2 to 3 mm in size. Circular colonies; flat, waxy appearing, metallic blue, pseudomycelium visible on rim. 1-2 mm in size. enticular colonies; shinny nauve with lighter edges. Dry cottony colonies with frosted glass appearance 1 to 2 mm in size. ROMaga Candida hromAga Candida Pink halo present to 2 mm in size. White circular waxy flat colonies, with visible nies; blue, dry, and slightly granular in the cen light blue edges with deep radially arranged grooves. translucent appear 1 to 2 mm in size HiCrome¹ Candida HiCrome pseudomycelia. 1 to 2 mm in size. glutini , my convex colonies; bright orange, no halo present. 1 to 3 mm in size. Convex colonies; shiny blue with whitish edge and presence blue halo. 2 to 3 mm in size. HROMagai CHROMagar Candida Plu Requires 4 to 5 days of Tiny convex colonies; bright orange, no halo present. 1 to 3 mm in size. Lenticular colonies; bright pink translucent appearant 2 to 3 mm in size. Candin Requires 4 to 5 days of incubation. Tiny convex colonies; bright orange, no halo present. 1 to 3 mm in size. lighter edges and p of blue halo. 1 to 2 mm in size HiCrome[®] Candida Requires 4 to 5 days of incubation.

Figure 3. Colony appearance of five non-*Candida species* on CHROMagar[™] Candida Plus, CHROMagar[™] Candida, and HiCromeTM Candida after 48 h incubation at 35 °C.

In the evaluation of mixed cultures, on CHROMagar Candida plus fluconazole, *C. auris* showed smooth, convex, shiny, light purple colonies with an egg yolk morphology (white periphery and purple center). After 48 h of incubation, in cultures including *C. auris* in combination with *C. albicans*, *C. metapsilosis*, and *C. glabrata*, identification was difficult in the three tested media due to the whitish to pinkish coloration of the colonies. However, after 48 h incubation, *C. auris* colonies acquired a turquoise blue halo in CHROMagarTM Candida Plus that encouraged differentiation. All *C. albicans*, *C. metapsilosis*, and *C. glabrata* isolates were completely inhibited on this medium (Figure 4).



Figure 4. Appearance of *Candida species* and other yeast species colonies from mixed cultures on CHROMagar[™] Candida Plus, CHROMagar[™] Candida, CHROMagar[™] Candida plus fluconazole, and HiCromeTM Candida after 48 h incubation at 35 °C.

In mixed cultures, the number of colonies was similar for CHROMagarTM Candida Plus and CHROMagarTM Candida with some differences depending on the species assayed. CHROMagar Candida plus fluconazole was able to recover 97–100% of *C. auris* colonies before 48 h of incubation. Differentiation was not possible with HiCromeTM Candida in mixed cultures containing *C. auris* and *C. guillermondii* (Table 1).

Table 1. Comparison of CHROMagarTM Candida Plus, HiCromeTM Candida, CHROMagarTM Candida, and CHROMagar Candida supplemented with 32 mg/L of fluconazole for the recovery (CFU) of mixed cultures after 48 h of incubation at 35 °C.

Mixed Culture/Organism	Recovered Colonies (n)				
	CHROMagar TM Candida Plus	HiCrome TM Candida	CHROMagar TM Candida	CHROMagar Candida Plus Fluconazole	
Mix A					
C. auris	142	74	144	112	
C. glabrata	137	134	137	NG	
C. metapsilosis	32	85	30	NG	
C. albicans	65	65	74	NG	
Total	376	358	385	112	
<i>C. auris</i> recovered vs. total (%)	37.77	20.67	37.40	100	
Mix B					
C. auris	100	* 325	81	124	
C. guillermondii	213		166	NG	
Č. tropicalis	21	37	13	NG	
C. krusei	89	77	67	3	
Total	423	114	327	127	
<i>C. auris</i> recovered vs. total (%)	23.64	ND	24.77	97.64	

Mixed Culture/Organism	Recovered Colonies (n)			
	CHROMagar TM Candida Plus	HiCrome TM Candida	CHROMagar TM Candida	CHROMagar Candida Plus Fluconazole
Mix C				
C. auris	105	122	107	103
C. lusitaniae	108	107	108	NG
S. cerevisiae	39	35	28	NG
C. orthopsilosis	85	20	70	NG
Total	337	284	313	103
<i>C. auris</i> recovered vs. total (%)	31.16	42.96	34.19	100
Mix D				
C. auris	105	85	72	74
C. rugosa	57	57	35	NG
C. lusitaniae	96	49	130	NG
G. capitatum	37	24	19	NG
Total	295	215	256	74
<i>C. auris</i> recovered vs. total (%)	35.59	39.53	28.13	100

Table 1. Cont.

NG = Non-growth; ND = not determined; * uncountable.

4. Discussion

Rapid and reliable diagnosis of *C. auris* infection or colonization is essential for prompt and accurate antifungal therapy as well as for transmission control. Classical yeast identification techniques still available such as API-20C, AuxaColorTM 2, VITEK-2 YST, BD-PhoenixTM, and MicroScan misidentify *C. auris* [20,21]. To avoid the identification problems of available commercial techniques, the use of commercial media such as CHROMagar Candida[®] supplemented with Pal agar (sunflower seed extract) has been suggested as a simple and low-cost method to distinguish between *C. auris* isolates and the *C. haemulonii* complex [22].

New media such as HiCromeTM Candida and CHROMagarTM Candida Plus have been previously studied with identification rates between 0% and 33%, respectively [21]. Our results using CHROMagarTM Candida Plus agree with those described by other authors who reported good results in the detection of *C. auris* in both artificial cultures and surveillance samples after 36 h of incubation, respectively [17,21]. Identification of *C. auris* was difficult in the presence of white to mauve growing species at 24–36 h of incubation [17]. In our study, after 24 h of incubation, inhibition of susceptible species in CHROMagarTM Candida supplemented with 32 mg/L fluconazole allows a rapid recovery of *C. auris* colonies, although those were small and not colored enough.

In studies carried out with CHROMagarTM Candida Plus, it has been described that *C. auris, C. albicans, C. parapsilosis/orthopsilosis,* or *C. tropicalis* species can have a similar bluish color before 36 h [17,21,23]; in our study, we found similar results for *C. auris, C. orthopsilosis, C. metapsilosis, C. nivariensis,* and *C. glabrata* in this medium.

According to our results, in culture assays performed on CHROMagar Candida with 32 mg/L of fluconazole, the presence of *C. auris* in combination with other species of similar color was not limiting for its early identification, due to this medium selecting only fluconazole-resistant strains. Although *C. auris* isolation in this medium requires further confirmation, it allows for a prompt response that enables rapid isolation of patients and thus improved control of its transmission.

Currently, molecular methods or mass spectrophotometry (MALDI-TOF MS) are recommended for the definitive identification of *C. auris*. Among the molecular methods, sequencing of genetic loci (D1/D2, RPB1 and RPB2), rRNA internal transcribed spacers (ITS), and commercialized molecular technologies are the most commonly used [24,25]. These techniques are reliable but their availability outside reference laboratories is limited, especially in developing countries [26,27]. In these settings, the use of CHROMagarTM Candida agar with 32 mg/L of fluconazole as the primary culture medium in patients with suspected *C. auris* colonization allows fast presumptive identification of *C. auris* (24–48 h). This achievement is similar to that observed by other authors who used CHROMagarTM Candida Plus for the isolation and identification of *C. auris* from clinical samples in monoculture [17,21]. Although a presumptive identification can be achieved in 36 h, in all the media studied, a minimum incubation of 48 h is necessary to facilitate the differentiation of the species, and definitive identification must be done by spectrophotometric techniques such as MALDI-TOF or by molecular methodology sequencing the ITS regions [28–30].

In conclusion, the use of low-cost and accessible chromogenic media could allow the presumptive identification of *C. auris* infection or colonization in areas with a high prevalence of this species. This advantage is more evident with the use of CHROMagarTM Candida agar with 32 mg/L of fluconazole: this medium offers better performance for the identification of *C. auris* by inhibiting other species susceptible to this antifungal and selecting rare strains or those with high fluconazole MIC (>32 mg/L). Nevertheless, in all cases, definitive confirmation is required.

Author Contributions: Conceptualization, A.R.-G., I.S.-G. and J.P.; Literature search, V.G.-B., M.G.-H. and E.V.-G.; Writing—original draft preparation, A.R.-G., J.P. and V.G.-B.; editing, S.G.A. and P.W.J.d.G.; Elaboration of Tables and Figures, A.R.-G., I.S.-G. and J.M.P.-R.; Supervision, A.R.-G. All authors have read and agreed to the published version of the manuscript.

Funding: This research received no external funding.

Institutional Review Board Statement: Not applicable.

Informed Consent Statement: Not applicable.

Data Availability Statement: Not applicable.

Acknowledgments: We thank Jacques Meis (Centre of Expertise in Mycology Radboudumc/Canisius-Wilhelmina Hospital, Nijmegen, The Netherlands), and Indira Berrio (Medical and Experimental Mycology Group, Biological Research Corporation, Medellín, Colombia) for the supply of clinical isolates recovered from Indian, Korean, Japanese, Venezuelan and Colombian isolates.

Conflicts of Interest: The authors declare no conflict of interest.

References

- Satoh, K.; Makimura, K.; Hasumi, Y.; Nishiyama, Y.; Uchida, K.; Yamaguchi, H. Candida aurissp. nov., a novel ascomycetous yeast isolated from the external ear canal of an inpatient in a Japanese hospital. *Microbiol. Immunol.* 2009, 53, 41–44. [CrossRef]
- Jeffery-Smith, A.; Taori, S.K.; Schelenz, S.; Jeffery, K.; Johnson, E.M.; Borman, A.; Manuel, R.; Brown, C.S. Candida auris: A review of the literature. *Clin. Microbiol. Rev.* 2017, 31, e00029-17. [CrossRef]
- Ruiz-Gaitán, A.; Moret, A.M.; Tasias-Pitarch, M.; Aleixandre-López, A.I.; Martínez-Morel, H.; Calabuig, E.; Salavert-Lletí, M.; Ramírez, P.; López-Hontangas, J.L.; Hagen, F.; et al. An outbreak due to *Candida auris* with prolonged colonisation and candidaemia in a tertiary care European hospital. *Mycoses* 2018, *61*, 498–505. [CrossRef] [PubMed]
- Welsh, R.M.; Bentz, M.L.; Shams, A.; Houston, H.; Lyons, A.; Rose, L.J.; Litvintseva, A.P. Survival, Persistence, and Isolation of the Emerging Multidrug-Resistant Pathogenic Yeast *Candida auris* on a Plastic Health Care Surface. *J. Clin. Microbiol.* 2017, 55, 2996–3005. [CrossRef] [PubMed]
- Garcia-Bustos, V.; Cabanero-Navalon, M.D.; Ruiz-Saurí, A.; Ruiz-Gaitán, A.C.; Salavert, M.; Tormo, M.Á.; Pemán, J. What Do We Know about *Candida auris*? State of the Art, Knowledge Gaps, and Future Directions. *Microorganisms* 2021, 9, 2177. [CrossRef] [PubMed]
- 6. Centers for Disease Control. Antibiotic Resistance Threats in the United States, 2019; CDC: Atlanta, Georgia, 2019.
- Wang, X.; Bing, J.; Zheng, Q.; Zhang, F.; Liu, J.; Yue, H.; Tao, L.; Du, H.; Wang, Y.; Wang, H.; et al. The first isolate of *Candida auris* in China: Clinical and biological aspects. *Emerg. Microbes Infect.* 2018, 7, 1–9. [CrossRef]
- Bidaud, A.; Chowdhary, A.; Dannaoui, E. *Candida auris*: An emerging drug resistant yeast—A mini-review. J. Med. Mycol. 2018, 28, 568–573. [CrossRef] [PubMed]
- Prakash, A.; Sharma, C.; Singh, A.; Singh, P.K.; Kumar, A.; Hagen, F.; Govender, N.; Colombo, A.; Meis, J.; Chowdhary, A. Evidence of genotypic diversity among *Candida auris* isolates by multilocus sequence typing, matrix-assisted laser desorption ionization time-of-flight mass spectrometry and amplified fragment length polymorphism. *Clin. Microbiol. Infect.* 2016, 22, 277.e1–277.e9. [CrossRef]

- 10. Perry, J.D. A Decade of Development of Chromogenic Culture Media for Clinical Microbiology in an Era of Molecular Diagnostics. *Clin. Microbiol. Rev. Am. Soc. Microbiol.* **2017**, *30*, 449–479. [CrossRef]
- Fasciana, T.; Cortegiani, A.; Ippolito, M.; Giarratano, A.; Di Quattro, O.; Lipari, D.; Graceffa, D.; Giammanco, A. *Candida auris*: An Overview of How to Screen, Detect, Test and Control This Emerging Pathogen. *Antibiotics* 2020, *9*, 778. [CrossRef]
- Das, S.; Singh, S.; Tawde, Y.; Chakrabarti, A.; Rudramurthy, S.M.; Kaur, H.; Shankarnarayan, S.A.; Ghosh, A. A Selective Medium for Isolation and Detection of *Candida auris*, an Emerging Pathogen. J. Clin. Microbiol. 2021, 59, e00326-20. [CrossRef] [PubMed]
- 13. Orenga, S.; James, A.L.; Manafi, M.; Perry, J.D.; Pincus, D.H. Enzymatic substrates in microbiology. J. Microbiol. Methods 2009, 79, 139–155. [CrossRef] [PubMed]
- 14. De Angelis, G.; Menchinelli, G.; Torelli, R.; De Carolis, E.; Posteraro, P.; Sanguinetti, M.; Posteraro, B. Different detection capabilities by mycological media for Candida isolates from mono- or dual-species cultures. *PLoS ONE* **2020**, *15*, e0226467. [CrossRef]
- 15. Tamura, T.; Alshahni, M.M.; Makimura, K. Evaluation of CHROMagarTM Candida Plus chromogenic agar for the presumptive identification of *Candida auris*. *Microbiol*. *Immunol*. **2022**, *66*, 292–298. [CrossRef] [PubMed]
- 16. Borman, A.M.; Fraser, M.; Johnson, E.M. CHROMagarTM Candida Plus: A novel chromogenic agar that permits the rapid identification of *Candida auris*. *Med. Mycol.* **2020**, *11*, 1–6. [CrossRef]
- Mulet Bayona, J.V.; Salvador García, C.; Tormo Palop, N.; Gimeno Cardona, C. Evaluation of a novel chromogenic medium for *Candida* spp. identification and comparison with CHROMagarTM Candida for the detection of *Candida auris* in surveillance samples. *Diagn. Microbiol. Infect. Dis.* 2020, *98*, 115168. [CrossRef]
- Bayona, J.V.M.; García, C.S.; Palop, N.T.; Martín, A.V.; Padrón, C.G.; Rodríguez, J.C.; Pemán, J.; Cardona, C.G. Novel Chromogenic Medium CHROMagarTM Candida Plus for Detection of *Candida auris* and Other *Candida* Species from Surveillance and Environmental Samples: A Multicenter Study. *J. Fungi* 2022, *8*, 281. [CrossRef]
- 19. Bellemain, E.; Carlsen, T.; Brochmann, C.; Coissac, E.; Taberlet, P.; Kauserud, H. ITS as an environmental DNA barcode for fungi: An in silico approach reveals potential PCR biases. *BMC Microbiol.* **2010**, *10*, 189. [CrossRef]
- Chowdhary, A.; Sharma, C.; Meis, J.F. Candida auris: A rapidly emerging cause of hospital-acquired multidrug-resistant fungal infections globally. PLoS Pathog. 2017, 13, e1006290. [CrossRef]
- de Jong, A.W.; Dieleman, C.; Carbia, M.; Tap, R.M.; Hagen, F. Performance of Two Novel Chromogenic Media for the Identification of Multidrug-Resistant *Candida auris* Compared with other Commercially Available Formulations. *J. Clin. Microbiol.* 2021, 59, e03220-20. [CrossRef]
- Kumar, A.; Sachu, A.; Mohan, K.; Vinod, V.; Dinesh, K.; Karim, S. Simple low cost differentiation of *Candida auris* from *Candida haemulonii* complex using CHROMagar Candida medium supplemented with Pal's medium. *Rev. Iberoam Micol.* 2017, 34, 109–111. [CrossRef] [PubMed]
- Sasoni, N.; Maidana, M.; Latorre-Rapela, M.G.; Morales-Lopez, S.; Berrio, I.; Gamarra, S.; Garcia-Effron, G. Candida auris and some Candida parapsilosis strains exhibit similar characteristics on CHROMagarTMCandida Plus. Med. Mycol. 2022, 60, myac062. [CrossRef] [PubMed]
- 24. Sharma, C.; Kumar, N.; Meis, J.F.; Pandey, R.; Chowdhary, A. Draft Genome Sequence of a Fluconazole-Resistant *Candida auris* Strain from a Candidemia Patient in India. *Genome Announc.* **2015**, *3*, e00722-15. [CrossRef] [PubMed]
- 25. Chowdhary, A.; Sharma, C.; Duggal, S.; Agarwal, K.; Prakash, A.; Singh, P.K.; Jain, S.; Kathuria, S.; Randhawa, H.S.; Hagen, F.; et al. New Clonal Strain of *Candida auris*, Delhi, India. *Emerg. Infect. Dis.* **2013**, *19*, 1670–1673. [CrossRef]
- Berkow, E.L.; Lockhart, S.R. Fluconazole resistance in Candida species: A current perspective. *Infect. Drug Resist.* 2017, 10, 237–245. [CrossRef]
- Lockhart, S.R.; Etienne, K.A.; Vallabhaneni, S.; Farooqi, J.; Chowdhary, A.; Govender, N.P.; Colombo, A.L.; Calvo, B.; Cuomo, C.A.; Desjardins, C.A.; et al. Simultaneous Emergence of Multidrug-Resistant *Candida auris* on 3 Continents Confirmed by Whole-Genome Sequencing and Epidemiological Analyses. *Clin. Infect. Dis.* 2017, 64, 134–140. [CrossRef]
- Pino-Calm, B.; de Artola, D.G.M.; Gil-Campesino, H.; Alcoba-Flórez, J. Restriction fragment length polymorphism de las regiones ITS1-ITS2 como método para identificar Candida auris. *Rev. Iberoam Micol.* 2018, 35, 167–168. [CrossRef]
- Bao, J.R.; Master, R.N.; Azad, K.N.; Schwab, D.A.; Clark, R.B.; Jones, R.S.; Moore, E.C.; Shier, K.L. Rapid, Accurate Identification of *Candida auris* by Using a Novel Matrix-Assisted Laser Desorption Ionization–Time of Flight Mass Spectrometry (MALDI-TOF MS) Database (Library). J. Clin. Microbiol. 2018, 56, e01700-17. [CrossRef]
- 30. Mahmoudi, S.; Afshari, S.A.K.; Gharehbolagh, S.A.; Mirhendi, H.; Makimura, K. Methods for identification of Candida auris, the yeast of global public health concern: A review. *J. Med. Mycol.* **2019**, *29*, 174–179. [CrossRef]

Disclaimer/Publisher's Note: The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.